Description

MultiCore® Flow is a dual-curing, fluoride-containing, radiopaque composite for core build-ups. It cures chemically without the use of light. Light-curing is optional. MultiCore Flow is available in a double-barrel cartridge or a double-push syringe with automixing tips for fabricating core build-ups using matrices.

Composition

The monomer matrix consists of dimethacrylate (29 wt %).

The inorganic fillers are barium glass, ytterbiumtrifluoride, Ba-Al-fluorosilicate glass and highly dispersed silicon dioxide (70 wt %).

Additional contents are catalysts, stabilizers and pigments (1 wt %).

The total content of inorganic fillers is 46 vol %. The particle size ranges from 0.04 to 25 µm.

Shades

MultiCore is available in four shades: The dentin shades Light (A1/B1), Medium (A2/A3) and White are particularly suitable for use under all-ceramic restorations, and Blue is recommended for use under metal or porcelain-fused-to-metal crowns.

Working time

90 to 120 s (at 37°C/98.6°F)

Mixing ratio

The base and catalyst pastes of MultiCore Flow are mixed at a ratio of 1:1.

MultiCore Flow is mixed by pressing the pastes through the static mixing tip.

Indications

- Core build-up of vital and non-vital teeth
- Adhesive cementation of glass fibre-reinforced endodontic posts

Contraindications

The use of MultiCore Flow is contraindicated

- if the operating field cannot be kept sufficiently dry, or if the stipulated application technique cannot be used;
- if the patient is known to be allergic to any of the ingredients in MultiCore Flow.

Side effects

- In individual cases, components of MultiCore Flow may lead to sensitization. MultiCore Flow should not be used in such cases.
- In vital teeth, dentinal areas close to the pulp should be covered with a suitable liner.

Interactions

Substances containing eugenol/oil of cloves may inhibit the polymerization of composite materials. Consequently, application of such materials together with MultiCore Flow must be avoided.

MultiCore® Flow

EN Instructions for Use	FI Käyttöohjeet
– Two-component core build-up material	– Kaksikomponenttinen pilarimateriaali
DE Gebrauchsinformation	NO Bruksanvisning
– Zweikomponenten Stumpfbaumaterial	– Tokomponent konusoppbygging materiale
FR Mode d'emploi	NL Gebruiksaanwijzing
– Matériau de reconstitution de faux-moignons	– Tweecomponenten stompopbouw materiaal
IT Istruzioni d'uso	EL Οδηγίες Χρήσεως
– Composito bicomponente per la ricostruzione del moncone	– Υλικό αναύστασης κολοβωμάτων δύο συστατικών
ES Instrucciones de uso	TR Kullanma Talimatı
– Material de muñones bicomponente	– İki bileşenli koruyarım materyali
PT Instruções de Uso	RU Инструкция по применению
– Material bi-componente para confecção de núcleos	– Двухкомпонентный материал для изготовления культовых вкладок

I. Fabrication of a reconstructive build-up

Apply a dentin bonding agent to the clean dentin surfaces. Ivoclar Vivadent recommends using either Adhese® Universal, Syntac® or Excite® F as the bonding agent when MultiCore Flow is used in the light-cure mode. If no light is used, the application of a dual-curing adhesive, e.g. Excite® F DSC, is recommended.

These Instructions describe the procedure with Adhese Universal. Please refer to the instructions of the respective bonding agents for further information.

1. Isolation

Adequate relative or absolute isolation using suitable auxiliaries, such as OptraGate® or OptraDam® Plus, is required.

2. Pulp protection / Base

In very deep cavities, areas close to the pulp should be selectively coated with a calcium hydroxide liner (e.g. ApexCal®) and subsequently covered with pressure-resistant cement (glass ionomer cement, e.g. Vivaglass® Liner).

3. Conditioning with phosphoric acid gel (optional)

Apply phosphoric acid gel (e.g. Total Etch) to the prepared enamel and then flow the etchant onto the prepared dentin. The etchant should be left to react on the enamel for 15–30 s and on the dentin for 10–15 s. Then rinse thoroughly with a vigorous stream of water for at least 5 s and dry with compressed air until the etched enamel surfaces appear chalky white.

4. Application of the adhesive

- Starting with the enamel, thoroughly coat the tooth surfaces to be treated with Adhese Universal.
- The adhesive must be scrubbed into the tooth surface for at least 20 s. This time must not be shortened. Applying the adhesive on the tooth surface without scrubbing is inadequate.
- Disperse Adhese Universal with oil- and moisture-free compressed air until a glossy, immobile film layer results.
- Light-cure Adhese Universal for 10 s using a light intensity of $\geq 500 \text{ mW/cm}^2$.
- Contamination by blood or saliva must be avoided after the application of the bonding agent.

Fig. 1

Fig. 2

Fig. 3

Fig. 4

5. Insertion of the cartridge

(If you are using an automix syringe, please proceed to point c)

- Press the black release lever located below the plunger on the back of the dispenser and pull plunger back as far as possible (Fig. 1).
- Lift the cartridge lock and insert the cartridge with the "V" shape on the cartridge base turned downwards. Lower the cartridge lock (Fig. 2).
- Remove the sealing cap by turning it 1/4 turn counter clockwise.
- It is **very important** to clear or bleed the cartridge or syringe **prior** to applying the mixing tip. Gently press the dispenser lever until both components (base and catalyst) begin to flow out evenly onto a mixing pad.
- For each application, place a new mixing tip on the cartridge/syringe. Push the tip completely down until the notch of the mixing tip is aligned with the notch on the housing of the cartridge/syringe. Secure the mixing tip in place by gripping the coloured base, and not the mixing tip, and by turning it 1/4 turn clockwise (Fig. 3). As MultiCore Flow will cure in the used mixing tip, it may serve as a seal for the contents of the cartridge/syringe until needed once again (replace with a new tip just before the next use).
- If necessary apply an intra-oral tip to the mixing tip (Fig. 4). Larger amounts of material can be extruded if the intra-oral tip is shortened by approx. 1 mm.

6. Application of MultiCore Flow

- Place matrix band around the prepared tooth.
- Apply MultiCore Flow directly to the cavity with the intra-oral tip in place.
- If no light is used, wait for 4–5 minutes until the chemical curing process is complete before removing the matrix. If the material is light-cured, it may be ground immediately after completing the curing cycle.
- Optionally, MultiCore Flow may be cured from the occlusal for approx. 10 s (approx. $1,100 \text{ mW/cm}^2$, e.g. Bluephase Style) ensuring a minimal distance between the light emission window and the build-up. Remove the matrix band. If necessary, cure incompletely polymerized areas for another 10 s.
- Use rotary instruments to complete the crown preparation.

II. Adhesive cementation of endodontic posts made of fibre-reinforced composite and core build-ups

1. Pretreatment of the root canal and the coronal dentin

- A dry operating field must be ensured. Ideally, a rubber dam should be used, e.g. OptraDam Plus.
- Prior to cementing endodontic posts, clean the root canal thoroughly to remove any residue of root canal sealers (Residue of eugenol-based sealers may inhibit the polymerization of the luting composite).
- Condition the root canal, the coronal dentin and the enamel surfaces according to the instructions of the manufacturer of the dentin bonding agent used (eg or Excite® F DSC).
- After applying the dentin bonding agent, make sure to remove all excess material from the root canal using paper points.
- Do not light-cure the bonding agent.
- Place matrix band around the tooth.

2. Adhesive cementation

- Clean and condition post surface according to the manufacturer's instructions.
- Coat the endodontic post, which has been prepared according to the instructions of the post manufacturer, with the mixed MultiCore

Flow. Additionally, MultiCore Flow may be applied to the root canal by means of the Intra Canal Tips. It is not recommended to use a lentulo spiral for loading MultiCore Flow into the root canal, as this will lead to the intermixing of bonding agent excess and MultiCore Flow. As a result, the working time of MultiCore Flow will be reduced, so that it may be impossible to push the post into the desired position within the root canal.

- Excess cement will be displaced when the post is inserted.
- Distribute excess cement across the occlusal preparation surface for total coverage.
- Subsequently, light-cure MultiCore Flow for 20 s (at a light intensity of approx. 1,100 mW/cm², e.g. Bluephase Style). In the process, hold the post in position with the curing light. If non-translucent, opaque endodontic posts are cemented, wait for the self-curing mechanism to take effect.

3. Core build-up

- Apply and cure MultiCore Flow as described above (see “I. Fabrication of a reconstructive build-up”).

Additional notes

MultiCore Flow should be at room temperature when processed. At lower temperatures (e.g. immediately after removal from the refrigerator) it may be difficult to extrude the material.

Warnings

- Unpolymerized MultiCore Flow should not come in contact with skin, mucous membrane, or eyes. Unpolymerized MultiCore Flow may have a slight irritating effect and may lead to a sensitization against methacrylates.
- Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

Storage and shelf life

- MultiCore Flow has to be stored at 2–8°C/36–46°F.
- Do not use MultiCore Flow after the indicated date of expiration.
- Leave the used mixing tip on the MultiCore Flow automix syringe/cartridge.
- Expiry date: See note on syringe, cartridge and/or packaging.

Keep out of the reach of children!

For use in dentistry only.

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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Ivoclar Vivadent S.A. de C.V.
Ivoclar Vivadent Ltd.
Ivoclar Vivadent LLC
Ivoclar Vivadent Marketing Ltd.
Ivoclar Vivadent AB