

Prime&Bond active™

Moisture levels vary. Our bond won't.



# Active Control for reliable restorations

HEMA-free  
BPA-free



Prime&Bond active™ introduces the new, patented Active-Guard™ Technology.


Our formulation provides well balanced hydrophobic and hydrophilic properties to ensure complete coverage and penetration at varying moisture levels. So even if dentin is overly wet or dry, you get a reliable bond that actively guards against the leading causes of adhesive failure. Prime&Bond active™ provides:

- All etching methods, all indications for more versatility
- Active moisture control
- Active spreading
- Low film thickness
- Active tubule penetration

Selective enamel etching – ideal for class II restorations.  
Get high bond strength in dentin, a tight seal at the enamel margin and minimal etchant exposure.

Prime&Bond active™ provides active control and consistent results – for virtually no post-operative sensitivity.

RESEARCH SHOWS THAT ON AVERAGE, **10%-15%** OF POSTERIOR COMPOSITE RESTORATIONS **RESULT IN POST-OPERATIVE SENSITIVITY.**<sup>1</sup>



no post-operative sensitivity  
post-operative sensitivity

Source: User survey 2015<sup>2</sup>

<sup>1</sup> Usman et al, Sensitivity in composite restorations, Pakistan Oral & Dental Journal Vol 34, No. 3 (September 2014); Berkowitz et al, Postoperative Hypersensitivity in class I Resin-based Composite Restorations, Compend Contin Educ Dent. 2009 ; 30(6): 356-363; Haller, Die Postoperative Hypersensibilität, zm 99, Nr. 6a, 13.03.2009 S. 44-51; Briso et al, Clinical Assessment of Postoperative Sensitivity in Posterior Composite Restorations, Operative Dentistry, 2007, 32-5, 421-426.

<sup>2</sup> Data on file.

<sup>3</sup> Not registered trademarks of Dentsply Sirona.

## ACTIVE CAVITY SPREADING

# COMPLETE COVERAGE

Inherent surface tension prevents other adhesives from spreading easily across the surface of the cavity. Dry spots that remain unbonded can weaken the restoration and may lead to post-operative sensitivity and microleakage.

Prime&Bond active™ is formulated to spread actively, minimizing the risk of inadvertently leaving dry spots.



Prime&Bond active™

A single drop, placed using a micro brush without further manipulation, naturally spreads across most of the surface.



Scotchbond Universal³

Other universal adhesives resist spreading.



iBOND Universal³


Other universal adhesives resist spreading.

## ACTIVE UNIFORMITY


# LOW FILM THICKNESS

Thick adhesives can have a tendency to pool in the corners of the proximal box of a class II restoration. These adhesive pools can show up as translucent areas on a radiograph, which can easily be misdiagnosed as a void, gap, or secondary decay, leading to unnecessary replacement. Thicker adhesives may also change crown prep geometry, impacting crown fit.


Prime&Bond active™ self-levels to actively create a thin, uniform adhesive layer.



Prime&Bond active™



Scotchbond Universal³



Clearfil™ Universal³


Source: "Film Thickness FE-SEM evaluation," Andre F. Reis, DDS, MS, PhD.

ACTIVE STRENGTH AND SIMPLICITY

# OPTIMIZED WITH CALIBRA® CERAM CEMENT

There's no need to apply an activator when Prime&Bond active™ is used with Calibra® Ceram cement. With high bond strength and easy gel cleanup, this is the ideal combination of all ceramic restorations.

## 6 Minute Dentin Shear Bond Strength



Test was performed following manufacturer's Instructions for Use. Scotchbond Universal³ was light cured as indicated for best performance.



ACTIVE ADVANTAGE

## CONVENIENCE AND CONTROL

- Flip-top cap for one-handed opening
- Precise dispensing for minimal waste and mess
- Low viscosity for easy application
- Up to 30 minutes working time in a closed CliXdish™, reducing waste when filling multiple cavities



ACTIVE ANSWERS

To learn more about Prime&Bond active™, contact your Dentsply Sirona representative.


<sup>3</sup> Not registered trademarks of Dentsply Sirona.

# ACTIVE MOISTURE CONTROL TOO WET OR TOO DRY: NO PROBLEM.

Most adhesives are comparatively hydrophobic, tending to separate from water. With too much water on the prepared surface, the adhesive can't achieve uniform coverage, leaving gaps that may weaken the bond and may cause post-operative sensitivity. If dentin is overly dried, collagen fibers can collapse and prevent effective sealing.

Prime&Bond active™ overcomes the surface tension of water and spreads across dentin and into dentinal tubules to form an even, homogeneous layer. Upon air-drying, the solvent and water evaporate together, leaving a consistent coating of adhesive across the entire surface to ensure a strong bond.

## Shear Bond Strength



Source: Internal Data. Shear Bond Strength of Universal Adhesives on dentin after water contamination. Data on file.

## ACTIVE MIXING



### Prime&Bond active™

To ensure a homogeneous application, water and solvent evaporate together, leaving an even, consistent adhesive layer.



### Scotchbond Universal³

Cannot overcome the surface tension of water. Adhesive and water remain separated.



### iBOND Universal³

Cannot overcome the surface tension of water. Adhesive and water remain separated.

PRIME&BOND ACTIVE™

## Ordering Information

**Mini Refill** 606.67.340  
1 x 2.5 ml Bottle  
IFU+ITG

**Standard Refill** 606.67.341  
1 x 4 ml Bottle  
IFU+ITG

**SUD Refill** 606.67.343  
80 x Single Unit Dose  
50 x Applicator Tip  
IFU+ITG

**Eco Refill** 606.67.342  
3 x 4 ml Bottle  
IFU+ITG

**SUD Eco Refill** 606.67.345  
240 x Single Unit Dose  
100 x Applicator Tip  
IFU+ITG

**Intro Kit** 606.67.347  
1 x 4 ml Bottle  
1 x CliXdish™  
50 x Applicator Tip  
IFU+ITG

**CliXdish™ Refill** 606.67.346  
3 x CliXdish™



Standard Refill



SUD Refill



Eco Refill



Intro Kit



CliXdish™ Refill

For more information visit [www.dentsplysirona.com](http://www.dentsplysirona.com).

Find us on  [www.facebook.com/dentsplysirona.restorative](https://www.facebook.com/dentsplysirona.restorative)

K79200229-04  
© Dentsply Sirona 03/2018

THE DENTAL  
SOLUTIONS  
COMPANY™

 Dentsply  
Sirona