Tetric® plus Flow

[en] Instructions for Use

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Patient target group

- Patients with permanent teeth
- Patients with deciduous teeth

Intended users

- Dentists

Special training

None

Use

For dental use only.

Description

Tetric® plus Flow is a flowable, light-curing, radiopaque composite resin ($\geq 200\%$ Al*) for the direct restorative treatment of anterior and posterior teeth and for restoring occlusal surfaces (according to ISO 4049:2019 Type 1, Class 2, Group 1).

*Natural dentin has a radiopacity equivalent to that of 100% aluminium. Natural enamel has a radiopacity equivalent to that of 200% aluminium.

Indications

- Missing tooth structure in anterior and posterior teeth

Types of restorations:

- Initial layer up to 4 mm in Class I and II restorations*
- Small direct restorations in permanent teeth (not occlusion-bearing)*
- Restorations of deciduous teeth
- Liner under direct restorations of all cavity classes
- Class V restorations
- Extended fissure sealings
- Repair of composite / ceramic restorations
- Sealing of implant screw channels

*In case of occlusion-bearing restorations, apply a final layer using a sculptable composite material.

Contraindications

The use of this product is contraindicated if the patient is known to be allergic to any of its ingredients.

Limitations of use

The product must not be used in the following cases:

- Restorations of load-bearing cusps and incisal edges of permanent teeth
- After the indicated expiry date
- Insufficiently dry working field
- Failure to observe the stipulated working technique
-  Reuse of Cavifil and application tip for several patients, to avoid potential cross-contaminations

Side effects

- In rare cases, components of the product may lead to sensitization
- Irritation of the pulp. To avoid irritation of the pulp: Apply pulp / dentin protection in areas close to the pulp.

Interactions

- Phenolic substances (e.g. eugenol) and disinfectants with an oxidative effect impair the curing process.
- Discolouration may occur when the final restoration comes in contact with cationic mouthwashes, plaque disclosing agents or chlorhexidine.

Clinical benefit

- Restoration of esthetics
- Reconstruction of chewing function

Composition

Barium glass, co-polymer, mixed oxide ($\text{SiO}_2/\text{ZrO}_2$), BPEMA, MOMA, Bis-GMA, Bis-EMA, DCP, UDMA, ytterbium fluoride, AFCT.
Total content of inorganic fillers: 65 wt% / 50–51 vol%
Particle size of the inorganic fillers: 0.02–3.0 μm

2 Application

- ✓ The product should be at least at ambient temperature.
- ✓ Application tips for the product line are ready-to-hand.
When using other application tips, the material may be difficult to extrude.

2.1 Placing the direct restoration

Determining the tooth shade

1. Clean the teeth professionally.
2. Determine the tooth shade with the tooth still moist using the Tetric plus shade guide.

Table 1:

Shade correspondence				
Tetric® plus	A2 plus	A3 plus	A3.5 plus	Bleach plus
VITA® classical shades	A1, A2, B2, C1, D2	A3, C2, C3, D3, D4	A3.5, A4, B3, B4, C4	B1, Bleach

* No registered trademark of Ivoclar.

Preparing the cavity according to the principles of the adhesive technique

1. Remove caries / defective restoration, preserving as much of the healthy tooth structure as possible.
For caries-free cervical defects after optional preparation: Clean the cavity with fluoride-free cleaning paste, a rubber cup or a rotary brush.
2. Round out internal edges.
3. Anterior region: Bevel the enamel margins.
Posterior region: Round out the enamel margins.
4. Rinse the cavity with water spray.
5. Dry the cavity with moisture- and oil-free compressed air.

Ensuring isolation

1. Provide relative or absolute isolation.

Placing the matrix

1. Place the matrix and secure it with wedges.

Lining the cavity close to the pulp

1. Selectively coat areas close to the pulp with a calcium hydroxide liner. Observe the manufacturer's instructions.
2. Selectively apply a pressure resistant glass ionomer cement on the calcium hydroxide liner. Observe the manufacturer's instructions.

Conditioning the cavity

1. Apply the bonding agent. Observe the manufacturer's instructions.

Applying the composite resin

1. Attach a new application tip to the syringe. / Insert Cavifil into the applicator.
2. Cover the syringe with a protective sleeve. Observe the manufacturer's instructions.
3. Apply the material to the cavity while observing the maximum layer thickness (see table 2,) and keep the application tip in the material during application to avoid air entrapments.
4. Distribute the material in the cavity to ensure tight marginal adaptation.

Table 2:

Shade	Maximum layer thickness
A2 plus, A3 plus, A3.5 plus, Bleach plus	4 mm*

* Apply the shade A3.5 plus in a maximum layer thickness of 3.5 mm when light-curing for 3 s or 5 s.

* The recommended increment thickness is based on hardness profile measurements.

Curing the composite resin

NOTICE! Avoid direct exposure of the gingiva or mucous membrane to the light emitted by the curing light.

1. Completely cure the restoration with the blue light of a curing light (400–500 nm) while observing the curing times (table 3). Observe the manufacturer's instructions.
2. If the diameter of the light guide is smaller than the restoration: Perform an overlapping polymerization.
3. If the light guide could not be positioned ideally: Light-cure the restoration again.
4. When using a metal matrix and a third-party curing light: Repeat the polymerization of the buccal and the lingual/palatal surfaces.
5. In case of adjustments before polishing: Apply Tetric plus Flow directly to the polymerized material.

Table 3:

Light intensity	Curing time
500–900 mW / cm ²	20 s
1,000–1,300 mW / cm ² *	10 s
1,800–2,200 mW / cm ²	5 s
2,700–3,300 mW / cm ² **	3 s

* Polymerize restorations in deciduous teeth at light intensities of $\leq 1,300$ mW / cm².

** The 3sCure program of the Bluephase PowerCure must only be used for light-curing posterior Tetric plus Flow restorations from the occlusal surface.

** The 3sCure polymerization program of the Bluephase PowerCure must not be used in the case of caries profunda and for restorations close to the pulp.

** Light-cure Tetric plus Flow for 3 seconds only with the Bluephase PowerCure in the 3sCure polymerization program.

Finishing the restoration

1. Finish the restoration with tungsten carbide or diamond polishers.
2. Check the occlusion and articulation.
3. In case of premature contacts or undesired articulation paths: Adjust the occlusion and articulation with tungsten carbide or diamond polishers.
4. Polish the restoration.
5. In case of adjustments after polishing: Roughen the restoration and wet it with adhesive material.

2.2 Repairing composite / ceramic restorations

1. Roughen the surface to be repaired with a diamond bur.
2. Clean the surface with water spray.
3. Dry the surface with moisture- and oil-free compressed air.
4. In case of defective ceramic restorations: Apply ceramic primer to the surface. Observe the manufacturer's instructions.
5. Apply the bonding agent. Observe the manufacturer's instructions.
6. Apply the composite resin. (See section "Applying the composite resin")
7. Cure the composite resin. (See section "Curing the composite resin")
8. Finish the restoration. (See section "Finishing the restoration")

2.3 Sealing an implant screw channel

1. Prepare the implant screw channel. Observe the manufacturer's instructions.
2. Apply the ceramic primer. Observe the manufacturer's instructions.
3. Apply the bonding agent. Observe the manufacturer's instructions.
4. Apply the composite resin. (See section "Applying the composite resin")

5. Cure the composite resin. (See section "Curing the composite resin")
6. Finish the restoration. (See section "Finishing the restoration")

3 Safety information

- The product has been developed solely for use in dentistry. Process according to the Instructions for Use.
- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/ Liechtenstein, www.ivoclar.com and your responsible competent authority.
- Product contains methacrylates. Irritation of eyes, respiratory tract and skin as well as sensitization are possible. Contact of monomer and uncured material with skin / mucous membrane and eyes must be avoided (commercially available medical gloves do not offer sufficient protection).
- Wear a face mask, protective gloves and goggles.
- Syringes or Cavifils must not be disinfected with oxidizing disinfection agents.
- Exposure to light causes premature polymerization. Close product immediately after use.
- Before use, visually inspect the packaging and the product for damage. In case of any doubts, please contact Ivoclar Vivadent AG or your local dealer.

Supporting documents

Document	To be found in:
Current version of Instructions for Use	www.ivoclar.com/eIFU
Structure of Instructions for Use and Warnings	www.ivoclar.com/eIFU
Explanation of symbols	www.ivoclar.com/eIFU
Safety Data Sheet (SDS)	www.ivoclar.com
Summary of Safety and Clinical Performance – SSCP	https://ec.europa.eu/tools/eudamed
Basic UDI-DI: 76152082AFILL005JJ	

Disposal information

Dispose contaminated products, remaining stocks or removed restorations according to the corresponding national legal requirements.

Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks.

The following known residual clinical risks exist:

- Failure of the adhesive bond (loss of the filling)
- Postoperative sensitivity
- Air entrapments during placement of the filling
- Chipping, fractures
- Ingestion of material

4 Shelf life and storage

- Storage temperature 2–28 °C

5 Additional information

Keep material out of the reach of children!

The manufacturer accepts no liability for damages resulting from a failure to observe the Instructions for Use. Furthermore, the user is responsible for testing the products for their suitability for any purpose not explicitly stated in the Instructions.

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