# Multilink® Automix

## [en] Instructions for Use

Polymer-based dental luting material (intraoral dual-curing)  
Self-curing and self-etching primer

Rx ONLY

CE 0123

Date information prepared:  
2023-08-31 / Rev. 1

Manufacturer:  
Ivoclar Vivadent AG  
Bendererstrasse 2  
9494 Schaan/Liechtenstein  
[www.ivoclar.com](http://www.ivoclar.com)

ivoclar

English

### 1 Intended use

#### Intended purpose

Adhesive cementation of indirect anterior and posterior restorations

#### Patient target group

Patients with permanent teeth

#### Intended users / Special training

- Dentists
- No special training required.

#### Use

For dental use only.

#### Description

Multilink® Automix is a polymer-based dental luting material (intraoral dual-curing) for the adhesive cementation of indirect restorations made of metal and metal-ceramic, all-ceramic and composite materials.

Multilink Automix is exclusively applied together with the self-etching and self-curing Multilink Primer. Multilink Primer is used to condition the dental hard tissue and any existing composite build-up fillings.

The special filler composition gives Multilink Automix a very high radiopacity ( $\geq 300\%$  Al).

Monobond Plus is recommended as a coupling agent to achieve a strong bond to precious and non-precious alloys, as well as to all-ceramics made of zirconium and aluminium oxide and silicate ceramics.

#### Technical specifications

Classification acc. to ISO 4049

- Type 2 / Class 3
- Dual-curing luting composite that is light-activated intraorally
- Optionally, Multilink Automix can be cured with light in the wavelength range of 400–500 nm.

#### Shades

Multilink Automix is available in four shades with different degrees of translucency:

- transparent (high translucency)
- yellow (high translucency)
- opaque (low translucency)
- white (medium translucency)

#### Working time

The working and setting times depend on the ambient temperature. Once Multilink Automix has been dispensed from the automix syringe, the following times apply:

When used purely self-curing	At room temperature 23 °C ± 1 °C (73 °F ± 2 °F)	Intraorally (in combination with Multilink Primer) 37 °C ± 1 °C (99 °F ± 2 °F)
Working time	approx. 3 min	approx. 2 min
Setting time (including working time)	approx. 8 min	approx. 5 min

#### Mixing ratio

Multilink Automix is always dispensed from the automix syringe in the optimum ratio. Multilink Primer A and Multilink Primer B are mixed in a 1:1 ratio (e.g. 1 drop Primer A + 1 drop Primer B).

#### Indications

- Missing tooth structure in anterior and posterior teeth
- Partial edentulism in the anterior and posterior region

#### Areas of application:

- Permanent cementation of glass-ceramic, lithium disilicate glass-ceramic and composite resin restorations (inlays, onlays, partial crowns, crowns, bridges)
- Permanent cementation of opaque restorations made of e.g. oxide ceramics or metal
- Permanent cementation of root canal posts

#### Contraindications

The use of the product is contraindicated if the patient is known to be allergic to any of its ingredients.

#### Limitations of use

If a dry working field cannot be established or the prescribed application technique cannot be applied.

#### Side effects

Systemic side effects are not known to date. In individual cases, allergic reactions to individual components have been reported.

#### Interactions

- Phenolic substances (e.g. eugenol, wintergreen oil) inhibit polymerization. Consequently, the application of products containing these components, e.g. mouth rinses and temporary cements, must be avoided.
- Disinfectants with an oxidative effect (e.g. hydrogen peroxide) may interact with the initiator system, which in turn may impair the curing process. Therefore, do not disinfect the preparation using oxidative agents.
- Alkaline jet media (e.g. Airflow) on dentin can compromise the effect of the Multilink Primer.
- Hemostatic substances may inhibit polymerization and/or lead to discolouration. Therefore, the instructions for use of these substances must be observed.

#### Clinical benefit

- Reconstruction of chewing function
- Restoration of esthetics

#### Composition

##### Multilink Automix

Barium glass, ytterbium trifluoride, Bis-EMA, HEMA, Bis-GMA, Si-Zr mixed oxide, barium aluminium fluorosilicate glass, UDMA, highly dispersed silicon dioxide

The total content of inorganic filler is approx. 40 vol%.

Particle size of inorganic fillers is 0.15 µm–15.5 µm.

##### Multilink Primer A and B

Water, phosphonic acid acrylate, HEMA, sulfinate, p-Toluidine, methacrylate-modified polyacrylic acid

## 2 Usage

For more detailed information also refer to the Instructions for Use of the products used in conjunction with Multilink Automix.

### 2.1 Removal of the temporary restoration and cleaning of the cavity

Remove possible residue of the temporary luting cement from the cavity or the preparation with a polishing brush and an oil- and fluoride-free cleaning paste (e.g. Proxyl<sup>®</sup> fluoride-free prophyl paste). Rinse with water spray. Subsequently, dry with oil- and moisture-free air. Avoid overdrying.

**Note:** Cleaning with alcohol can lead to dehydration of dentin.

### 2.2 Try-in of the restoration and isolation

Seat the restoration using the desired Multilink Try-In Paste and check the shade, fit and occlusion of the restoration.

Care should be taken when checking the occlusion of fragile and brittle ceramic objects, as there is a risk of fracture before they are permanently cemented. If necessary, make adjustments with fine diamonds at medium speed and with slight pressure and adequate water cooling. Polish ground surfaces. Adequate relative or absolute isolation using suitable auxiliaries, such as OptraGate<sup>®</sup> or OptraDam Plus, is required. If the dental hard tissue is contaminated with blood or saliva during the try-in procedure, it needs to be cleaned again as described in section 2.1.

### 2.3 Pre-treatment of the restoration

2.3.1 If a restoration which has been pre-treated in the dental laboratory is tried in, the bonding surface of the restoration needs to be prepared as follows after the try-in:

- Following try-in, thoroughly rinse the restoration with water spray and dry with oil-free air.
- Shake Ivoclean before use and completely coat the bonding surface of the restoration with the product using a microbrush or brush.
- Leave Ivoclean to react for 20 s, then thoroughly rinse with water spray and dry with oil-free air.
- Next, prime the bonding surface of the restoration with a suitable bonding agent (e.g. Monobond<sup>®</sup> Plus).

Make sure to observe the Instructions for Use of the bonding agent used.

2.3.2 Restorations that have not been pre-treated by the dental laboratory should be rinsed with water spray and dried after the try-in procedure.

2.3.2.1 Subsequently, condition the bonding surface of the restoration as follows:

- a) Glass-ceramic restorations (e.g. IPS Empress<sup>®</sup>)
  - Etch the restoration with 5% hydrofluoric acid (e.g. IPS<sup>®</sup> Ceramic etching gel) for 60 s or according to the instructions of the manufacturer of the restorative material being used.
  - Thoroughly rinse the restoration with water spray and dry with oil-free air.
- b) Lithium disilicate glass-ceramic restorations (e.g. IPS e.max<sup>®</sup> Press, IPS e.max CAD)
  - Etch the restoration with 5% hydrofluoric acid (e.g. IPS Ceramic etching gel) for 20 s or according to the instructions of the manufacturer of the restorative material being used.
  - Thoroughly rinse the restoration with water spray and dry with oil-free air.
- c) Restorations made from zirconium oxide (e.g. IPS e.max ZirCAD) or aluminium oxide ceramics
  - Sandblast the bonding surface (sandblasting parameters according to the instructions of the manufacturer of the restorative material being used).
  - If necessary, clean the restoration in an ultrasonic unit for about one minute.
  - Thoroughly rinse the restoration with water spray and dry with oil-free air.
  - **IMPORTANT!** In order to achieve an optimum bond, do not clean zirconium oxide surfaces with phosphoric acid.
- d) Metal or metal-supported restorations
  - Sandblast the inner restoration surfaces (sandblasting parameters according to the instructions of the manufacturer of the restorative material) until a uniformly mat surface is achieved.
  - If necessary, clean the restoration in an ultrasonic unit for about 1 minute.
  - Thoroughly rinse the restoration with water spray and dry it with oil-free air.
  - **IMPORTANT!** In order to create a strong bond, do not clean the metal surfaces with phosphoric acid.

2.3.2.2 Apply Monobond Plus with a brush or a microbrush to the pre-treated surfaces, allow to react for 60 s and then disperse with a strong stream of air.

2.3.3 Alternatively, restorations made of glass-ceramics and lithium disilicate glass-ceramics (e.g. IPS Empress, IPS e.max Press, IPS e.max CAD) may be conditioned with Monobond Etch & Prime<sup>®</sup> as follows:

- Following try-in, thoroughly rinse the restoration with water spray and dry with oil- and moisture-free air.
- Apply Monobond Etch & Prime on the bonding surface using a microbrush and scrub it into the surface for 20 s. Allow to react for another 40 s.
- Thoroughly rinse off Monobond Etch & Prime with water until the green colour has disappeared. If any residue remains in the microporosities, the restoration may be cleaned with water in an ultrasonic unit for up to 5 min.
- Dry the restoration with a strong stream of oil- and moisture-free air for approx. 10 s.

2.3.4 Composite and fibre-reinforced composite restorations should always be conditioned according to the instructions of the manufacturer of the restorative material being used.

With Tetric<sup>®</sup> CAD restorations, please observe the instructions stated in the Tetric CAD Instructions for Use.

### 2.4 Mixing of Multilink Primer A and Multilink Primer B

Mix the two Multilink Primer liquids A and B at a 1:1 mixing ratio (e.g. 1 drop of Primer A and 1 drop of Primer B). The mixed Primer A/B is purely self-curing and does not need to be protected against light. However, it must be applied within 10 min.

### 2.5 Application of the mixed Multilink Primer A/B to the enamel, the dentin and the build-up

Apply the mixed Multilink Primer A/B onto the entire bonding surface using a microbrush, starting from the enamel surface. Subsequently, scrub it in for another 30 s. Disperse excess with blown air until the mobile liquid film is no longer visible.

As the Primer is purely self-curing, no light-curing is necessary! **Note:** Previous conditioning of enamel with phosphoric acid is required for Maryland bridges.

### 2.6 Application of Multilink Automix on the restoration

- For each application place a new mixing tip on the syringe. Dispense Multilink Automix from the automix syringe and apply the desired amount onto the restoration. As the luting material will cure in the used mixing tip, it serves as a seal for the remaining contents of the syringe until needed again (replace with a new mixing tip before the next use).
- Multilink Automix should be processed quickly after having been dispensed from the automix syringe and the restoration seated into place.

Apply Multilink Automix directly to the inner surface of the restoration.

**Note:** The direct application of Multilink Automix to the preparation or cavity pretreated with Multilink Primer is contraindicated, as this would significantly accelerate the curing process and could adversely affect the accuracy of fit of the restoration.

### 2.7 Placement of the restoration and removal of excess cement

a) Purely self-curing

Seat the restoration in place and secure it. Remove excess material immediately with a micro-brush/brush/foam pellet/dental floss or scaler. Make sure to remove excess material in time, especially in areas that are difficult to reach (proximal areas, gingival margins, pontics). Due to the reaction between Multilink Automix and Multilink Primer A/B, a high bond strength and a high degree of cure are achieved within a few minutes after the restoration has been seated.

b) Self-curing with additional light-curing

(quarter technique, indicated for cases with up to 2 bridge abutments = 3- to 4-unit bridges)  
Light-curing of all excess cement and subsequent removal: Seat the restoration in place and secure it. Light-cure excess material in quarter segments (mesio-oral, disto-oral, mesio-buccal, disto-buccal) using a curing light at a distance of max. 10 mm.

For this procedure, the following parameters must be observed:

Light intensity	Exposure time per quarter segment
approx. 650 mW/cm <sup>2</sup>	3 s
approx. 1,000 mW/cm <sup>2</sup>	1-2 s e.g. Bluephase G4, PreCure mode

Following this procedure, excess material can be easily removed with a scaler. Make sure to remove excess material in time, especially in areas that are difficult to reach (proximal areas, gingival margins, pontics). Subsequently, light-cure all the margins again for 20 s (approx. 1,000 mW/cm<sup>2</sup>, e.g. using Bluephase G4 in High Power mode).

- c) Self-curing with additional light-curing (quarter technique, indicated for cases with up to 6 bridge abutments = circular bridges)  
Light-curing of all excess cement and subsequent removal. Seat the restoration in place and secure it. Successively light-cure excess material on all the abutments in quarter segments (mesio-oral, disto-oral, mesio-buccal, disto-buccal) using a curing light at a distance of max. 10 mm. For this procedure, the following parameters must be observed:

Light intensity	Exposure time per quarter segment
approx. 650 mW/cm <sup>2</sup>	1 s

After this procedure, excess cement is easy to remove with a scaler. Make sure to remove excess material in time, especially in areas that are difficult to reach (proximal areas, gingival margins, pontics). Subsequently, light-cure all the margins again for 20 s (approx. 1,000 mW/cm<sup>2</sup>, e.g. using Bluephase G4 in High Power mode).

#### Note

As with all composites, Multilink Automix is subject to oxygen inhibition. This means that the surface layer (approx. 100 µm) does not polymerize during curing, as it comes in contact with atmospheric oxygen. To prevent this, it is advisable to cover restoration margins with glycerine gel/air block (e.g. Liquid Strip) immediately after excess removal and rinse off the protective film after complete polymerization.

### 2.8 Finishing of the completed restoration

- Check occlusion and functional movements and make adjustments if necessary.
- Finish the cement lines with finishing diamonds if necessary.
- Smooth out the cement lines using finishing and polishing strips and polish them with suitable polishing instruments (e.g. OptraGloss®).
- If necessary, finish the restoration margins with suitable polishers (e.g. OptraGloss).

### Special instructions for the cementation of endodontic posts

- For the cementation of endodontic posts, carefully clean the root canal to remove any residue of root canal filling materials. (Residue of eugenol-based sealers may inhibit the polymerization of the luting composite.)
- Ideally, the mixed Multilink Primer A/B should be applied to the root canal and the occlusal surface of the preparation for about 15 s each using the thin (violet) microbrushes.
- Remove excess material from the root canal using paper points.
- Coat the endodontic post, which you have prepared according to the instructions of the endodontic post manufacturer, with the dispensed Multilink Automix cement.  
**IMPORTANT!** Do not apply Multilink Automix cement with the lentulo into the root canal that has been wetted with Primer A/B. Otherwise, premature curing might occur, which could possibly prevent the post from being properly positioned.
- Place the endodontic post. Excess cement will be displaced when the post is inserted. Uniformly distribute excess cement across the occlusal preparation surface. Then, light-cure Multilink Automix for 20 s.
- Core buildup material can then be directly applied onto Multilink Automix, which serves as bonding agent, and cured according to the instructions of the manufacturer (light-curing).

### 3 Safety notes

- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: [www.ivoclar.com](http://www.ivoclar.com), and your responsible competent authority.
- The current Instructions for Use are available in the download section of the Ivoclar Vivadent AG website ([www.ivoclar.com](http://www.ivoclar.com)).
- Explanation of the symbols: [www.ivoclar.com/eIFU](http://www.ivoclar.com/eIFU)
- The Summary of Safety and Clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at <https://ec.europa.eu/tools/eudamed>.
- Basic UDI-DI: 76152082ACEME001DV

### Warnings

- Observe the Safety Data Sheet (SDS) (available in the download section of the Ivoclar Vivadent AG website [www.ivoclar.com](http://www.ivoclar.com)).
- Contact of the uncured material with the skin, mucous membrane or eyes must be avoided.
- Uncured material may cause slight irritation and may lead to a sensitization against methacrylates.
- Commercial medical gloves do not provide protection against the sensitizing effects of methacrylates.

### Disposal information

Remaining stocks must be disposed of according to the corresponding national legal requirements.

### Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks.

The following clinical residual risks are known:

- Failure of the adhesive bond
- Postoperative sensitivity

### 4 Shelf life and storage

- Storage temperature: 2 - 28 °C (36 - 82 °F)
- Store automix syringe with mixing tip attached after use.
- Multilink Primer bottles must be tightly closed with the correct cap after use.
- Do not use the product after the indicated date of expiration.
- Expiry date: see note on bottle, automix syringe and package.

Before use, visually inspect the packaging and the product for damage. In case of any doubts, please contact Ivoclar Vivadent AG or your local dealer.

### 5 Additional information

Keep material out of the reach of children!  
Not all products are available in all countries.

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from misuse or failure to observe the Instructions. The user is responsible for testing the product for its suitability and use for any purpose not explicitly stated in the Instructions.

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d) Metal or metal-supported restorations

- Sandblast the inner restoration surfaces (sandblasting parameters according to the instructions of the manufacturer of the restorative material) until a uniformly mat surface is achieved.
- If necessary, clean the restoration in an ultrasonic unit for about 1 minute.
- Thoroughly rinse the restoration with water spray and dry it with oil-free air.
- **IMPORTANT!** In order to create a strong bond, do not clean the metal surfaces with phosphoric acid.

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