Calibra® Abutment

Resin Cement

INSTRUCTIONS FOR USE **ENGLISH** 2

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Cemento composito

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Cemento de resina

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Cimento resinoso

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Calibra® Abutment Resin Cement

INSTRUCTIONS FOR USE – ENGLISH

Caution: This is a medical device. For dental use only by dental professionals.
USA: Rx only.

1. PRODUCT DESCRIPTION

Calibra® Abutment Resin Cement is a two-component, self-cure, self-adhesive cement.

Intended Purpose

Calibra® Abutment Resin Cement is a dental cement which is intended for the extraoral cementation assembly of prosthetic components.

User Profile

Calibra® Abutment Resin Cement is a prescription device, used by or under the supervision of dental professionals licensed by law, with high technical knowledge, specialized education and training. Calibra® Abutment Resin Cement is to be used only in a laboratory, clinical or hospital environment, following good dental laboratory practice, by qualified dental professionals such as laboratory technicians, general practitioners as well as specialists and Dental Assistants.

Clinical Benefits to be expected

Calibra® Abutment Resin Cement has the ability to bond to titanium and is chemically compatible with silanated ceramic materials, making it suitable for extraoral assembly of prosthetic components. This self-cure only, opaque shade has been optimized for extraoral assembly when blocking of metal gray-ing or show-through is desired.

Cured Calibra® Abutment Resin Cement is essentially hydrophobic, minimizing post-cure water sorp-tion, solubility and hygroscopic expansion.

Intended Patient Population and Medical Conditions

Calibra® Abutment Resin Cement is intended for use in fabricating custom implant hybrid abutments or screw-retained hybrid prostheses for the following patient population/target groups and applications: Age-appropriate patients with dental implant(s) undergoing implant-supported dental prosthetic restor-ative procedures.

1.1 Indications

Calibra® Abutment Resin Cement is indicated for the cementation of: Extraoral assembly of prosthetic components, e.g., structures made of ceramic or zirconia onto custom or pre-fabricated titanium/tita-nium alloy or zirconia bases/implant inserts/abutments.

1.2 Contraindications

Calibra® Abutment Resin Cement is contraindicated for use with patients who have a known hypersen-sitivity to methacrylate resins.

1.3 Delivery forms (Some delivery forms may not be available in all countries)

Calibra® Abutment Resin Cement is available in:

- a convenient dual-barreled syringe
- 1 self cure, extraoral use only, opaque shade, for blocking undesirable substrate show-through

1.4 Composition

Silanated barium boron-fluoro-aluminosilicate-glass < 70%, Urethane dimethacrylate < 10%, Trimeth-acrylate resin < 10%, Phosphoric acid modified acrylate resin < 10%, Triethyleneglycol dimethacrylate < 10%, Titanium dioxide < 10%, Hydroxyethylmethacrylate < 5%, Highly dispersed silicon dioxide < 5%, Hydroxypropyl methacrylate < 5%, Initiator (Hydroperoxide) < 1%, Acrylic acid < 1%, Co-initiator (Thio-urea) < 1%, UV stabilizer < 1%, Butylated hydroxy toluene < 1%, Iron oxide pigments < 1%.

1.5 Compatible Substrates

Calibra® Abutment Resin Cement is compatible with all high-strength ceramic and zirconia substrates for cementation to titanium/titanium alloy or zirconia bases.

2. GENERAL SAFETY NOTES

Be aware of the following general safety notes and the special safety notes in other chapter of these Instructions for Use.

Safety alert symbol.

This is the safety alert symbol. It is used to alert you to potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury.

2.1 Warnings

1. Calibra® Abutment Resin Cement is acidic in nature and contains polymerizable acrylate and meth-acrylate monomers which may be irritating to skin, eyes and oral mucosa, and may cause allergic contact dermatitis in susceptible persons.
 - **Avoid eye contact** to prevent irritation and possible corneal damage. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.
 - **Avoid skin contact** to prevent irritation and possible allergic response. In case of contact, red-dish rashes may be seen on the skin. If contact with skin occurs, immediately remove material with cotton and wash thoroughly with water and soap. In case of skin sensitization or rash, dis-continue use and seek medical attention.
 - **Avoid contact with oral soft tissues/mucosa** to prevent inflammation. If accidental contact occurs, immediately remove material from the tissues. Flush mucosa with plenty of water after the restoration is completed and expectorate/evacuate the water. If sensitization of mucosa per-sists, seek medical attention.
2. Calibra® Abutment Resin Cement should not be used with patients who have a known hypersensi-tivity to any of the components.
3. Do not use adhesive in conjunction with Calibra® Abutment Resin Cement. Cement should be placed directly on prepared substrates (see step-by-step). Adhesives may degrade during pre-delivery cleaning/disinfection/sterilization processes, leading to debonding.

2.2 Precautions

1. This product is intended to be used only as specifically outlined in the Instructions for Use. Any use of this product inconsistent with the Instructions for Use is at the discretion and sole responsibility of the practitioner.
2. Insufficient data exist to support the use of Calibra® Abutment Resin Cement for intraoral cemen-tation of restorations.

3. Wear suitable protective eyewear, clothing and gloves. Protective eyewear is recommended for patients.
4. Devices marked "single use" on the labeling are intended for single use only. Discard after use. Do not reuse in other patients in order to prevent cross-contamination.
5. The syringe cannot be reprocessed. To prevent the syringe and from exposure to spatter or spray of body fluids or contaminated hands it is mandatory that the syringe is handled with clean/disin-fected gloves.
 - As additional precautionary measure, syringe may be protected from gross debris but not from all contamination by applying a protective barrier.
6. Do not insert or use prostheses which have not been properly processed. Before delivery, assem-bly must be properly cleaned/disinfected/autoclaved according to manufacturer's Instructions to reduce the risk of cross contamination.
7. Remove excess cement while cement is in the gel stage. Do not allow excess to harden, making cleanup difficult. Remove excess cement before restoration delivery. Residual excess cement may lead to inflammation.
8. Cement at the margins may appear set before cement under the restoration is set. Do not move, torque or disturb restoration until final set of the cement (6 minutes from the start of mix).
9. Exercise caution when placing multiple units not to exceed available working time. Load and ensure complete seating of each unit prior to loading successive units. The use of a fresh mix of cement is recommended if available working time is approached.
10. Containers should be tightly closed immediately after use. Syringe should be tightly closed by replacing the original cap immediately after use.
11. Calibra® Abutment Resin Cement should extrude easily. **DO NOT USE EXCESSIVE FORCE.** Excessive pressure may result in unanticipated extrusion of the material or cause syringe rupture.
12. Material should be mixed in equal volumes for optimum performance. Dispense and discard a small amount of material (bleed syringe ports) prior to attaching mixing tip. Uneven, non-homoge-neous dispensing may lead to compromised physical properties and performance.
13. Numerous factors can affect restoration longevity. Marginal leakage can result in restoration loss. Check restoration marginal integrity at every periodic examination.

Interactions

- Eugenol containing materials should not be used in conjunction with this product because they may interfere with hardening and cause softening of the polymeric components of the material.
- Apply silane (available separately) only to ceramic restorations designed to be silanated. Do not use pre-treatment agents, e.g., primers or bonding agents on any surface when cementing with Calibra® Abutment Resin Cement. Use of any such pre-treatments may compromise the bonding of the restoration subjected to cleaning/disinfecting/autoclaving. See step-by-step instructions below for recommended restoration treatment.

2.3 Adverse reactions/Undesirable side effects

1. Product may irritate the eyes and skin. Eye contact: irritation and possible corneal damage. Skin contact: irritation or possible allergic response. Reddish rashes may be seen on the skin. Mucous Membranes: inflammation (see Warnings).
2. Residual hardened excess cement may lead to soft tissue injury or irritation (see Precautions, Step-by-Step instructions).

2.4 Storage

Inadequate storage conditions may shorten the shelf life and may lead to malfunction of the product.

- Keep out of direct sunlight
- Store in a well-ventilated place
- Store at temperatures between 2 °C-24 °C (35 °F-75 °F).
- Use the product at room temperature. Allow material to reach room temperature prior to use.
- Protect from moisture.
- Do not freeze.
- Do not use after expiration date.
- Do not use if the lot number and/or expiration date are missing or illegible.

3. STEP-BY-STEP INSTRUCTIONS

3.1 Extraoral Cementation of Mesostructure abutment/Full-contour Prosthesis to TiBase abutment

3.1.1 Preparation of TiBase

- The diameter of the TiBase must not be reduced, modified or adjusted e.g., by grinding.
- Shortening the TiBase is not permitted.
- The underside contact/connection surfaces of the TiBase to the implant should not be sand blasted, polished or otherwise modified or processed. Use of an appropriate analogue to protect connection surfaces is strongly recommended. Wax or impression material may also be placed on exposed emergence profile surfaces for protection.

1. Attach the TiBase to the lab implant analogue or polishing tool using screw.
2. Cover the hex head of the abutment screw and fill the screw channel with easily removable mate-rial, e.g., soft wax, impression material, cotton or plumber's tape.
3. Ensure that the ceramic/zirconia Mesostructure/Prosthesis can be fully inserted onto the TiBase. The shoulder of the Mesostructure/Prosthesis must contact the TiBase flange without perceptible gap. Place small alignment marks on the Mesostructure/Prosthesis and the analogue using a felt-tip marker to assist aligning the correct Bu-Li position.
4. Sandblast only the surfaces of the TiBase intended for cementation (50 µm aluminum oxide, max. 2.0 bar). Do not sandblast the polished tissue surface of the flange.
5. Clean the TiBase with alcohol, ultrasonic or steam. If necessary, reseal the abutment screw and screw channel. Dry with clean, oil-free air.
6. Once cleaned, do not touch or contaminate the bonding surface. Resin cement may be placed directly on the sandblasted bonding surface. Do not place adhesive on the bonding surface of the TiBase restoration.

3.1.2 Preparation of Mesostructure/Prosthesis

A. Glass Ceramic

1. Treat restoration according to manufacturer's Instructions for Use, e.g., clean the restoration with an ultrasonic cleaner, steam cleaner or with alcohol.
2. Apply 5% to 9% hydrofluoric acid etching gel (available separately, see complete Instructions for Use) to the bonding surface of the restoration only. Do not allow etchant into the screw channel.

Etching time: 30 seconds.

Adicid material – To reduce the risk of injury.

- Observe manufacturer's Warnings, Instructions for Use, Material Safety Data Sheets.
- Do not allow acid to come into contact with tissue or the eyes.
- Remove the hydrofluoric acid according to the manufacturer's Instructions for Use.

3. It is recommended to immediately silanate the etched surfaces. Apply silane only to surfaces that are required for adhesive cementation. Allow to react for 60 seconds. If the silane layer is then no longer liquid, apply silane again. Blow off with a strong stream of air (Recommended material: Calibra® Silane Coupling Agent – sold separately, see full instructions for use).
4. Once silanated, resin cement may be placed directly on the silanated surface. Do not place adhe-sive on the bonding surface of silanated restorations.

B. Zirconia

1. Treat restoration according to manufacturer's Instructions for Use.
2. Internal surface microetching (sandblasting with 50 µm alumina) of bonding surfaces of the restoration is recommended. Do not sandblast the screw channel.
3. Clean the restoration with an ultrasonic cleaner, steam cleaner or with alcohol.

Debonding – To reduce risk of inadequate retention due to pre-delivery processing.

- Do not apply primers or adhesives to zirconia bonding surfaces.
- Apply cement directly to sandblasted and cleaned zirconia surfaces.

3.1.3 Cementation of Mesostructure/Prosthesis to TiBase

1. Prepare and use syringe – Follow dispensing instructions outlined below.

Material Extrusion – To reduce risk of injury due to excessive force.

- Apply slow and steady pressure on the syringe.
- Do not use excessive force – syringe rupture may result.

2. Dual Barreled syringe dispensed – Extraoral Cementation.

- Remove syringe cap. Dispense and discard a small amount of material from the dual-barreled syringe. Be sure material is flowing freely from both ports. Holding syringe vertically, carefully wipe away excess so base and catalyst do not cross contaminate and cause obstruction of the ports. Save syringe cap for replacement following use.

Mixing Discrepancies – To reduce risk of compromised performance.

- Store with original cap only.
- Bleed both ports prior to attaching mixing tip.
- Wipe excess from port openings after bleeding and before replacing storage cap. Do not allow base/catalyst cross contamination.

- Install a mixing tip on the cartridge by lining up the v-shaped notch on the outside of the mixing tip with the v-shape notch on the syringe flange. Turn colored mixing tip cap 90 degrees in a clockwise direction to lock in place on syringe.
- Gently depress syringe plungers to begin the flow of material. DO NOT USE EXCESSIVE FORCE. If force is encountered, remove syringe from operating field, remove and discard mixing tip. Check for obstruction and confirm material is flowing from both syringe barrels. Wipe barrels and install new mixing tip as outlined above. Dispense a small amount through the mixing tip onto a mixing pad and discard.

Technique Tip: The opaque shade has been optimized for extraoral assembly and blocking of metal show-through. Material is designed to be self-cure only.

3. Apply the cement to the bonding surfaces of the TiBase and the inside of the restoration. Minimize the cement allowed to enter the restoration screw channel.
 - At room temperature, Calibra® Abutment Resin Cement offers a minimum work time of 2 minutes.
 - Immediately seat the restoration with gradual pressure. Verify complete seating. A gentle rocking or vibratory motion may be helpful to insure optimal seating.
4. Align the pen marks and push the Mesostructure/Prosthesis over the TiBase bonding chimney up to the end stop. Note the locking of the rotation and position lock.
 - Ensure margin between TiBase and Mesostructure/Prosthesis is closed and sealed.

To reduce risk of incomplete cementation.

- Assure restoration is fully seated within 2 minutes.
- Placement of multiple units may require multiple mixes of cement.

5. Protect restoration from contamination and movement until the final set of the cement (6 minutes from start of mix).

Debonding – To reduce risk of inadequate retention due to insufficient stabilization.

- Stabilize restoration during cleanup.
- Protect restoration from movement during setting.
- Allow cement to self-cure without disturbing for 6 minutes.

3.1.4 Cleaning Marginal Excess

Self-cure cleanup

1. Remove coarse cement excess immediately by wiping with a soft tissue, cotton or similar material. Ensure there is no cement excess in the screw channel. Remove excess while cement is still fluid using a flocked applicator tip or brush.
2. Do not allow any movement between the Mesostructure/Prosthesis and the TiBase for three minutes. A clamp may be useful for stabilizing the Mesostructure/Prosthesis on the TiBase.
3. The excess cement will reach the "gelled" state after approximately 3 minutes, allowing easy removal. Excess cement will remain in the "gelled" state for approximately 1 minute. Complete excess cement removal using an instrument such as a rubber tip, a scaler or an explorer.

Note: Cement within the restoration has not yet set. Do not move, torque, or disturb the restoration during cleanup.

4. Apply water-soluble air barrier to the cement margin.
5. Protect assembly from contamination and movement until the final set of the cement (6 minutes from start of mix).

3.1.5 Finishing and Polishing

1. Rinse air barrier material under running water.
2. Remove wax or other sealing material and any residual cement from the screw channel.
3. Remove any cement excess at the margin, finish and polish margin using rubber rotary tools. Proceed to Hygiene for pre-use prosthesis preparation.

Residual excess cement – To reduce the risk of soft tissue irritation.

- Ensure all excess cement is removed.
- Ensure all surfaces are polished and free of contamination.

Calibra® Abutment Resin Cement is radio-opaque, with a 1 mm radio-opacity equivalent to 2 mm radio-opacity of aluminum. Aluminum has a radio-opacity equivalent to that of dentin. Thus 1 mm of material having a radio-opacity equivalent to 1 mm of aluminum has a radio-opacity equivalent to that of dentin.

4. HYGIENE, PROCESSING AND DISPOSAL

4.1 Assembled Prosthesis

Cross-contamination – To reduce risk of infection.

- Prostheses must be polished before cleaning, disinfection and/or sterilization.
- Prostheses must be processed according to manufacturer's validated processing Instructions prior to delivery.

Calibra® Abutment Resin Cement has been shown to tolerate two cycles of the processes listed below without significant effect on physical properties.

See prosthesis manufacturer's Instructions for Use for detailed instructions and validated cleaning/disinfection/sterilizing process parameters as applicable.

- Wiping with or immersion in 70% isopropyl alcohol for up to 10 minutes.
- Immersion in boiling water for up to 5 minutes. (Recommended only for non-US markets)
- Automated Washer-Disinfector meeting the requirements of ISO 15883, run program with A0 value ≥ 3000 (e.g., 5 min at $\geq 90^\circ\text{C}$ (194 $^\circ\text{F}$)) (e.g., Miele Vario TD) using appropriate detergents (e.g., neodisher® MediClean and neodisher® Z, both from Dr. Weigert, Germany or equivalent).
- Dynamic-Air-Removal (Pre-Vacuum) Steam autoclaved at:
 - 132 $^\circ\text{C}$ (270 $^\circ\text{F}$) for 4 minutes; Dry Time = 20 minutes
 - 134 $^\circ\text{C}$ (270 $^\circ\text{F}$) for 3 minutes; Dry Time = 16 minutes (Recommended only for non-US markets)
 - 135 $^\circ\text{C}$ (275 $^\circ\text{F}$) for 3 minutes; Dry Time = 16 minutes

To reduce risk of inadequate retention due to wrong processing method – Do not use ultrasonic cleaner.

- Cleaning/Disinfection of cemented assembly using an ultrasonic unit is not validated and should not be used.
- Physical properties may be adversely affected.
- Manual wiping or static immersion only for manual cleaning and manual disinfection.

- Follow Prosthesis manufacturer's Instructions for Use for pre-use disinfection and/or sterilization.
- Locally applicable legal regulations and the hygiene standards applicable for a dental practice must be observed.
- Use only the validated processing procedures specified by the prosthesis manufacturer.
- The responsibility for the microbial processing of the prosthesis lies with the user.
- The equipment and devices must be properly maintained and serviced at regular intervals.
- The fabricator (laboratory or in-office dental technician) of the assembled prosthesis must inform the dentist of the need to process the device before using or inserting it in the patient's mouth.

4.2 Cement Syringes

Cross-contamination – To reduce the risk of infection.

- Do not reuse single use products. Dispose of in accordance with local regulations.
- Syringe cannot be reprocessed. Dispose of contaminated syringe in accordance with local regulations.
- When not in use, all supplies and instruments should be in covered storage, such as drawers and cabinets, and away from potential contamination.
- During treatment clinicians with patient contact should not handle the syringe.
- To prevent the syringes from exposure to spatter or spray of body fluids or contaminated hands it is mandatory that the syringes are handled with clean/disinfected gloves, outside of the patient treatment operatory.
- Assembling components with clean/disinfected gloves in a separate room, bringing only clean, disinfected and sterilized assembly into the operatory is strongly recommended.
- Dispensing materials with clean/disinfected gloves in a separate room, bringing only what will be used into the operatory is strongly recommended.
- Do not reuse syringes if contaminated.

For dual-barreled syringe, remove used mixing tip and discard appropriately. Replace original syringe cap prior to storing.

To prevent syringes from exposure to spatter or spray of body fluids or contaminated hands, or oral tissues, use of a protective barrier is recommended. The use of protective barriers is an additional precautionary measure against gross debris but not against all contamination. The barrier is designed for single use and must be disposed of after each use in accordance with local regulations.

Incidental contact of the syringe with water, soap or a water-based hospital-level disinfection solution will not damage syringe body. Do not allow any solution contact with contained material. Discard material that has been in contact with any fluid or non-sterile instrument.

NOTE: Vigorous wiping can destroy the label. Wipe syringe gently.

Repeated liquid contact may damage label. Dry syringes with a lint-free single-use cloth.

4.3 Disposal

- Based on the information in the SDS, the waste generated by the disposal of the product does not present characteristics of physical danger, for the person or the environment; It can therefore be managed, according to local regulations, as a non-hazardous waste.
- Contaminated containers and accessories that have been in contact with the patient must be cleaned and disinfected as outlined in Hygiene above prior to disposal or disposed of as hazardous waste at risk of biological contamination.
- Always wear gloves when handling packaging and accessories.

5. LOT NUMBER, EXPIRATION DATE AND CORRESPONDENCE

1. Do not use after expiration date. Do not use if the lot number and/or expiration date are missing or illegible. ISO standard uses: "YYYY-MM-DD"
2. The following numbers should be quoted in all correspondences:
 - Reorder Number (REF)
 - Lot number
 - Expiration Date
3. Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.
4. A summary of safety and clinical performance (SSCP) for this product can be found (upon activation) at <https://ec.europa.eu/tools/eudamed> by searching using the Basic UDI-DI number ++D002DENTALCEMENTSC6 and at <https://www.dentsplysirona.com/ifu> using the reference number (REF).

150-DV0

© 2024 Dentsply Sirona Inc. **507702WEB** (R 9/30/24)

Manufactured by
Dentsply Caulk
38 West Clarke Avenue
Milford, DE 19963, USA
Made in USA
1-302-422-4511
www.dentsplysirona.com

CE/REP
Dentsply DeTrey GmbH
De-Trey-Strasse 1
78467 Konstanz, Germany
49-7531-583-0
www.dentsplysirona.com

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Silaneret barium-bortrifluorid-aluminosilicaglas < 70%, Urethandimethacrylat < 10%, Trimethacrylatplast < 10%, Fosforsyremodificeret acrylatplast < 10%, Triethylenglycoldimethacrylat < 10%, Titaniumdioxid < 10%, Hydroxyethylmethacrylat < 5%, Højt dispergeret siliciumdioxid < 5%, Hydroxypropylmethacrylat < 5%, Initiator (hydroperoxid) < 1%, Akrylsyre < 1%, Co-initiator (thiourea) < 1%, UV-stabilisator < 1%, Butyleret hydroxytoluen < 1%, Jernoxidpigmenter < 1%.

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© 2024 Dentsply Sirona Inc. 507702WEB (R 9/30/24)