INSTRUCTIONS FOR USE**PARKELL PRODUCT INFORMATION:****Blu-Mousse® &
Green-Mousse®****VPS Bite Registration / Impression Materials**

A00775revE0318

Product Name Setting Time (mouth temperature) Final Shore A Surface Hardness	2 x 50ml Cartridges "Split" (Total vol. 100ml) (MixPac DS-50 1:1/2:1 gun)	2 x 50ml Cartridges "Standard Non-Split" (Total vol. 100ml) (MixPac Original 1:1 gun)	2 x 80ml Tubes (Total vol. 160ml) (Hand mixed)
Blu-Mousse® Super-Fast VPS (30 second set / 85 Durometer)	S457S	S457	n/a
Blu-Mousse® Super-Fast Peppermint VPS (30 second set / 85 Durometer)	S448S	n/a	n/a
Blu-Mousse® 60 VPS (60 second set / 85 Durometer)	S440S	n/a	n/a
Blu-Mousse® Classic VPS (120 second set / 85 Durometer)	S438S	S438	S436
Blu-Mousse® Classic Peppermint VPS (120 second set / 85 Durometer)	S437S	n/a	n/a
Green-Mousse® VPS (120 second set / 57 Durometer)	S455S	n/a	n/a

PRODUCT DESCRIPTION AND INDICATIONS FOR USE

Parkell's **Blu-Mousse®** and **Green-Mousse®** are addition-reaction silicones, also known as A-Silicones, polyvinyl siloxanes (PVS), or vinyl polysiloxanes (VPS). Mousse products display a light, fluffy consistency after mixing with virtually no resistance to closure, as well as a rigid, hard final set. They are available in setting times of 30 seconds, 60 seconds or 120 seconds. Mousse is indicated for use in occlusal bite registrations, transfer impressions, instant bite rims, and double-arch, one-stage or two-stage impressions (as the tray material). They are capable of producing accurate impressions with excellent tear strength. Mousse materials set to a dull surface that is scannable without powder, using CAD-CAM digital impression scanners.

Parkell's Mousse materials are provided as listed in the table above. Each kit contains the appropriate number and type of MixPac mixing tips or tube squeezers for that product. Accessory products such as VPS Adhesive (Stock No. HS446), the DS-50 mixing gun (Stock No. S343), replacement pink "split" mixing tips (Stock No. S445) or replacement older-style "standard" mixing tips (Stock No. S439K) are available separately from Parkell or many dental dealers.

STORAGE, GENERAL USE INSTRUCTIONS AND PRECAUTIONS FOR MOUSSE IMPRESSION MATERIALS

- Avoid excessive heat, light, moisture or freezing.
- Store and use Mousse products at room temperature (18°-24° C / 64°-75° F). If extended time is required, material may be chilled. Cold material may be difficult to express, and take longer to cure. Excess heat will shorten shelf life.
- After initial use, store in a cool, dark place with used mixing tip tightly screwed on. Do not replace mixing tip until ready for next use, and do not re-use the original cap.
- Excessive or rapid gun force or mixing tip blockage may cause cartridge rupture. To avoid injury or property damage, slow and steady pressure should be applied to the dispensing gun.
- Avoid contact with latex or sulfur-containing "rubber" products, (latex rubber gloves or dams), which may inhibit the set of VPS materials. Nitrile or vinyl gloves are acceptable
- Avoid hemostatic agents with sulfates (e.g. Aluminum Sulfate, Ferric Subsulfate). If field is contaminated with these products, scrub with 5% hydrogen peroxide on gauze or cotton for 15 seconds, and rinse with copious water before impressioning. Chloride hemostatic agents (e.g. Aluminum Chloride, Ferric Chloride) are an acceptable substitute.
- Vinyl polysiloxanes like Mousse are not compatible with polysulfide or polyether materials.
- Bis-acryl resins, core composites and certain bonding agents may leave uncured residue on the surface of restorations or impressions, which may inhibit the setting of VPS products. Remove uncured resin from any surface before impressioning.
- BLU-MOUSSE IMPRESSIONS THAT WILL BE USED FOR WASH IMPRESSIONS LATER SHOULD NOT BE USED AS A MATRIX FOR BIS-ACRYL TEMPORARIES. THE UNCURED BIS-ACRYL SLICK LAYER INSIDE THE BLU-MOUSSE PRIMARY IMPRESSION WILL CONTAMINATE THE SURFACE AND INTERFERE WITH THE BOND TO THE SECONDARY WASH IMPRESSION.**
- Set Blu-Mousse may be soaked or spray-disinfected using readily-available dental surface disinfectants as per the manufacturer's instructions. Consult www.CDC.gov for the most recent version of the "Guidelines for Infection Control in Dental Health-Care Settings".
- Allergic reactions may occur in susceptible individuals. Flush affected tissues with copious water and consult a physician.
- Blu-Mousse sets almost as hard as plaster. It should not be used wherever undercuts may cause locking of the impression, such as full-arch tooth impressions, periodontally-involved teeth, or impressions involving pontics. Flexible Green-Mousse may be preferable.

(Continued)

13. Avoid use of rim-locking metal trays with Blu-Mousse, as these trays will be extremely difficult to remove if the Blu-Mousse locks into dental undercuts. Block out or lubricate undercuts before impressioning, using wax, petroleum jelly or other materials.
14. For best results when using Mousse with other materials, match the setting times to be the same or similar. This will eliminate the error of removing the impression from the mouth too soon, when the outside is set, but the inside is not.
15. **USE ONLY GENUINE PARKELL/MixPac MIXING TIPS, NOT COUNTERFEIT MIXING TIPS! PARKELL IS NOT RESPONSIBLE FOR DAMAGES CAUSED BY THE USE OF UNAUTHORIZED MIXING TIPS.**

USING IMPRESSION GUNS WITH MOUSSE MATERIALS ("THE DOUBLE-BLEED TECHNIQUE")

1. Insert Mousse cartridge into impression gun. Twist off and discard the sealing cap. Bleed about 1/4" of material from both orifices, ensuring that it is extruding evenly from both holes. Wipe the ends clean, avoiding cross-contamination. Attach appropriate mixing tip and tighten 1/4 turn.
2. When ready for use, bleed a small amount of Mousse thru the mixer onto a pad and immediately dispense the material for the desired application.
3. Wait the appropriate time and remove the set impression with a quick, firm motion. Use the chart above to determine the correct setting time.
4. Do not reuse the shipping cap. Leave the used mixing tip on the cartridge as the new sealing cap.
5. The cartridge, used mixing tip and dispensing gun should not be directly sprayed with or soaked in disinfectant. They can be cleaned of debris with a damp paper towel, disinfected by wiping with a paper towel that has been saturated with an EPA-registered low-level (HIV/HBV claim) to intermediate-level (tuberculoocidal claim) hospital disinfectant, wiped of chemical residue with a water-damp paper towel, and dried. Utilize the surface disinfecting protocol of the disinfectant manufacturer. Consult www.CDC.gov for the most recent version of the "Guidelines for Infection Control in Dental Health-Care Settings".

6. Store the Mousse cartridge horizontally until the next use.
7. Parkell's Quick-Pour formula allows Mousse impressions to be poured in any dental stone in 15 minutes, up to 30 days later.

USING BLU-MOUSSE CLASSIC IN TUBES

1. Dispense equal lengths of base and catalyst (1:1) onto a coated mixing pad. Altering the ratio will impair final set. Recap tubes immediately.
2. Mix quickly by hand with clean spatula until color is uniform and streak-free, with no air bubbles. Load the material into the appropriate dispensing syringe and dispense it immediately.

BITE REGISTRATION PROCEDURE WITH MOUSSE MATERIALS

1. Mousse material may be dispensed directly from a mixing tip or syringe onto the biting surfaces of teeth, or into a conventional disposable plastic tray, or a side-less tray, or a "3 in 1" tray.
2. All bites should be protected with padding during shipping to the lab to preserve their accuracy.

WARRANTY AND TERMS OF USE:

Parkell will replace defective material. **This warranty is in lieu of all warranties of merchantability, fitness for purpose or other warranties, express or implied.** Parkell does not accept liability for any loss or damage, direct, consequential or otherwise, arising out of the use of or the inability to use the product herein described. Before using, the user shall determine the suitability of the product for its intended use and the user assumes all risk and liability whatsoever in connection herewith.

For full Warranty and Terms of Use information, please see www.parkell.com.

Safety Data Sheets (SDS) are available at www.parkell.com.

Parkell's Quality System is certified to ISO 13485.

European Authorized Representative (Not a dealer/distributor):

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Directa AB - P.O. Box 723, Finvids väg 8