Inhalt
Contents

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Carbide instruments

General information

Instruments are delivered not sterile unless labelled otherwise and must therefore be processed appropriately prior to use.

According to KRINKO BfArm (which is the Commission for Hospital Hygiene and Infection Prevention of the German Federal Institute for Drugs and Medical Devices) a classification in the Critical A risk group is recommended; country-specific requirements must be observed.

The frequency of repeat deployment of the device is at the user's discretion and should be based on the respective application. Damaged and worn instruments should be promptly discarded.

Only use cleaning and disinfection equipment that is compliant with EN ISO 15883 and an autoclave with a fractionated vacuum cycle that is compliant with DIN EN 13060/DIN EN 285/ANSI AAMI ST79 and validated as per DIN EN ISO 17665.

Only use suitable and approved cleaning agents and disinfectants according to the manufacturer's instructions.

Machine processing is preferable. A standardised procedure is specified for manual processing.

Storage/Pretreatment:

Remove visible residues with a soft nylon brush under running cold water for at least 1 minute to avoid their build-up.

Then place the instruments in a drill bath (preferably with an ultrasonic insert) with a non-protein-fixing and non-aldehyde-containing disinfecting and cleaning solution for the specified time.

Remove instruments from the bath and rinse thoroughly under running water. Check the instruments for any visible residues. If any residues still persist (particularly bone or dentine particles) either repeat the pre-cleaning process or dispose of the instrument.

Machine cleaning and disinfection:

Instruments should be rinsed under running water for at least 1 minute immediately before machine cleaning.

Use a suitable instrument stand to load instruments into the cleaning and disinfection machine to avoid damaging the devices. Run the following programme:

- 2 min pre-cleaning with cold tap water (10°C - 25°C)
- Drain
- 5-minute cleaning at 55°C with tap water and an alkaline cleaner
- Drain
- 3 min additional rinse with cold deionised water (10°C - 25°C)
- Drain
- 2 min additional rinse with cold deionised water (10°C - 25°C)
- Drain
- Thermal disinfection: at least 5 minutes at 90°C

When finished, remove instruments and dry (preferably with compressed air, as recommended by the Robert Koch Institute).

The basic compatibility of devices for effective machine cleaning was demonstrated in a Miele G 7836 CD washer/disinfector using an alkaline detergent (0.2% Neodisher FA, Dr. Weigert, Hamburg; 0.5%, Neodisher MediZym Dr. Weigert, Hamburg).

Manual cleaning and disinfection:

Rinse instruments under running cold water.

Heavily soiled instruments should be pre-cleaned manually by soaking in cold water (for at least 5 minutes) and then removing any visible debris with a soft nylon brush. Visually inspect at 10x magnification. If instruments are still visibly soiled, repeat the process until all traces of contamination have been removed.

Place instruments in an appropriate basket in an ultrasonic bath containing a cleaning and disinfecting agent.

Instruments should be thoroughly rinsed under running water (preferably distilled water) upon completion of the specified exposure time and then immediately dried (preferably with compressed air, as recommended by the Robert Koch Institute).

The basic compatibility of devices for effective machine cleaning was demonstrated in an ultrasonic bath using the cleaning and disinfecting agent (0.5%, Neodisher MediZym, Dr. Weigert, Hamburg) at 40°C for 5 minutes in an ultrasonic bath.

Visual inspection:

Visually inspect for cleanliness and integrity; we recommend using a 10x magnifying glass. Repeat the cleaning process to remove any remaining traces of contamination. Please dispose of the instruments if contamination is still visible after repeated cleaning and disinfection. Damaged and worn instruments should be discarded.

Packaging:

Use suitable steam sterilisation packaging (temperature resistant up to at least 138 °C (280 °F) and good steam permeability) with dimensions which ensure that the seal is not under tension, and that complies with DIN EN ISO/ANSI AAMI ISO 11607-1.

Sterilisation:

Autoclave using the following parameters:

- Fractionated prevacuum (3x)
- Sterilisation temperature: 134°C
- Holding temperature: 5 min (full cycle)
- Drying time: 10 min*

*The actual drying time required is directly dependent on the parameters for which the user bears sole responsibility (e.g. load configuration and density, condition of the autoclave) and must therefore be determined by the user. Nevertheless, drying times of 20 minutes should not be exceeded.

The basic supportive proof to determine that the devices can be effectively steam-sterilised was established using a Webeco A65-1 autoclave.

Storage:

Prior to initial use, devices should be stored in their original packaging at room temperature, and in an environment protected from dust and moisture.

Instruments should be stored and transported in hygienically maintained trays, racks, or other suitable containers. Stored items must be protected from dust, moisture, and recontamination.

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