










serres

ACCESSORIES

EN	ACCESSORIES.....	4
DE	ZUBEHÖR	6
FR	ACCESSOIRES.....	8
ES	ACCESORIOS	10
IT	ACCESSORI	12
SV	TILLBEHÖR.....	14
DA	TILBEHØR.....	16
NO	TILBEHØR	18
FI	LISÄVARUSTEET	20
NL	ACCESSOIRES.....	22
CS	RŘÍSLUŘENSTVÍ.....	24
SK	RRÍSLUŘENSTVO.....	26
TR	AKSESUARLAR.....	28
HU	KIEGÉSZÍTŐK.....	30
PL	AKCESORIA.....	32
PT	ACESSÓRIOS	34
SL	DODATKI	36
EL	ΠΑΡΕΛΚΟΜΕΝΑ	38
TH	สํานประกอบตําง ๑.....	40
LT	PRIEDAI	42

EN ACCESSORIES

(Availability must be checked from the reseller)

Serres Collection Cup ☒

Installed in a suction bag, a collection cup **Figure 7** is used to collect particles from the suctioned liquid. Disconnect the white angle connector of the suction bag and attach the collection cup to the patient connector. Connect the patient tube to the angle connector of the cup lid. During the procedure, the particles are collected into the basket of the cup.

After the particle collection, detach the patient tube and angle connector and close the connection (smaller plug on the cup lid). Detach the collection cup from the suction bag and turn the cup upside down.

Specimen jars are recommended to use to retain or transport specimens.

Serres Measuring Cup ☒ CE 0598

A measuring cup is used with suction bags when the precise volume of the suctioned liquid is needed. The device is intended to measure quantity of liquids and secretions removed from the human body. The scale on the measuring cup displays the volume of the suctioned substance in milliliters. The accuracy of the measuring cup is ± 5 ml up to 50 ml and then ± 10 ml. Hang the measuring cup in an upright position on the lid of a suction bag and connect the tubes according to **Figure 8**.

Connect the patient tube to the angle connector in the middle of the measuring cup. During the procedure, the cup can be emptied into the suction bag without disconnecting the vacuum. To do this, the cup is lifted from the bracket and tilted.

Serres Splash-Vac ☒

The Serres Splash-Vac **Figure 9** can be used to suction liquids from the floor. Connect the tubing to the patient connector of the suction bag.

NOTE! THE TUBING OF SPLASH-VAC IS MADE OF PVC PLASTIC. The product is not a medical device and thus does not have a CE marking.

Straight connector ☒

A straight connector **Figure 10** is used when there is a need to use wider-bore patient tubing. The patient connector in the suction bag can be replaced by the straight connector.

Twin Connector ☒

A twin connector **Figure 11** can be used to connect two patient tubes to the suction bag simultaneously. During suction, one patient tube is kept closed. The patient connector on the suction bag is replaced with the twin connector.

Smoke filter ☒

The smoke filter **Figure 12** can be used to prevent the premature blockage of the suction bag overflow protection in surgical operations during which there is very heavy smoke production.

NOTE! THE SMOKE FILTER IS NOT INTENDED FOR SUCTIONING SMOKE AND IT DOES NOT PREVENT SMOKE FROM RE-ENTERING THE ROOM AIR.

Replace the suction canister grey, angled connector with a T-connector or Serres valve **Figure 16** and open the suction bag serial port. Press the smoke filter into the serial port so that the filter locks in place.

Connect the smoke filter tubing to the T-connector or valve now on the rear of the canister. Different positions of the Serres valve:

12A Suction is closed

12B Suction only through the bag filter (used when no smoke is created)

12C Suction through bag filter and smoke filter (used when smoke is created)

The smoke filter can also be used in serial connection. In such a case, install the smoke filter in the last suction bag of the serial connection.

WARNING! REPLACE T-CONNECTOR WITH ANGLED CONNECTOR OR TURN THE SERRES VALVE TO POSITION AS SHOWN IN PICTURE **12B** IF THE SMOKE FILTER IS NO LONGER USED.

Serres Vacuum Shift

The Serres Vacuum Shift is mounted to the device bracket between two suction canisters. The Vacuum Shift is designed for shifting the vacuum between two suction lines or for the simultaneous use of two suction bag systems. The suction canisters may also be easily connected on each side in a serial chain.

1. Mount the Vacuum Shift on the device bracket, using the connector on the Vacuum Shift

2. Connect the tube of the suction source to the tube connector of the Vacuum Shift, (**figure 13, connector A**)

3. Connect the tubing from the suction canisters to the tube connectors on the right side (**figure 13, connector B**) and left side (**figure 13, connector C**) of the Vacuum Shift.

4. The Vacuum Shift has four different positions.

NOTE! ENSURE THAT THE TIP OF THE LATCH POINTS DIRECTLY TO THE SIDE OR UPWARDS OR DOWNWARDS, AS INDICATED IN THE IMAGES. AN INCORRECT ALIGNMENT MAY WEAKEN THE SUCTION EFFECT.

- figure 13.1**, suction closed
- figure 13.2**, suction from the right
- figure 13.3**, suction from the left
- figure 13.4**, suction from the right and left

Serres recommends the use of silicone tubing between the suction canister/source and the Vacuum Shift.

Serres Vacuum source protector

The vacuum source protector **figure 14** prevents liquid entering the vacuum source in case of an installation error of a suction bag. Connect the vacuum source protector between the suction canister and the vacuum source. Ensure the arrow is pointing to the vacuum source. The protector is replaced every 30 days or if the suction power reduces or if contaminated.

Serres on-off valve

The on-off valve is used to open and shut off suction. Install the on-off valve in the suction canister connection to replace the grey, angled connector. Connect the tubing from the suction source to the on-off valve connector. To open (**figure 15 A**) or shut off (**figure 15 B**) the suction, turn the on-off valve lever.

Serres Valve

The Serres valve **figure 16** is used to open and shut off suction in connection to using a smoke filter. For further instructions, please refer to the smoke filter section.

Serres T-valve

The T-valve **figure 17** is used to open and shut off suction in serial connection. Install the T-valve in the suction canister connection to replace the grey, angled connector.

Serres Control valve

Serres control valve facilitates the speedy individual use of a Serres suction bag. The Serres control valve is installed inside the canister holder of a Serres trolley (ref: 57940, 57941).

Installation:

- Detach the canister holder from the trolley.
- Open the holes (6 pcs) of on-off valves on the canister holder's top surface by driving the cover plates in. **Figure 18.1**.
- Install the control valve inside the canister holder and rotate the control valve until it reaches the bottom.
- Use three screws to attach the control valve to the canister holder.
- Re-install the canister holder to the trolley, adjust the height of the holder and tighten the canister holder, **figure 18.2**
- Install as many on-off valves to the holes opened on top of the canister holder as there are suction canisters. If you commission less than 6 suction canisters the extra holes must be closed with a control valve plug (ref. 57943). In order for the suction to be created correctly, all six holes must either be fitted with an on-off valve or a control valve plug.
- Install the silicone tubing from the on-off valve to the angle connector of the suction canister below. Install the tubing so that it is not bent and does not prevent air from flowing.
- Install the tube of the suction device to the tube connector on the bottom surface of the control valve. **Figure 18.2. C**
- Test the functionality of the product by ensuring one by one that each on-off valve is tight.

Use:

- Start the suction source
- Install the required number of suction canisters and suction bags.
- Open one on-off valve to position I
- Connect the patient tube to the suction bag below the on-off valve
- When the suction bag is full, remove the patient tube, plug the full suction bag and turn the on-off valve to position O.
- Turn the next on-off valve to position I and connect the patient tube to the suction bag below it. In order to achieve sufficient suction power, it is recommended that only one suction canister at a time should be used, one on-off valve is in the I position and the other on-off valves are in position O.

Disposal of used products

Serres suction bags, collection cups, measuring cups and serial tubes are disposable and must be replaced patient-specifically. If the same patient undergoes long-term treatment, it is recommended that the suction bag be replaced at least every 24 hours. Reusable products must be decommissioned if they have become damaged or no longer meet the manufacturer's intended

operational specifications. Used or decommissioned products must be disposed of in accordance with recycling labelling and hospital-specific instructions. Suction bags containing solidifying agent must not be emptied into a sewer.

WARNING! REUSE OF DISPOSABLE PRODUCTS IS STRICTLY FORBIDDEN. REUSE REDUCES THE PERFORMANCE OF THE PRODUCT AND CAN CAUSE AN INFECTION RISK. WHEN HANDLING A USED SUCTION BAG, ATTENTION SHOULD BE PAID TO THE FACT THAT IT MAY CONTAIN INFECTIOUS WASTE.

Cleaning of re-usable products

The suction canister and angle connector can be washed (95°C) and autoclaved (121 °C). Remove the grey, angled connector before washing or autoclaving. Before autoclaving, make sure that the product does not contain any detergent remnants. Do not autoclave

the product more than 30 cycles. Adapters and brackets (excluding 57815) can be washed (95°C). Other reusable products can be wiped with disinfectant.

Storage

Protect the packages from humidity, dirt and dust. Disposable products can be used for 5 years after the date on the label, excluding pre-gelled suction bags and solidifying agents that can be used for 2 years after the date on the label.

Reporting of serious incidents

Any serious incident related to the use of this product, should be reported to both the manufacturer and the health authority/competent authority where the product is installed.

Symbols used

 Batch code

 Catalogue number

 Do not re-use

 Keep away from rain

 Keep away from sunlight

 Operating instruction

 Medical Device

MFGDT (Manufacturing date)

 CE-mark

 DEHP-Free symbol

 Latexfree


 Quantity

 Refer to instruction manual booklet

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Serres Oy