

# Zeta 7 Solution



## Special applications

A concentrated disinfectant for impressions.

**Active ingredients:** quaternary ammonium salts, phenoxyethanol.

**1% DILUTION**



## PRODUCT BENEFITS

- ▶ High protection: bactericidal, yeasticidal, tuberculocidal and virucidal against enveloped viruses.
- ▶ Compatibility with impression materials (addition and condensation silicones, alginate, polyether, polysulphide and polyvinyl).
- ▶ High efficiency concentrated formula.
- ▶ Respects the dimensional stability of impressions and their compatibility with gypsum.
- ▶ Safe to use: free from aldehydes and phenols.

## Use

**1% dilution:** add 10 ml of **Zeta 7 Solution** to every litre of water. Rinse the impression under running water. Immerse the impression in the disinfectant solution for 10 minutes. Remove the impression from the bath and rinse it thoroughly.

**ONE LITRE OF ZETA 7 SOLUTION GIVES 100 LITRES OF DISINFECTANT SOLUTION. (UP TO 40 IMPRESSIONS CAN BE DISINFECTED WITH ONE LITER OF SOLUTION\*).**

\*The solution must be replaced every day and immediately in case of visible contamination (blood, saliva or organic tissue).

## Spectrum



**Bactericidal:** EN 13727 (*S. aureus*, *P. aeruginosa*, *E. hirae*)

**Yeasticidal:** EN 13624 (*C. albicans*)

**Tuberculocidal:** EN 14348 (*M. terrae*)

**Virucidal against enveloped viruses:** EN 14476 (Vaccinia Virus including HIV, HBV, HCV)

*Tests carried out in dirty conditions.*