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RECOMMENDED INDICATIONS
GC Fuji PLUS CONDITIONER is a tooth conditioning agent designed for use with GC Fuji PLUS glass ionomer luting cement. It helps to enhance the bond between GC Fuji PLUS and tooth structure.

CONTRAINDICATIONS
In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

DIRECTIONS FOR USE
1. After tooth preparation, apply GC Fuji PLUS CONDITIONER to the bonding surface for 20 seconds using a cotton pellet or sponge.
2. Rinse thoroughly with water and dry. DO NOT DESICCATE.
3. Proceed according to GC Fuji PLUS directions.
4. Close the bottle immediately after use.

STORAGE
Store at temperature of 4-25°C (39.2-77.0°F).
(Shelf life: 2 years from date of manufacture)

PACKAGES
Bottle of 7g (6.5mL)

CAUTION
1. In case of contact with oral tissue or skin, flush immediately with water.
2. In case of contact with eyes, flush immediately with water and seek medical attention.
3. GC Fuji PLUS CONDITIONER contains 2% Ferric Chloride. Not recommended for use with products that contain Epinephrine, which is known to react with iron containing products and may stain oral tissue.

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RECOMMENDED INDICATIONS

GC Fuji PLUS CONDITIONER is a tooth conditioning agent designed for use with GC Fuji PLUS glass ionomer luting cement. It helps to enhance the bond between GC Fuji PLUS and tooth structure.

CONTRAINDICATIONS

In rare cases the product may cause sensitivity to some persons. If such reactions are experienced, discontinue the use of the product and refer to a physician.

DIRECTIONS FOR USE

- After tooth preparation, apply GC Fuji PLUS CONDITIONER to the bonding surface for 20 seconds using a cotton pellet or sponge.
- Rinse thoroughly with water and dry. DO NOT DESICCATE.
- Proceed according to GC Fuji PLUS directions.
- Close the bottle immediately after use.

STORAGE

Store at temperature of 4-25°C (39,2-77,0°F).
(Shelf life: 2 years from date of manufacture)

PACKAGES

Bottle of 7g (6.5mL)

CAUTION

- In case of contact with oral tissue or skin, flush immediately with water.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- GC Fuji PLUS CONDITIONER contains 2% Ferric Chloride. Not recommended for use with products that contain Epinephrine, which is known to react with iron containing products and may stain oral tissue.

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Viimati parandatud : 05/2008

GC Fuji PLUS EWT

RADIOPAQUE REINFORCED GLASS IONOMER LUTING CEMENT WITH EXTENDED WORKING TIME

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DIRECTIONS FOR USE

Powder / Liquid ratio (g/g)	2,0/1,0
Mixing time (sec.)	20"
Working time (min., sec.) (at 23 °C) (from start of mixing)	3' 30"
Final finishing commencing time (after seating the restoration)	5' 45"

- TOOTH PREPARATION
 - Prepare tooth in usual manner. For pulp capping, use calcium hydroxide.
 - Clean the prepared teeth with pumice and water. Note: Use GC Fuji PLUS CONDITIONER on the prepared tooth surface for 20 seconds for optimal adhesion.
 - Rinse thoroughly with water. Remove excess moisture by blotting with a cotton pellet or gently blowing with an oil-free air. DO NOT DESICCATE. Prepared surface should appear moist (glistening).
- RESTORATION PREPARATIONMake sure that the restoration is pretreated and handled according to the manufacturer's instructions.
- POWDER AND LIQUID DISPENSINGNote : Use GC Fuji PLUS Liquid. Refer to the instructions for use of GC Fuji PLUS.
 - The standard powder to liquid ratio is 2.0g / 1.0g (1 level large scoop of powder to 3 drops of liquid, or 1 level small scoop of powder to 1 drop of liquid).
 - For accurate dispensing of powder, lightly tap the bottle against the hand. Do not shake or invert.
 - Hold the liquid bottle vertically and squeeze gently.
 - Close bottles immediately after use.
- MIXINGDispense powder and liquid onto the pad. Using the plastic spatula, add all the powder to the liquid. Mix rapidly for 20 seconds. Note : When mixing larger amounts, divide the powder into 2 equal parts. Mix the first portion with all of the liquid for 15 seconds. Incorporate the remaining portion and mix the whole thoroughly for 15-25 seconds (total 30-40 seconds).
- CEMENTATION
 - Coat the internal surface of the restoration with sufficient cement and seat immediately. The working time is 3 minutes 30 seconds after start of mixing at 23°C (73.4°F). Higher temperatures will shorten working time.
 - Maintain moderate pressure.
 - Start removing excess cement when the excess cement feels rubbery.
 - Finishing can be started 5 minutes 45 seconds after seating the restoration.

STORAGE

Store in a cool and dark place (4-25°C) (39.2-77.0°F). (Shelf life: 3 years from date of manufacture)

SHADE

A3 (Vita® shade)
Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.

PACKAGE

Powder 15g with scoop

CAUTION

- In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water. To avoid contact, a rubber dam and/or cocoa butter can be used to isolate the operation field from the skin or oral tissue.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- DO NOT mix powder or liquid with components of other glass ionomer cements.
- GC Fuji PLUS CONDITIONER contains 2% Ferric Chloride. It is not recommended to use GC Fuji PLUS CONDITIONER in combination with products containing Epinephrine as this may stain oral tissue.
- Do not use GC Fuji PLUS EWT in combination with eugenol containing materials as eugenol may hinder GC Fuji PLUS EWT from setting or bonding properly.
- This product is not indicated for filling or core build-up.

Last revised : 05/2008

0086

'GC.'

MANUFACTURED BY
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

DISTRIBUTED BY
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

 EU: GC EUROPE N.V.
Research-Park, Interleuvenlaan 13, B-3001 Leuven, Belgium
TEL: +32. 16. 39. 80. 50.

GC AMERICA INC.
3737 West 127th Street, Alsip, IL 60803 U.S.A.
TEL: +1-708-597-0900

GC ASIA DENTAL PTE. LTD.
19 Loyang Way, #06-27 Singapore 508724
TEL: +65 6546 7588

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Prior to use, carefully read the instructions for use.

EN

GC Fuji PLUS

RADIOPAQUE REINFORCED GLASS IONOMER LUTING CEMENT

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

- Cementation of metal-based inlays, onlays, crowns and bridges
- Cementation of resin inlays, onlays, crowns and bridges
- Cementation of all ceramic inlays
- Cementation of high strength (zirconia based) all ceramic crowns and bridges
- Cementation of metal, ceramic and fiber posts

CONTRAINDICATIONS

- Pulp capping.
- In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

DIRECTIONS FOR USE

Powder / Liquid ratio (g/g)	2.0/1.0
Mixing time (sec.)	20"
Working time (min., sec.) (at 23 °C) (from start of mixing)	2' 30"
Final finishing commencing time (after seating the restoration)	4' 30"

1. TOOTH PREPARATION

- Prepare tooth in usual manner. For pulp capping, use calcium hydroxide.
- Clean the prepared tooth with pumice and water.

Note:
Use GC Fuji PLUS CONDITIONER on the prepared tooth surface for 20 seconds for optimal adhesion.
- Rinse thoroughly with water. Remove excess moisture by blotting with a cotton pellet or gently blowing with an oil-free air. DO NOT DESICCATE. Prepared surfaces should appear moist (glistening).

2. RESTORATION PREPARATION

Make sure that the restoration is pretreated and handled according to the manufacturer's instructions.

3. POWDER AND LIQUID DISPENSING

- The standard powder to liquid ratio is 2.0g / 1.0g. (1 level large scoop of powder to 3 drops of liquid, or 1 level small scoop of powder to 1 drop of liquid).
- For accurate dispensing of powder, lightly tap the bottle against the hand. Do not shake or invert.
- Hold the liquid bottle vertically and squeeze gently. Close bottles immediately after use.

4. MIXING

Dispense powder and liquid onto the pad. Using the plastic spatula, add all the powder to the liquid. Mix rapidly for 20 seconds.

Note :

When mixing larger amounts, divide the powder into 2 equal parts. Mix the first portion with all of the liquid for 15 seconds. Incorporate the remaining portion and mix the whole thoroughly for 15-25 seconds (total 30-40 seconds).

5. CEMENTATION

- Coat the internal surface of the restoration with sufficient cement and seat immediately. The working time is 2 minutes 30 seconds from start of mixing at 23°C (73.4°F). Higher temperature will shorten working time.

Note:
For more complex procedures requiring longer working time, it is recommended to mix powder and liquid on a cold mixing slab and or use GC Fuji PLUS EWT.
- Maintain moderate pressure.
- Start removing excess cement when the excess cement feels rubbery.
- Finishing can be started 4 minutes 30 seconds after seating the restoration.

STORAGE

Store in a cool and dark place (4-25°C) (39.2-77.0°F). (Shelf life: powder 3 years, liquid 2 years from date of manufacture)

SHADES

Yellow
A3 (Vita® shade)
Translucent

Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.

PACKAGES

- 1-1 package : 15g powder, 8g (7mL) liquid, GC Fuji PLUS CONDITIONER 7g (6.5mL), a double-ended powder scoop, a plastic spatula, mixing pad (#23)
- Powder 15g with scoop
- Liquid 8g (7mL)
- GC Fuji PLUS CONDITIONER 7g (6.5mL)

CAUTION

- In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water. To avoid contact, a rubber dam and/or cocoa butter can be used to isolate the operation field from the skin or oral tissue.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- DO NOT mix powder or liquid with components of other glass ionomer cements.
- GC Fuji PLUS CONDITIONER contains 2% Ferric Chloride. It is not recommended to use GC Fuji PLUS CONDITIONER in combination with products containing Epinephrine as this may stain oral tissue.
- Do not use GC Fuji PLUS in combination with eugenol containing materials as eugenol may hinder GC Fuji PLUS from setting or bonding properly.
- This product is not indicated for filling or core build-up.

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