VOCO Grandio® Seal

EN Instructions for use

MD EU Medical device

In accordance with DIN EN ISO 6874

Product description:

Grandio Seal is a light-curing, highly filled, fluoride-containing nano-hybrid fissure sealant with outstanding mechanical properties.

Grandio Seal contains 70 % w/w inorganic fillers in a methacrylate matrix and cures under blue light (halogen or LED).

Indications:

- Sealing/filling of pits and fissures
- Sealing/facing of damaged enamel surfaces
- Covering of caries predilection sites during orthodontic treatments
- Sealing of deciduous teeth

Contraindications:

Grandio Seal contains aromatic and aliphatic dimethacrylates as well as BHT. **Grandio Seal** should therefore not be used in patients with a known hypersensitivity (allergy) to these ingredients.

Patient target group:

Grandio Seal is suitable for use in all patients without any age or gender restrictions.

Performance features:

The product's performance features satisfy the requirements of the intended use and the relevant product standards.

User:

Grandio Seal should only be used by a professionally trained dental practitioner.

Use:

Fissures should be sealed as early as possible, as soon as the eruption of the tooth allows treatment and moisture control is possible. After isolating the operating site, thoroughly clean the surfaces to be sealed (rotating brush/cleaning paste), spray and dry.

In case of incipient carious lesions or unclear diagnosis, extend the fissure with a fine instrument, if necessary.

Etching:

Apply **Vococid** etching gel (or liquid) (blue colour for control), work in with a short-haired disposable brush and let it act for 15 - 30 s according to the condition of the enamel. Thoroughly spray and dry. The etched enamel should have an opaque and matt-white appearance (avoid contamination at all costs!). Follow the **Vococid** instructions for use.

Application:

Remove the protective cap from the syringe. Screw a new cannula onto the **Grandio Seal** syringe.

Grandio Seal flows homogeneously from the cannula upon medium pressure and this is very easy to dispense accurately.

Apply **Grandio Seal** from the margin into the fissure or cavity and allow it to penetrate for 15 - 20 s. Remove any air bubbles with the fine cannula.

In order to ensure the functioning of the non-dripping syringe do not pull back the plunger during or after use.

The attached application tip is intended for single use only. Remove it after use and seal the syringe again tightly with its original cap. The use of a protective barrier is recommended to protect the syringe against contamination with bodily fluids, dirt from hands or oral tissue.

Light-curing:

Polymerise **Grandio Seal** for 20 s per surface and check the sealing. Conventional polymerisation devices are suitable for light-curing this material. The light output should be at least 500 mW/cm² on halogen polymerisation devices and 300 mW/cm² on LED devices. Re-seal undercuts, remove excess cautiously (e.g. finishing diamond), polish if necessary (rubber polisher). Then apply a fluoride varnish on the treated and neighbouring teeth.

Warnings, precautionary measures:

- An operating site which cannot be kept sufficiently dry involves the risk of a (partial) premature loss of the sealing. In such cases, the application of a fluoride varnish or sealing with a glass ionomer cement is recommended until perfect sealing is possible.
- If the etched surface is contaminated (e.g. by saliva), the adhesion of **Grandio Seal** may be reduced. In case of contamination, clean and etch again for approx. 10 - 20 s, spray and dry.
- Adhesion can be increased by using a bonding agent after etching (e.g. in case of difficult moisture control/extended sealing). Follow the bonding agent manufacturer's instructions for use. Then apply **Grandio Seal**.
- Only expose **Grandio Seal** to the operating light for a short time as otherwise polymerisation will start (viscosity increases).
- Check the sulcus after application and remove residues of material.
- Avoid contact with the mucosa as far as possible. Do not use cleaning or drying agents containing phenolic substances (e.g. eugenol, thymol) as they inhibit complete polymerisation.
- Carious lesions (e.g. around fillings) can be stopped by appropriate sealing if renewal of the filling is not desired or possible. For this purpose, clean the areas affected, prepare as minimal as possible and etch in the enamel area. In such cases, sealing of the complete filling is recommended.
- In the case of orthodontic treatment in areas of reduced oral hygiene, sealing these areas can provide good and permanent protection. Fluoridation should be continued in all cases to protect the unsealed tooth areas from caries.
- Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for their intended purpose.

Constituents (in descending order according to content):

Barium aluminium borosilicate glass, TEGDMA, silicon dioxide, BisGMA, fumed silica, PEGDMA, initiators, stabilisers, pigments

Storage:

Store **Grandio Seal** away from light at temperatures between 4 °C and 23 °C. Do not use after the expiry date.

Disposal:

Dispose of the product in accordance with local regulations.

Reporting obligation:

Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's health condition and a serious risk to public health that arise or could have arisen in association with the use of **Grandio Seal** must be reported to VOCO GmbH and the responsible authority.

Note:

The Summary of Safety and Clinical Performance of **Grandio Seal** can be found in the European database on medical devices (EUDAMED – <https://ec.europa.eu/tools/eudamed>). Detailed information can also be found at www.voco.dental.

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VC 60 AA1061 E1 1025 99 © by VOCO