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# VOCO Grandio®SO Flow Grandio®SO Heavy Flow

<b>EN</b>	<b>Instructions for use</b>
<b>MD</b>	<b>EU Medical device</b>
	<b>In accordance with DIN EN ISO 4049</b>

**Product description:**

**GrandioSO Flow / Heavy Flow** is a flowable composite. In contrast to conventional flows, it is characterised by its increased filler content and stability. The radiopaque restorative contains 81 (**GrandioSO Flow**), 82 (**GrandioSO Flow WO**) and 83 % w/w (**GrandioSO Heavy Flow**) inorganic fillers in a methacrylate matrix and cures under blue light. In **GrandioSO Flow WO**, an extremely radiopaque version is at your disposal. **GrandioSO Flow / Heavy Flow** is used with a dentine-enamel bond.

**Shades:**

A1, A2, A3, A3.5, A4, B1, B2, C2, OA2, BL, WO, Incisal, GA5

**Indications:**

- Fillings of class I - V cavities
- Minimally invasive filling therapy
- Extended fissure sealing
- Blocking out undercuts
- Lining or coating cavities
- Repairing fillings and veneers
- Luting translucent prosthetic pieces (e.g., full ceramic crowns, etc.)
- Interlocking and splinting teeth (e. g., with glass fibre strands)
- Use as a base material in combination with glass fibres or similar products for fabricating semi-permanent crowns and bridges

**Contraindications:**

**GrandioSO Flow / Heavy Flow** contains methacrylate and BHT. **GrandioSO Flow / Heavy Flow** should therefore not be used in patients with a known hypersensitivity (allergy) to these ingredients.

**Patient target group:**

**GrandioSO Flow / Heavy Flow** is suitable for use in all patients without any age or gender restrictions.

**Performance features:**

The product’s performance features satisfy the requirements of the intended use and the relevant product standards.

**User:**

**GrandioSO Flow / Heavy Flow** should only be used by a professionally trained dental practitioner.

**Use:**

**Preparation:**

Clean the teeth being treated. If applicable, mark occlusal contact points. Bring material to room temperature before use.

**Shade selection:**

To select a shade, compare the enclosed shade guide with the cleaned tooth before anaesthesia. Moisten the shade sample and tooth and select the shade, preferably in daylight.

**Creating a dry working field:**

Ensure that the working field is sufficiently dry. Use of a rubber dam is recommended.

**Cavity preparation:**

In principle, cavity preparation should be carried out according to the rules of adhesive filling therapy and should be minimally invasive to conserve healthy tooth substance. Bevel the enamel margins on anteriors. Play attention to special preparation shapes for extended fissure sealing and deciduous teeth! Clean and dry the cavity afterwards. Non-carious cervical lesions do not have to be prepared; thorough cleaning is sufficient here.

**Lining:**

A suitable pulp-protective liner that can be covered with a stable cement, if necessary, should be applied in close proximity to the pulp.

**Matrix placement:**

Take appropriate measures for areas where the applied **GrandioSO Flow / Heavy Flow** can flow out. Translucent matrices that are wedged in the approximal area are advantageous. Minimal separation facilitates the shaping of the approximal contact and placement of the matrix.

**Bonding material:**

**GrandioSO Flow / Heavy Flow** is used in the adhesive technique with a dentine-enamel bond. All light-curing bond materials can be used. Follow the respective instructions for use with regard to preparation (etching technique) and application.

**Application of GrandioSO Flow / Heavy Flow:**

**Syringe:** Lock the enclosed application tip on the syringe with a clockwise turn and directly apply **GrandioSO Flow / Heavy Flow**. For fillings with a thickness of over 2 mm, apply and cure in layers. Do not pull the plunger back on the syringe during or after use to guarantee the function of the non-dripping syringe. The attached application tip is intended for single use only. Remove it after use and seal the syringe again tightly with its original cap. The use of a protective barrier is recommended to protect the syringe against contamination with bodily fluids, dirt from hands or oral tissue.

**Caps:** Remove protective cap. Insert Caps into the opening of the dispenser (follow the instructions for use). Bend the metal cannula according to the indication into the desired direction. Apply **GrandioSO Flow / Heavy Flow** directly into the cavity by slightly and evenly pressing the levers of the dispenser together.

Use each **GrandioSO Flow / Heavy Flow** Caps for one patient only!

**Light-curing:**

Conventional polymerisation devices are suitable for light-curing the material. Curing time if an LED/halogen lamp with a light output of at least 500 mW/cm² is used is

- 20 s for shades A1, A2, A3, A3.5, B1, B2, C2, BL and Incisal
- 40 s for shades A4, OA2, WO and GA5

Bring the light aperture as close as possible to the light-curing unit as otherwise the curing will be less effective. Inadequate curing can cause discolouration and complaints.

**Finishing:**

The restoration can be finished and polished immediately after removing the matrices (e.g., fine or extra-fine diamond burs, polisher). The tooth should be fluoridated as a final step.

**Warnings, precautionary measures:**

Phenolic substances, especially products containing eugenol and thymol, lead to curing failures in restorative composites. Consequently, the use of zinc oxide/eugenol cements or other materials containing eugenol should be avoided in combination with restorative composites.

– Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for their intended purpose.

**Constituents (in descending order according to content):**

Barium aluminium borosilicate glass, silicon dioxide, HEDMA, BisGMA, TEGDMA, BisEMA, fumed silica, initiators, stabilisers, pigments

**Storage instructions and application method:**

Storage at 4°C-23°C. Seal the syringes immediately after removing the material in order to avoid exposure to light and subsequent polymerisation. Do not use after the expiry date.

**Disposal:**

Dispose of the product in accordance with local regulations.

**Reporting obligation:**

Serious events such as death, temporary or permanent serious deterioration of a patient’s, user’s or other person’s health condition and a serious risk to public health that arise or could have arisen in association with the use of **GrandioSO Flow / Heavy Flow** must be reported to VOCO GmbH and the responsible authority.

**Note:**

The Summary of Safety and Clinical Performance of **GrandioSO Flow / Heavy Flow** can be found in the European database on medical devices (EUDAMED – https:// ec.europa.eu/tools/eudamed). Detailed information can also be found at www.voco.dental.

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