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Lyostypt®

Description

Lyostypt® is a wet-stable, felt-like collagen haemostatic and is composed of native, absorbable collagen fibrils of bovine origin. The product keeps its shape, even when it is soaked with blood, and, thus, can be completely removed after bleeding has stopped, if so desired.

The felt structure being rich in surface offers an ideal framework for the adhesion of platelets and, thus, provides an additional impetus to clotting.

Lyostypt® is radiation-sterilized and is available in the form of rectangular pieces of various sizes.

Composition:

10 cm² contain 100 mg collagen of bovine origin.

Indications

For local haemostasis of capillary bleeding, bleeding from parenchymal organs, oozing haemorrhages and as a supportive measure for other techniques of haemostasis.

Mode of action

The contact between collagen and blood leads to platelet aggregation. The collagen structure leads to additional stabilization of the clot.

Contra-indications

- Lyostypt® must not be applied to infected areas.
- Lyostypt® should not be employed in connection with cemented endoprosthesis since this would reduce the adhesive capacity of the bone cement.
- Known hypersensitivity to proteins of bovine origin.

Warning notes

None.

Precautions in use

See under mode of application.

Interactions with other drugs

Reduces the bonding strength of bone cement.

The efficacy of Lyostypt® is reduced by aggregation inhibitors and coagulation-inhibiting drugs.

Mode of application

Lay a Lyostypt® compress of the appropriate size on the wound or introduce it into the wound. The compress may be cut to the desired size with scissors. The affected wound should be cleansed as far as possible of blood. Then the compress should be lightly pressed onto the wound using forceps or clamp until fibrin formation assures adhesion. Excess compress material should be removed after haemostasis is achieved. In the case of skin wounds Lyostypt® must not be allowed to penetrate the skin borders.

The Lyostypt® compress is designed to remain, in general, at the site of application, but be removed after bleeding has stopped. Absorption takes place within ca. 3 weeks as a result of phagocytosis and enzymatic degradation. The rate of absorption also depends on the amount of Lyostypt® applied.

Side effects

Adhesions can occasionally occur at the site of application after intra-abdominal use. Allergic reactions to collagen of bovine origin can occur.

Note and details of shelf-life

Lyostypt® should not be used after the expiry date given.

Lyostypt® should be used directly after opening the pack.

Lyostypt® should be stored at 25 ± 5 °C.

General information

In addition to its good haemostyptic properties Lyostypt® is characterized by good shaping properties, softness and absorbency. Lyostypt® possesses a high affinity for bloody surfaces. The adhesion of Lyostypt® to Instruments and gauze strips is slight.

Lyostypt® must not be resterilized.

Symbols used on labelling

Do not reuse!

Use until Year + Month

Temperature range for storage

Sterile unless package is opened or damaged.
Method of Sterilization: Gamma-Irradiation

CE-mark and identification number of notified body. Product conforms to the essential requirements of the Medical Device Directive 93/42/EEC.

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B. Braun Surgical, S.A.