**EN** Instructions for use

**MD** EU Medical device

In accordance with DIN EN ISO 4049

**Product description:**

**Grandio Flow** is a light-curing flowable restorative combining proven composite technology with innovative nano-technology. The radiopaque **Grandio Flow** contains 80% w/w inorganic fillers ( $\leq 65.6\text{ % Vol.}$ ) in a methacrylate matrix and cures under blue light (halogen/LED). It is suited for restorative therapy in anterior and posterior areas. Because of its low viscosity, the material exhibits very good wetting and adaptation to cavity walls. **Grandio Flow** is used with a dentine enamel bond.

**Shades:**

A1, A2, A3, A3.5, A4, B1, B2, C2, D2, OA2, BL, WO

**Indications:**

- Filling minimally invasive cavities of all classes
- Filling small class I cavities and extended fissure sealing
- Filling class II - V cavities including V-shaped defects and cervical caries
- Blocking out undercutts
- Lining or coating cavities
- Repairing fillings and veneers
- Luting translucent prosthetic pieces (e.g. full ceramic crowns, etc.)

**Contraindications:**

**Grandio Flow** contains methacrylate and BHT. **Grandio Flow** should therefore not be used in patients with a known hypersensitivity (allergy) to these ingredients.

**Patient target group:**

**Grandio Flow** is suitable for use in all patients without any age or gender restrictions.

**Performance features:**

The product's performance features satisfy the requirements of the intended use and the relevant product standards.

**User:**

**Grandio Flow** should only be used by a professionally trained dental practitioner.

**Use:**

**Grandio Flow** is a material which adapts very well to the cavity, helps to avoid air enclosures and accelerates restorations.

**Preparation:**

Clean the teeth being treated. If applicable, mark occlusal contact points. Bring material to room temperature before use.

**Shade selection:**

Clean the teeth prior to shade selection. The shade is selected in comparison with the tooth while it is still moist.

**Cavity preparation:**

In principle, the cavities should be prepared according to the rules for adhesive restorative treatment; preparation should be minimally invasive in order to protect healthy dental hard tissue.

Bevel enamel margins. Observe special preparation forms in extended fissure sealing and deciduous teeth! Then clean and dry the cavity. Caries-free lesions in the cervical area do not need to be prepared; thorough cleaning is sufficient in this case.

**Creating a dry working field:**

Ensure that the work area is sufficiently dry. Use of a rubber dam is recommended.

**Matrix placement:**

Take appropriate measures for areas where the applied **Grandio Flow** can flow out. Translucent matrices that are wedged in the approximal area are advantageous. Minimal separation facilitates the shaping of the approximal contact and placement of the matrix.

**Pulp protection:**

A suitable pulp-protective liner that can be covered with a stable cement, if necessary, should be applied in close proximity to the pulp.

**Bonding materials:**

**Grandio Flow** is applied with a dentine enamel bond, using the adhesive technique. Any light-curing bonding materials can be used. Observe the relevant instructions for use during preparation (etching technique) and processing.

**Application:**

Lock the application cannula (type 41) onto the syringe with a clockwise turn and apply **Grandio Flow** directly. For fillings of more than 2 mm, apply and polymerise in layers. In order to ensure the function of the non-dripping syringe do not pull back the plunger during or after use.

The attached application tip is intended for single use only. Remove it after use and seal the syringe again tightly with its original cap. The use of a protective barrier is recommended to protect the syringe against contamination with bodily fluids, dirt from hands or oral tissue. Insert Caps into the opening of the dispenser following the respective instructions for use. Turn Caps to the desired direction and remove the protective cover. Apply **Grandio Flow** directly into the cavity by slowly and evenly pressing the levers of the handle together.

**Light-curing:**

Conventional polymerisation devices are suitable for light-curing the material. Halogen light devices should have a minimum light output of 500 mW/cm<sup>2</sup>, LED lamps should have a minimum output of 300 mW/cm<sup>2</sup>.

The curing time is 20 s minimum per layer, for opaque shades A4, OA2 and WO a minimum of 40 s. Hold the tip of the curing light as close as possible to the surface of the filling. If the distance is more than 5 mm, the curing depth will be reduced. Insufficient curing can lead to discolouration and pulp irritation.

**Finishing:**

The filling can be finished and polished immediately after removing the matrices (e.g. fine and extra-fine grain diamonds, polishing discs). During finishing and polishing, the filling must be cooled. The margin of the filling, or even better, the entire tooth should be fluoridated as a final step.

**Warnings, precautionary measures:**

- No known side effects. Hypersensitive patients may experience sensitivity.
- Each **Grandio Flow** Caps unit should be used for one patient only and should be applied exerting uniform but not excessive pressure.
- Phenolic substances, especially preparations containing eugenol or thymol, inhibit the curing of **Grandio Flow**. Therefore, do not use zinc oxide eugenol cements or other eugenol-containing substances in combination with **Grandio Flow**.
- Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for their intended purpose.

**Constituents (in descending order according to content):**

Strontium aluminium borosilicate glass, Barium aluminium borosilicate glass, silicon dioxide, HEDMA, BisGMA, TEGDMA, fumed silica, initiators, stabilisers, pigments

**Storage instructions and application method:**

Storage at 4°C-23°C. Seal the syringes immediately after removing the material in order to avoid exposure to light and subsequent polymerisation. Do not use after the expiry date.

**Disposal:**

Dispose of the product in accordance with local regulations.

**Reporting obligation:**

Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's health condition and a serious risk to public health that arise or could have arisen in association with the use of **Grandio Flow** must be reported to VOCO GmbH and the responsible authority.

**Note:**

The Summary of Safety and Clinical Performance of **Grandio Flow** can be found in the European database on medical devices (EUDAMED – <https://ec.europa.eu/tools/eudamed>). Detailed information can also be found at [www.voco.dental](http://www.voco.dental).

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