**VOCO****Clip/Clip F**

Instructions for use

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Mode d'emploi

**Product description:**

**Clip** is a light-curing material for temporary filling and sealing. The elastic consistency assures an easy removal, even in case of undercuts, and finishing the cavity is not necessary. Thus the fitting of the temporary filling will be improved.

**Clip** is available in syringes and Caps for direct application into the cavity.

**Clip F** contains fluoride and continuously releases fluoride ions.

**Indications:**

All types of temporary fillings, temporary protection in inlay and onlay technique.

**Application:**

Preparation, lining and impression to be carried out as usual. Rinse and dry the cavity.

**Clip syringe:**

Apply **Clip** into the cavity with a suitable instrument (round condenser, spatula, etc.).

**Clip Caps:**

Insert Caps into the opening of the dispenser (follow the instructions for use). Turn Caps into the required direction and remove the protective cap. Apply **Clip** directly into the cavity by slightly and evenly pressing the levers of the dispenser together.

**Further handling:**

Contour with slight excess. For use in the inlay and onlay technique, keep the approximate contact points free of resin with dental floss or a matrix before curing to maintain tooth mobility.

Cure for 20 - 40 s (depending on layer thickness) with a conventional polymerisation device. The light output should be a minimum of 500 mW/cm<sup>2</sup> on halogen polymerisation devices and 300 mW/cm<sup>2</sup> on LED devices. As a rule, 5 mm deep layers are sufficiently cured in 40 s. Increased exposure times are required when the light source cannot be brought close enough to the material.

Check occlusion. Excess can be removed with a scalpel. Silicone points are suitable for grinding.

**Removal of the temporary filling:**

Insert a rigid instrument into the temporary filling and extract it.

**Information, precautionary measures:**

Adhesion to dental resins or linings (e.g. light-curing materials) may occur. In such cases apply a suitable insulating agent (e.g. glycerine, vaseline). Phenolic substances, especially preparations containing eugenol or thymol, inhibit the curing of **Clip/Clip F**. Therefore, do not use zinc oxide eugenol cements or other eugenol-containing substances in combination with **Clip/Clip F**.

In case of more than two weeks' stay in the mouth, and in occlusal or marginal areas, reinforcement with composite is recommended. Apply a small quantity of a composite material onto the corresponding areas and polymerise after contouring together with **Clip**.

The slight expansion of **Clip** ensures a tight seal of the temporary restoration. When treating adjacent teeth simultaneously, remove excess material in approximal areas prior to polymerisation to maintain tooth mobility and contact points.

**Clip** contains acrylates, methacrylates and BHT.

**Clip F** contains methacrylates, BHT and fluorides.

**Clip/Clip F** should not be used in case of known hypersensitivities (allergies) to any of these ingredients. Avoid contact of unpolymerised paste with mucous membranes to avoid sensitivity reactions.

**Storage:**

Store at temperatures of 4 °C to 23 °C. Refrigeration is advised if the material is not used regularly. Do not use after expiry date.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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VC 60 BB1284 E1 0818 99 © by VOCO
VOCO GmbH