(DE) Gebrauchsanweisung	3
(GB) Instructions for use	11
(FR) Mode d'emploi	19
(ES) Instrucciones de uso	28
(IT) Istruzioni per l'uso	37
(PT) Instruções de uso.....	45
(NL) Gebruiksaanwijzing	53
(SE) Bruksanvisning	61
(DK) Brugervejledning	69
(NO) Bruksanvisning	77
(FI) Käyttöohjeet	85

MD

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Medical device

Dispositif médical

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Medisch hulpmiddel

Medicinteknisk produkt

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i | BOND® Universal

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Stand: 2021-11

Medical Device – for use by dental health care professionals only.

Intended purpose

iBOND Universal is a light-curing, self-conditioning all-in-one adhesive for use in adhesive, restorative dentistry.

iBOND Universal was developed for the adhesive bonding and repair of direct filling materials and indirect restorations. Use iBOND Ceramic Primer additionally for the adhesive bonding and repair of glass ceramic.

As a universal adhesive, it may be used as follows depending on the indication:

- SELF-ETCH TECHNIQUE for etching, priming, bonding and desensitisation in a single step.
- SELECTIVE ETCH TECHNIQUE: adhesion to the tooth enamel can be maximised by applying selective enamel etching with phosphoric acid.
- ETCH & RINSE TECHNIQUE (total-etch) with an additional etching step with phosphoric acid on enamel and dentine.

This dental material is suitable for patients requiring dental treatment for the following indications with consideration of the contraindications. For pregnant and nursing women, the treatment risks must be weighed carefully against the benefits taking into consideration the unborn child or infant. Kulzer medical devices ensure the rehabilitation of oral functions such as chewing, speaking and aesthetics. They stabilize the remaining dentition and/or the alveolar ridge restoratively or prosthetically.

Indications

1. Bonding of direct restorations for all cavity classes (Black) using light-curing, dual-curing or self-curing methacrylate based composites / compomers.
2. Bonding of light-curing, dual-curing or self-curing core build-up materials.
3. Sealing of hypersensitive tooth areas.
4. Sealing of cavities prior to amalgam restorations.
5. Bonding of fissure sealants.
6. Sealing of cavities and core preparations prior to temporary cementation of indirect restorations (according to the immediate dentine sealing technique).
7. Cementation of indirect restorations with light-curing, dual-curing or self-curing adhesive resin cements.
8. Intraoral repair of composite and compomer restorations, porcelain fused to metal, all ceramic as well as metal restorations.

Contraindications

The use of this product is contraindicated in case of known or suspected allergies against (meth)acrylate compounds. Do not allow iBOND Universal to come into direct contact with exposed pulp tissue.

Potential side effects

The use of this product is contraindicated in case of known or suspected allergies against (meth)acrylate compounds.

Composition

Aceton, UDMA, Water, 4-META, MDP

Application steps

Preparation of the tooth surface (except fissure sealants)

- Prepare cavity according to the rules of the adhesive technique.
- Clean with oil-free and fluoride-free paste.
- Clean with water spray and dry with water-free and oil-free air.
- Apply pulp protection and/or base selectively.

Application options

iBOND Universal may be used as a self-etching adhesive in the self-etch technique.

For selective enamel etching, apply 35 % phosphoric acid etching gel (e.g. iBOND Etch 35 Gel) to the prepared enamel and allow to work for 20 – 30 seconds. Then rinse thoroughly with water and dry with an oil-free air flow. Do not overdry.

In the etch & rinse technique, apply 35 % phosphoric acid etching gel (e.g. iBOND Etch 35 Gel) to the prepared enamel and dentine and allow to work (20 – 30 seconds on enamel/15 seconds on dentine). Then rinse thoroughly with water and dry with an oil-free air flow. Do not overdry.

- Fissure sealing: after pretreating with polishing paste or powder jets and water spray, etch with 35 % phosphoric acid etching gel (e.g. iBOND Etch 35 Gel) for 20 – 30 seconds, rinse thoroughly, dry and apply adhesive.

Caution: Pay attention to the directions for using the etching gel.

Preparation of iBOND Universal

From the bottle:

- Drop iBOND Universal into the mixing well and use immediately (within 3 minutes). When product has been removed, reseal carefully.

From the single dose:

- Open the single dose just before use.
- When removing the liquid, briefly rotate the applicator tip or brush in the single dose so that it is completely wetted.

Application to tooth surface

- Use the applicator tip or brush to apply iBOND Universal to the entire surface of the cavity.
- Gently rub in the adhesive for 20 seconds.
- Air dry carefully with a gentle oil-free air flow until the adhesive film no longer moves. Move the air flow from outside to inside with increasing intensity at the same time. Caution: If the air flow at the start of blowing is too strong, this leads to the bonding agent thinning out and may result in poor adhesion and post-operative discomfort. Do not leave any excess bonding material/puddles, e.g. on the cavity base.
- The surface must be visibly glossy both after applying. If the cavity surface does not appear universally glossy, apply iBOND Universal again as described.
- Polymerise iBOND Universal for 10 seconds with a conventional curing lamp (wavelength range 440 – 480 nm; light intensity min. 600 mW/cm²).
- Introduce restoration material immediately in accordance with the manufacturer's instructions.
- Caution: If the light power is too low this leads to insufficient adhesion. Curing lamps should be inspected at regular intervals with calibrated test equipment. The light emission window should be placed as close as possible to the cavity surface during polymerisation.
- Caution: Carefully remove oxygen inhibited layer with an alcohol soaked pellet after light curing.

Indirect restorations

Preparation of the material surfaces

Composite and metal restorations in addition to zirconium and aluminium oxide ceramic:

Prepare the restoration in accordance with the manufacturer's instructions. If there is no information available, clean thoroughly with water spray and dry with oil-free air flow.

Caution: Do not clean aluminium and zirconium oxide ceramic and metal surfaces with phosphoric acid.

Silicate/glass ceramic and lithium disilicates:

Ceramic materials require different surface treatments; ask the manufacturer or laboratory about the correct procedure. Always etch Silicate/glass (incl. glass ceramic) with hydrofluoric acid in accordance with the manufacturer's instructions.

Always apply iBOND Ceramic Primer after pre-treating the ceramic surfaces.

- Measure out the desired quantity of drops of iBOND Ceramic Primer into the mixing well.
- Use the applicator tip or brush to apply to the whole surface of the material without delay and allow to evaporate for 20 seconds.
- Then dry briefly with an oil-free air flow

Prepare and apply iBOND Universal immediately as previously described in "Application to tooth surface". Ensure complete coverage of the entire restoration surface. At the same time take care to ensure that no puddles of adhesive remain before light curing to ensure an optimal fit.

After curing of iBOND Universal (see above), apply the adhesive resin cement to the restoration in accordance with the manufacturer's instructions, place on the tooth surface prepared with adhesive and cure. In the case of dual-curing cements, additional light-curing of the cement is recommended after removal of the remaining cement.

Intraoral repair

Preparation of the material surfaces

Composite, compomer and metal restorations in addition to zirconium and aluminium oxide ceramic:

Roughen the restorations to be repaired. Rinse thoroughly and dry with an oil-free air flow.

Silicate / glass ceramic and lithium disilicate:

Roughen the surface to be repaired with an intraoral powder jet device (aluminum oxide, 50 µm, 1.5 bar) or a diamond approx. 1 mm beyond the margin of the fracture site. Rinse and dry with oil-free air. Apply iBOND Ceramic Primer as previously described under „Indirect Restorations“.

Use iBOND Universal as described in “Application to tooth surface”. After curing of the adhesive, process appropriate restoration material in accordance with the manufacturer’s instructions.

General information

Please pay attention to the following safety instructions and all processing instructions in other sections. This product must be used according to the instruction for use only. Any use not complying with these instructions is at the discretion and sole responsibility of the dental professional itself.

- Visually inspect for damage before use. Damaged products must not be used.
- Materials containing eugenol may impair the polymerisation of iBOND Universal.
- iBOND Universal is light curing and should be protected against ambient light.
- Do not treat the surface of zirconium and aluminium oxide ceramic and metal with phosphoric acid.
- The use of disinfectants, astringents and rinsing solutions containing EDTA or hydrogen peroxide is not recommended. Residues of these materials may impair the adhesion and curing reaction of the adhesive. But GLUMA Desensitizer and chlorhexidine solutions may be used.

Warnings and safety instructions

Highly flammable liquid and vapour. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources - no smoking. May cause skin and eye irritation. Avoid contact with skin. Skin contact may cause sensitization. Wear suitable protective gloves. In case of contact with skin, immediately wash with plenty of water and soap. Avoid contact with eyes. In case of contact with eyes, rinse with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If symptoms persist, consult an eye specialist. May cause respiratory irritation. Avoid breathing dust/fume/gas/mist/vapours/spray. May cause drowsiness or dizziness. Harmful to aquatic life with long lasting effects. The adhesive must not be swallowed by the patient, otherwise incompatibilities may arise.

- Do not apply adhesive directly to the exposed pulp.
- Eye protection is recommended for the patient. The use of a rubber dam is strongly recommended.
- Seal iBOND Universal and iBOND Ceramic Primer bottles carefully after each use.
- Packages marked with the symbol ☒ are intended for single use. Dispose of after using once and do not use on other patients.
- Do not reuse bottles after contamination with body fluids in order to prevent cross-contamination.
- Saliva or blood in the field of application during treatment may lead to failure. Use rubber dam (e.g. Ivory) or other suitable means of isolation.

Special advice

For SSCP please see EUDAMED database (<https://ec.europa.eu/tools/eudamed>) or contact sscp-service@kulzer-dental.com

To find the SSCP in the EUDAMED use the UDI code attached at the packaging of the product.

Storage conditions

Do not expose to direct sunlight. Store in a sufficiently well ventilated dry place. Do not store above room temperature 25°C (77°F). Incorrect storage may shorten the shelf-life and the product may prematurely lose its effectiveness. Do not use after the expiry date

Disposal information

Recommendation: Dispose of in accordance with official regulations. Do not dispose of contents or partially emptied packaging in the household waste or allow them to enter the sewage system. European Waste Catalogue: 180106 Chemicals consisting of or containing hazardous substances.

Complaint notice

Please quote batch number and article number in all correspondence about the product.

Article number **REF** and batch number **LOT**: See note on syringe/packaging.

According to EU Medical Device Regulation users/ patients are obliged to report serious incidents with a medical device to the manufacturer and to the competent authority of the country, the incident occurred.

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Kulzer, LLC