Instructions for use

- EN more information on www.fkg.ch
- FI Lisä tietoja www.fkg.ch
- DE weitere Informationen auf www.fkg.ch

D-Race

Endodontic instruments
Endodontiainstrumentit
Endodontische Instrumente

FKG REF. - 99.7AA.10.13A.XX - notice n°102 - 2015/06

D-Race

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EN Golden rules

- > Let the instrument work
- > Work 3-4 seconds in a row, get out, clean the blade and irrigate the canal

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D-Race 1 blister includes/sisältää/enthält: 2 × DR1 + DR2

«Desobturation» prior to retreatment, because access must be made, the entrance of the canal «desobtured» to reach inside the canal itself so that it could be properly retreated. Retreatment must be achieved because the first RCT failed or needs improvement.

D-Race will desobturate most canals, whatever the filling used – gutta-percha, carriers, paste, resin-based materials. D-Race removes a large part of the ancient filling material. With only 2 instruments time is saved and safety of treatment is increased.

Once desobtured, the canal will be retreated.

General advice to take into consideration for removal of filling material

1. Take 2 x-rays from different angles.
2. A precise appreciation of canal anatomy is very important.
3. Find entrance of canal and access the obturation material.
4. Prepare, if needed, a cavity 1-2 mm with DR1 and drop solvent to soften obturation material. If necessary, use a heat plugger or ultrasonics.

Coronal third preparation

1. Gently engage DR1 rotating at 1000 rpm (recom. torque 1.5 Ncm) into the obturation material.

The active tip of DR1 facilitates the initial penetration.

Median and Apical thirds preparation

1. Insert DR2 rotating at 600 rpm (1 Ncm).
2. Don't force the way, clean and check the blade regularly.
3. Continue to advance apically as long as filling material is visible on the flute.
4. If needed, drop solvent to help removal of obturation material.
5. Establish final WL and finalize shaping with standard NiTi Race files.

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EN

Sterilization protocol for reusable medical devices

Non sterile products

First use

Carefully follow the rules herebelow
Cleaning and Sterilization

Sterile products

First use

1. Open the alveolus
2. Extract the instrument (use of gloves imperative)
3. Use the instrument according to good dental practice

Second use:

Carefully follow the rules herebelow
Second and following uses.

Cleaning

1. Manual cleaning with/without Ultrasonic assistance
 2. Through cleaning allows efficient disinfection & sterilization.
 2. Rinse
- Under distilled/demineralized water at least 1 min.
Then dry the instruments.

3. Inspection - Check
- Sort out damaged or worn out instruments.
4. Packing
- Place the instruments in a suitable device/container and pack the devices in sterilisation pouches ISO 11607-1.

Sterilization

5. According EN/ISO 17664 protocol
- Autoclave: 134°C/273°F, 2.2 bar during 18 min.

6. Storage
- Keep devices in sterilization wrap/pouch in a dry and clean environment.

 Always refer to manufacturer's instructions for cycle & duration

Second and following uses

7. Pre-disinfection
- Plunge instruments after use in detergent solution and eventually brush them manually

8. Rinse
- Under running water at least 1 min. Then dry the instruments.
9. Follow steps 1 - 6 above

 Complete operating protocol and warnings on www.fkg.ch

General information

Pre-/Disinfection DO NOT USE

- Solutions containing:
- phenol (corrosion)
 - aldehyde (blood fixation)
 - di-/triethanolamines (corrosion)

Sterilization with chemiclav or dry air devices have not been validated by FKG Dentaire. Refer to manufacturer's instructions for cycle & duration. Nevertheless our instruments withstand such methods.

Caution:

This product contains nickel and should not be used for individuals with known allergic sensitivity to this metal.

A single patient use is recommended to avoid cross-contamination.

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