VOCO

Ceramic Bond

GB**CE** 0482

Instructions for use

Product description:

Ceramic Bond is a one-component silanisation agent for durable adhesive bonding between metal-free indirect restorations (also zirconium dioxide) and luting composites, as well as for the pretreatment of glass fibre reinforced composite posts.

Indication:

Ceramic Bond is used as a coupling agent to achieve an optimal bond strength between indirect restorations made of either silicate ceramics, aluminium oxide, zirconium dioxide or composite, as well as glass fibre-reinforced composite posts, on the one side and methacrylate-based luting composites on the other side.

Application:

Luting of indirect restorations:

Silanising the surface of indirect restorations, including those made from zirconium dioxide, creates optimum conditions for longterm bonding of restoration and luting composite.

1. Once the fit of the restoration has been checked, etch or sandblast its bonding surfaces to clean them, rinse carefully with water and dry with oil-free air.

2. When using the bottle:

Dispense **Ceramic Bond** onto a **Single Tim** or a disposable brush. Alternatively, dispense **Ceramic Bond** onto a mixing pallet.

When using the SingleDose:

Detach a **SingleDose** at the perforated line and turn the printed side upwards. Pierce the foil with a **Single Tim**, enlarge the opening and moisten the brush with liquid through circular motion. Apply **Ceramic Bond** onto the conditioned surfaces with a **Single Tim** or disposable brush, allow to take effect for 60 s and dry carefully with oil-free air.

3. Apply the luting composite as described in the respective instructions for use and place the restoration.

Luting of glassfibre-reinforced posts:

1. After preparation of the root canal, try the fit of the glassfibre-reinforced post, shorten the post extraorally, then clean with alcohol and dry with oil-free air.

2. Coat the entire post with **Ceramic Bond** (for details on dispensing from bottle or **SingleDose**: see above) and allow to dry for 60 s.

3. Cement the post in the root canal according to the instructions of the manufacturer.

Information, precautionary measures:

Ceramic Bond contains organic acid, 3-methacryloxypropyltrimethoxysilane and acetone. **Ceramic Bond** should not be used in case of known hypersensitivities (allergies) to any of these ingredients.

Contamination of the silanised bonding surface (e.g. saliva or blood) may reduce the bonding strength of **Ceramic Bond**. Contaminated surfaces must be cleaned again with etchant, rinsed and dried. Close the bottle after use.

Do not allow liquid **Ceramic Bond** to contact the skin/mucosa or eyes as this may have a slightly irritating effect.

Storage:

Store **Ceramic Bond** at temperatures of 4 °C to 23 °C.

Refrigeration is recommended.

Do not use after expiry date.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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VOCO GmbH