

## Registrado® X-tra

Addition-curing silicone for bite registration, extra hard, fast setting



### Indications

Occlusal bite registration  
 Fixing of face bow registrations  
 Location material for intraoral pin tracing templates

### Advantages

- Optimised final hardness
- Allows application without time pressure
- Highly thixotropic qualities
- Good grinding properties
- Very precise impressions
- With economic mixing tips

### Application



Excellent surface affinity



No resistance on occlusion



The bite registration can be easily trimmed with a scalpel



**SC**  
self-curing

**auto mix**

**Shore-A 90**

REF 2035	Cartridge 2 x 50 ml, accessories
REF 2036	Cartridge 8 x 50 ml, accessories
REF 2158	Dispenser – type 2
REF 2196	Mixing tips type 16, 50 pcs.
REF 2140	Intraoral tips type 2, 50 pcs.

## Fit Test® C & B

Addition-curing silicone material for checking of prosthetic reconstructions

### Indications

Control of precision of fit, marginal adaptation and marginal gaps of

- full or partial crowns as well as bridges
- inlays / onlays
- cast post abutments / post crowns
- metal denture bases
- implant-supported crowns and bridges
- attachments and other secondary parts of dentures

### Advantages

- Fast diagnostic help for all restorative pieces
- QuickMix syringe for simple handling – automatic bubble-free mixing
- Direct application
- Optimal colour for metal and ceramic restorations
- Fast results due to short working and setting time

**auto mix**

REF 2095	QuickMix syringe 5 ml, accessories
REF 2261	Dispenser – QuickMix
REF 2188	Mixing tips type 9, 50 pcs.
REF 2191	Mixing tips type 11, only with intraoral tips type 1, 50 pcs.
REF 2139	Intraoral tips type 1, 50 pcs.



CAD / CAM material

Luting

Core build-up

Temporaries

Denture relining

Impression

Repair