

V-Posil® Putty Fast · Putty Soft Fast · Heavy Soft Fast · Mono Fast

Precision impression material, VPS

Indications

V-Posil Putty Fast / Soft Fast is used as preliminary material for

- two-step putty-wash impression technique
- one-step putty-wash impression technique
- two-step putty-wash impression technique using a foil (plastic putty spacer)
- one-step putty impression technique for forming functional peripheries

V-Posil Heavy Soft Fast is to be used as heavy bodied material for

- one-step impression technique (simultaneous technique) using dual viscosities
- two-step impression technique using dual viscosities
- functional impressions

V-Posil Mono Fast is to be used as a medium bodied tray or syringeable impression material for

- taking impressions over fixed / removable restorations and implants (i.e., transferring impression posts and bridge components)
- functional impressions
- fabricating crown and bridgework or inlays
- fabricating full or partial dentures
- relining impressions
- transferring root posts when fabricating posts and cores indirectly
- multi tray technique
- use in the simultaneous mixing technique as well as the putty-wash technique

Advantages

- High precision – strongly hydrophilic addition-curing silicone for good wetting of oral structures
- Treatment comfort – long working time combined with a short intraoral setting time*
- Safe removal – high toughness allows removal without damaging filigree structures
- Perfect prosthetics – good hydrophilicity and high recovery after deformation of the set impression makes the pouring out and thus the restoration perfect

* V-Posil Putty Fast 380 ml: Extraoral working time ≤ 01:30 min., time in mouth ≥ 02:30 min.

Monophase impression



Double mixing impression



Correction impression



V-Posil Putty Fast

- REF 2561 Jar 2 × 450 ml, accessories
- REF 2562 Jar 8 × 450 ml, accessories
- REF 2564 Cartridge 2 × 380 ml
- REF 2204 Mixing tips type 31, 50 pcs.

V-Posil Putty Soft Fast

- REF 2579 Jar 2 × 450 ml, accessories

V-Posil Heavy Soft Fast

- REF 2565 Cartridge 380 ml, accessories
- REF 2566 Cartridge 2 × 380 ml
- REF 2567 Cartridge 2 × 50 ml, accessories
- REF 2568 Cartridge 10 × 50 ml
- REF 2196 Mixing tips type 16, 50 pcs.
- REF 2204 Mixing tips type 31, 50 pcs.

V-Posil Mono Fast

- REF 2569 Cartridge 380 ml, accessories
- REF 2570 Cartridge 2 × 380 ml
- REF 2571 Cartridge 2 × 50 ml, accessories
- REF 2572 Cartridge 10 × 50 ml
- REF 2196 Mixing tips type 16, 50 pcs.
- REF 2204 Mixing tips type 31, 50 pcs.
- REF 2133 Intraoral tips type 6, 50 pcs.

CAD / CAM material

Luting

Core build-up

Temporaries

Denture relining

Impression

Repair