

Cotton pellets

Cotton Pellets are cotton fiber spun into round ball, used to dry/absorb or apply substances during dental procedures.

PD Cotton pellets are highly absorbent and uniform in shape. The cotton used is a purified medical grade cotton with TCF (Total Chlorine Free) certification. The purification process is done with "oxygen purification" (through hydrogen peroxide), leading to a cotton **without any trace of chlorine**.

The product is available in different diameters. The cotton pellets are packaged either in bulk (18g or 32g) or in refill (6x3g).

*EU Class I medical device.

** Product availability may depend on your local regulation and product registration status.



		REF
N°1 Refill (model A), 6x3 g (~630 pcs)	Ø 8 mm	25908
N°2 Refill (model A), 6x3 g (~1470 pcs)	Ø 6 mm	25901
N°3 Refill (model A), 6x3 g (~2220 pcs)	Ø 4.3 mm	25902
N°4 Refill (model A), 6x3 g (~4050 pcs)	Ø 3.6 mm	25903
N°1 Bulk (model B), 32 g	Ø 8 mm	25904
N°2 Bulk (model B), 32 g	Ø 6 mm	25905
N°3 Bulk (model B), 18 g	Ø 4.3 mm	25906
N°4 Bulk (model B), 18 g	Ø 3.6 mm	25907

