



Chisels

Chisels are designed to dislodge bridges of calculus, using a push stroke, on anterior mandibular teeth.



Cushing #61 Chisel, Single-End

#24

SCU61



Cushing #62 Chisel, Single-End

#24

SCU62



Darby-Perry 5K/6K Chisel

#2

SDP5/6K



Darby-Perry 9K/11K Chisel

#2

SDP9/11K

DiamondTec™ File Scalers

DiamondTec™ File Scalers offer 360° of diamond coated coverage. The absence of blade edges allows use of all sides of the working end. Ideal for furcations, root depressions, mesial/distal line angles and developmental grooves.



DiamondTec™ Mesial Distal File Scaler

#2

SDCM/D7



DiamondTec™ Nabors File Scaler

#2

SDCN7