

# Konservierende Behandlung







## Restorative Dentistry

### Composite Instrumente

Composite instruments

mit blauer TITAN-OXID Beschichtung zur Verringerung der Klebewirkung des Composites

blue titanium oxide coating to reduce the sticking of the composite



### Composite-SET

#### SI-1051-SET6

Inhalt/Content: je/each 1x SI-1051/14-HBL,  
SI-1051/95-GE, SI-1051/104-OR,  
SI-1051/108-GR, SI-1051SF/1-BL,  
SI-1051/Carve-PK

