



Euronda[®]
Eyewear

Euronda Eyewear

Protection at first sight

Euronda Eyewear offers a wide range of protective glasses and visors for operators and patients, designed and manufactured to ensure maximum safety in every situation while guaranteeing comfort and an appealing design.

Versatile solutions for every need: can also be used over corrective glasses, available with orange or smoked lenses, with anti-fog and anti-scratch treatments, and protection against impact, UV400 rays and blue light (up to 525 nm).



Classification

Personal Protective Equipment Category 2, compliant with EU Regulation 2016/425



Regulations

Compliant with ISO 16321: Contemporary, Cube, Cube Orange, FitUp, Clear, Light Orange, Zero, Cube Smoke

Compliant with EN 166 CE: Evolution, Evolution Orange, Baby Orange, Small Visor, Operator Visor



Variants

Available in different models and colours, with clear, orange, or smoked lenses



Packaging

Individual pack
Box of 10 pieces (except Contemporary sold individually)

Treatments

Anti-fog

When used in combination with surgical masks, lenses may fog up, impairing visibility. Euronda Eyewear protective glasses offer different types of treatment to prevent the formation of condensation:

- **basic:** complies with the requirements of current legislation and prevents the formation of condensation for up to 8 seconds.
- **plus:** superior performance, prevents the formation of condensation for up to 40 seconds.
- **ultra:** advanced protection, prevents the formation of condensation for over 40 seconds.

Anti-scratch

Scratches can compromise lens transparency and thus visibility, exposing the operator to various risks. Euronda Eyewear protective glasses offer different levels of anti-scratch treatment:

- **basic:** compliant with current regulations, this treatment guarantees scratch resistance in accordance with the required standard tests.
- **plus:** ensures a longer lifespan and consistent clarity, with a reduced luminance factor of $\geq 2 \text{ Cd/M}^2 \times \text{lux}$.
- **ultra:** offers advanced performance even in challenging conditions, with high scratch resistance and a reduced luminance factor of $1.2 \text{ Cd/M}^2 \times \text{lux}$.

Lenses

The protective glasses in the Euronda Eyewear range are available with different lens types:

- **clear:** UV protection (UV400).
- **orange:** UV400 protection and blue light protection (up to 525 nm); improved contrast and reduced eye fatigue.
- **smoked:** UV400 protection; blocks light emitted by the dental chair, enhancing visual comfort.

Frame technologies

Euronda Eyewear protective glasses are designed to ensure comfort and ergonomics through advanced technical solutions:

- adjustable nose pad;
- integrated side protectors;
- tilting movement of the temples for a perfect fit and stability;
- adjustable tilt and length of the temples for a perfect fit.

	Anti-fog	Anti-scratch	Lenses	Technologies	Weight
Contemporary	Plus	Plus	Clear	Integrated side protectors	25 g
Cube	Plus	Plus	Clear and orange	Integrated side protectors Tilting movement Length-adjustable arms	37 g
Evolution	Ultra	Ultra	Clear and orange	Adjustable nose pad Tilting movement Length-adjustable arms	37 g
Fit Up	Plus	Plus	Clear	Adjustable nose pad Adjustable length and tilt of arms	25 g
Clear	Basic	Basic	Clear	Integrated side protectors Length-adjustable arms	30 g
Light Orange	Basic	Basic	Orange	Integrated side protectors Adjustable length and tilt of arms	43 g
Zero	Basic	Basic	Clear		19 g
Cube Smoke	/	Basic	Smoked	Integrated side protectors Tilting movement Length-adjustable arms	37 g
Baby Orange	/	Basic	Orange		23 g
Small Visor	/	Basic	Clear		9 g
Operator visor	/	Basic	Clear		25 g



Contemporary



Lightweight and incredibly robust protective glasses with an elegant and sophisticated design, ideal for doctors and operators.

The transparent lenses feature an anti-fog treatment that makes them twice as effective against condensation, and an anti-scratch treatment that provides five times the standard abrasion resistance.

The glossy all-black frame also allows the glasses to be worn over prescription glasses, while the transparent side guards ensure safety and maximum visibility.



Colours

Black and transparent



Weight

25 g



Treatments

Anti-scratch and
Anti-fog plus

Cube

Cube Orange

Protective glasses with a modern, square, wraparound design, engineered to ensure comfort and stability even during prolonged use.

The temples, which are adjustable in length and tilt, and the adaptable nose bridge ensure a personalized, ergonomic fit.

Available in two versions:

- Cube: Clear polycarbonate lenses for wide, clear vision;
- Cube Orange: Orange lenses that protect against impact, UV400 rays and UV525 blue light up to 525 nm.



Colours

Cube: white and light blue
Cube Orange: black and orange



Weight

37 g



Treatments

Anti-scratch
and Anti-fog plus

Evolution

Evolution Orange

Protective glasses with a lightweight and modern design, created to ensure comfort and an optimal fit.

Arms made of overmoulded material maintain the desired shape, while the rubber non-slip nose bridges ensure stability during use.

Available in two versions:

- Evolution: Clear lenses made of high-strength polycarbonate;
- Evolution Orange: Orange lenses that protect against impact, UV400 rays and UV525 blue light up to 525 nm.



Colours

Evolution: white and grey
Evolution Orange: black and orange



Weight

37 g



Treatments

Anti-scratch
and Anti-fog ultra

FitUp

Protective glasses with a wraparound design, developed to ensure maximum protection and visibility while combining comfort and ergonomics.

The temples, which are adjustable in length and tilt, and the soft, non-slip nose bridge ensure a stable and comfortable fit.

High-strength lenses with anti-fog and anti-scratch treatments for reliable performance over time. Available in two colour variants with a transparent frame.



Colours

Transparent and pink and transparent and green



Weight

25 g



Treatments

Anti-scratch and Anti-fog plus

Light Orange



Protective glasses for both operators and patients, featuring side protection and a single lens made of high-impact polycarbonate.

The temples, with adjustable length and tilt, ensure high comfort and a secure fit.

The lenses, which feature anti-fog and anti-scratch coatings, protect against impacts, UV400 rays, and blue light (up to 525 nm).



Colours

Black and orange



Weight

43 g



Treatments

Anti-scratch and Anti-fog basic

Clear



Protective glasses for operators and patients in a contemporary design.

Made entirely of polycarbonate, they offer the perfect balance between lightness and resilience.

The temples, which are adjustable in length and tilt, ensure an optimal fit and a high level of comfort during use.



Colours

Transparent and black



Weight

30 g



Treatments

Anti-scratch
and Anti-fog basic

Zero



Protective glasses with a minimalist design – simple and ultra-lightweight.

Made entirely of polycarbonate, they are resilient and easy to handle, ensuring comfort during use.

The slim-profile temples optimize comfort. They are the lightest model in the whole range.



Colour

Transparent



Weight

19 g



Treatments

Anti-scratch
and Anti-fog basic

Cube Smoke



Protective glasses with a square wraparound design, engineered to ensure comfort and stability even during prolonged use.

The adjustable temples ensure a customized fit, while the upper front section provides a shock-absorbing effect and enhanced protection.

The UV400-resistant, smoked polycarbonate lenses shield against light from the dental chair, promoting visual comfort.



Colours

Black and green



Weight

37 g



Treatments

Anti-scratch basic

Baby Orange



Lightweight, easy-to-handle protective glasses for children, designed to ensure comfort and safety.

Made of polycarbonate with black arms and orange lenses that protect against impact, UV400, and blue light up to 525 nm.



Colours

Black and orange



Weight

23 g



Treatments

Anti-scratch basic

Small visor



Short disposable operator visor designed to protect the eyes and cheekbones.

Ultralight and comfortable, equipped with interchangeable clear shields and a frame available in a variety of colours.



Colours

Black, blue, green, pink, yellow



Weight

9 g



Treatments

Anti-fog basic



Packaging

5 interchangeable coloured frames and 20 visors

Operator visor



Ultralight disposable visor for operators, consisting of a clear frame and interchangeable shield.

Designed to provide complete protection from splashes and potentially dangerous fluids, they can also be worn over corrective glasses and magnifying systems (loupes).

Available with replacement lenses.



Colour

Transparent



Weight

25 g



Treatments

Anti-scratch basic



Packaging

1 frame and 10 visors

