

OptiView[™]

Lip and Cheek Retractor





Enjoy the unobstructed view.

Access to the working area is among the first steps of dental procedures. Ranging from surgical implant placement to cosmetic teeth whitening, the use of a lip and cheek retraction system is a great way to help the patient keep his or her mouth open while you're working. Dental cheek retractors lock in place to hold the patient's mouth open throughout the time you are working. Most conventional retractors provide either lip or cheek retraction and are uncomfortable to wear due to their high structural tension.




OptiView is an innovative lip and cheek retractor that **provides retraction of lips and cheeks** for improved buccal and gingival access **at the same time**. Thanks to its anatomic shape, OptiView minimizes tension, offering patient comfort even during longer dental treatments.



OptiView™


Lip and Cheek Retractor

- **Dual-purpose:** permits both lip and cheek retraction with one device, without additional hands or instruments
- **Uniform perioral retraction:** maximizes visibility and access to the working field
- **Passive placement:** stable, flexible 3D design, ensures quick and gentle placement and removal, facilitating the placement of the saliva ejector
- **Cushioned lip-supports:** enhance patient comfort during longer procedures
- **Latex free and autoclavable** at 134° C for 3 minutes
- **Available in two sizes:** standard (white) and small (light blue)




Standard

104 mm



Small

97 mm



Indications

Prophylaxis, diagnostic, orthodontic, bleaching procedures, therapeutic dental applications. In general for children, elderly and handicapped patients. For those reacting with panic and claustrophobic symptoms under rubber dams; for those suffering from strong nausea or for patients with strong lip pressure. Convenience offered while taking pictures.



OptiView in orthodontic procedures
© Dr. Jessy Askar, France

OptiView - Kit

CONTENT: 2 Lip-and Cheek Retractors with mounted cushions

Part No.	Description
5500	OptiView Standard Kit
5502	OptiView Small Kit

OptiView - Refill

CONTENT: 10 refill cushions

Part No.	Description
5501	OptiView Refill Cushions

Kerr™
PREVENTION

prevention.KerrDental.eu

KerrHawe SA - P.O. BOX 272, 6934 Bioggio, Switzerland

Phone +41 91 610 05 05 - Fax +41 91 610 05 14 - info@kerrdental.eu