



Scaler and Curette Sets

We offer a wide variety of Scaler and Curette Sets including McCall, Goldman-Fox, and Cattoni. All scalers and curettes included in these sets are available in multiple handle options.

McCall Scalers and Curettes

McCall Scalers and Curettes are posterior instruments featuring dual cutting edges and available with either a rounded or pointed toe.




McCall 11/12 Scaler	
#6	SM11/126
#7	SM11/127
#8, Gray	SM11/128
#2	SM11/12




McCall 11/12A Scaler	
#6	SM11/12A6
#2	SM11/12A



McCall 11/12 Thin Scaler	
#6	S11/12S6



McCall 13/14 Curette	
Harmony	SM13/14SXE2
#9	SM13/149E2
#G9, Steel Grip	SM13/14G9E2
#C8, Lavender	SM13/14C8E2
#6	SM13/146
#7	SM13/147
#8, Gray	SM13/148
#2	SM13/14



McCall 13/14 Pointed Curette	
Harmony	SM13/14SXE2
#9	SM13/14S9E2
#G9, Steel Grip	SM13/14SG9E2
#C8, Lavender	SM13/14SC8E2
#6	SM13/14S6
#8	SM13/14S8
#7	SM13/14S7
#2	SM13/14S