

K-Flex™ Files

Stainless steel files

K-Flex™ Files feature a unique rhomboid design, which provides a short axis for flexibility while maintaining cross-sectional strength. K-Flex is more flexible than conventional “K” files of comparable size, particularly in sizes #30 and above, which facilitates the negotiation of curved canals and their preparation with minimal potential for ledging and perforating.



K-Files™

Stainless steel files

K-Files can bypass obstructions with ease. They are designed to provide a smooth tactile feel during instrumentation.



K-Flex files — Pack of 6

Size	21 mm	25 mm	30 mm
80		15174	15258
70		15170	15254
60	15334	15166	15250
55	15330	15162	15246
50	15326	15158	15242
45	15322	15154	15238
40	15318	15150	15234
35	15314	15146	15230
30	15310	15142	15226
25	15306	15138	15222
20	15302	15134	15218
15	15298	15130	15214
10	15294	15126	15210
08	15290	15122	15206
06	60120	60144	60193

K-Files — Pack of 6

Size	21 mm	25 mm	30 mm
80		06099	06170
70	06024	06085	06166
60	06023	06077	06163
55	06022	06076	06162
50	06021	06075	06161
45	06020	06074	06160
40	06019	06069	06159
35	06016	06067	06157
30	06015	06066	06156
25	06014	06065	06155
20	06005	06064	06153
15	06004	06063	06152
10	06003	06062	06151
08	14127	14128	14129
06	62249	62250	62251

K-Files — Assorted files — Pack of 6

Size	21 mm	25 mm	30 mm
15 — 40	821-8021	821-8025	821-8030
45 — 70	61363		
45 — 80	821-8121	821-8125	821-8130