


sterilisation pouches

Sterilisation pouches provide an easy sealing system that eliminates the need for expensive heat-sealing equipment. This ensures a uniform seal that withstands the sterilisation process and keeps instruments sterile.




6000

57 mm x 90 mm

(200 pcs/box)

*diamond and carbide burs,
scaler tips, orthodontic
appliances etc.*




6001

57 mm x 130 mm

(200 pcs/box)

*mirror heads, matrix retainers,
orthodontic appliances,
rubber dam clamps etc.*




6003

90 mm x 145 mm


(200 pcs/box)

*slow speed handpieces,
scissors, small
instruments etc.*




6005

70 mm x 230 mm
(200 pcs/box)
*probes, currettes, scalers,
tweezers and other hand
instruments etc.*




6007

90 mm x 230 mm
(200 pcs/box)
*forceps, elevators, scissors and
other hand instruments etc.*




6009

135 mm x 260 mm
(200 pcs/box)
*ruber dam forceps, rubber dam
frames and other large hand
instruments etc.*



6011


190 mm x 330 mm
(200 pcs/box)
*large instruments and
sterilisation cassettes*



6015

250 mm x 380 mm
(100 pcs/box)
large instruments and
sterilisation cassettes

- Manufactured from 70 gsm medical grade paper and 60 µm blue tinted film to help detect puncture holes.
- Easy to read steam indicators change from blue to dark grey to ensure sterile conditions have been met.
- Triple adhesive seal to protect contents.
- Internal and external indicators.
- Shelf life: 3 years before use, 6 months after processing.
- Each pouch is printed with LOT no. for traceability.
- Large range of sizes for all applications.
- Complies with ISO 11607-1 and EN 868-5.





6000
10 x boxes per carton



6001
10 x boxes per carton



6003
10 x boxes per carton



6005
10 x boxes per carton



6007
10 x boxes per carton



6009
10 x boxes per carton



6011
10 x boxes per carton



6015
10 x boxes per carton

Order	Product Description	Packaging	Carton QTY	Carton Dimensions L x W x H	Carton Weight	Cartons per EU Pallet
6000	6000 Pouches 200s 57 mm x 90 mm	200 pcs/box	10 boxes/carton	33 cm x 15 cm x 15 cm	2.35 kg	208 cartons
6001	6001 Pouches 200s 57 mm x 130 mm	200 pcs/box	10 boxes/carton	35 cm x 19 cm x 15 cm	3.05 kg	156 cartons
6003	6003 Pouches 200s 90 mm x 145 mm	200 pcs/box	10 boxes/carton	35 cm x 21 cm x 21 cm	4.95 kg	90 cartons
6005	6005 Pouches 200s 70 mm x 230 mm	200 pcs/box	10 boxes/carton	35 cm x 29 cm x 16 cm	5.80 kg	96 cartons
6007	6007 Pouches 200s 90 mm x 230 mm	200 pcs/box	10 boxes/carton	35 cm x 29 cm x 21 cm	7.30 kg	72 cartons
6009	6009 Pouches 200s 135 mm x 260 mm	200 pcs/box	10 boxes/carton	33 cm x 32 cm x 30 cm	11.70 kg	36 cartons
6011	6011 Pouches 200s 190 mm x 330 mm	200 pcs/box	10 boxes/carton	42 cm x 38 cm x 35 cm	20.60 kg	20 cartons
6015	6015 Pouches 100s 250 mm x 380 mm	100 pcs/box	10 boxes/carton	44 cm x 34 cm x 28 cm	15.75 kg	28 cartons