



## Thank you and goodbye iBOND Self Etch!

Thank you iBOND Self Etch for the great time during all these years. You can now take your well-deserved retirement. Your successor **iBOND Universal** is already in the starting blocks and is ready to take over!

Giving a hand to oral health.



**KULZER**  
MITSUI CHEMICALS GROUP

# One bottle, all benefits!

With iBOND Universal you can use your application routine from iBOND Self-Etch, but simultaneously benefit from all advantages universal adhesives have:


- Only one adhesive system needed for all bonding situations because it can be used also for luting of indirect restorations in combination with adhesive resin cement
- Compatible to all light, dual- and self-curing materials without additional dual-cure activator
- Bonds to all dental material surfaces (for silicate ceramic in combination with iBOND Ceramic Primer)
- Can be used with all etching techniques
- The usage of one adhesive minimizes the risk of application errors
- Dentine bond strength on gold standard level even after simulated ageing

\* Ehlers V, Ernst C-P, Kastrati A, Gerlach M, Eppinger R, Willershausen B: Bond strength comparison of a new universal adhesive to a gold-standard. J Dent Res 94 (Spec Iss B), # 0544 CED, 2015, (<http://www.iadr.org>).

*“After several years of thorough clinical experience and many completed doctoral theses I can state that universal adhesives such as iBOND Universal have reached the standard of conventional multi-step adhesives.”*



iBOND Universal has in each mode even after simulated ageing bond strength on gold standard level



Source: University of Mainz & Kulzer R&D, Germany, 2014



1. Total-etch: With an etching step for enamel and dentine.
2. Selective enamel etch: With an etching step for enamel.
3. Self-etch: No etching step is required.

Request your sample and contact your sales rep!



Art. code	Content
66061411	iBOND Universal Bottle Refill 1 x 4 ml Bottle
66061412	iBOND Universal Bottle Value Pack 3 x 4 ml Bottle
<b>Single Dose</b>	
66061419	iBOND Universal Single Dose Value Pack 100 Single Doses à 0,15 ml
66061416	iBOND Universal Ceramic Primer 1 x 4 ml Bottle