

Hi-Tray Light Dual Arch

Impression trays in plastic

Impression systems / Accessories

The **Hi-Tray Light Dual Arch** are impression trays in plastic for the dual arch technique in fixed prostheses.

Characteristics

- Plastic structure, of various colors
- Disposable
- Support webbing in white plastic
- Available in 5 types:
 - Full arch
 - Anterior
 - Quadrant
 - Posterior
 - Wide body posterior

Advantages

- The support webbing lets silicone through to form a single body between the upper and lower impressions
- The thickness of the support webbing is imperceptible to the patient (avoiding alterations during bite registration)
- Clear plastic lets you view the correct position of the impression tray on the preparation



Hi-Tray Light Dual Arch

Code	Packaging
D5DANT	Hi-Tray Light Dual Arch anterior 30 pcs
D5DFUL	Hi-Tray Light Dual Arch full arch 30 pcs
D5DPLE	Hi-Tray Light Dual Arch posterior sideless 48 pcs
D5DPOS	Hi-Tray Light Dual Arch posterior 48 pcs
D5DQUA	Hi-Tray Light Dual Arch quadrant 40 pcs