


## • APPLICATION •

### Standard

|                                                  |                                                                                   |                 |     |
|--------------------------------------------------|-----------------------------------------------------------------------------------|-----------------|-----|
| <b>Miraject P (Carpule)</b><br>REF 254 206       |  | 0.8 x 25 mm/21G | 21G |
| <b>Miraject P Super (Carpule)</b><br>REF 254 207 |  | 0.9 x 32 mm/20G | 20G |
| <b>Miraject PL Super (Luer)</b><br>REF 254 214   |  | 0.9 x 22 mm/20G | 20G |
| <b>Miraject Endo Luer</b><br>REF 254 222         |  | 0.3 x 40 mm/30G | 30G |
| <b>Miraject Endo Luer</b><br>REF 254 215         |  | 0.5 x 40 mm/25G | 25G |

### Endotec Luer

|                                             |                                                                                     |                 |     |
|---------------------------------------------|-------------------------------------------------------------------------------------|-----------------|-----|
| <b>Miraject Endotec Luer</b><br>REF 254 216 |  | 0.8 x 25 mm/21G | 21G |
| <b>Miraject Endotec Luer</b><br>REF 254 217 |  | 0.6 x 25 mm/23G | 23G |
| <b>Miraject Endotec Luer</b><br>REF 254 218 |  | 0.5 x 25 mm/25G | 25G |
| <b>Miraject Endotec Luer</b><br>REF 254 228 |  | 0.4 x 25 mm/27G | 27G |
| <b>Miraject Endotec Luer</b><br>REF 254 219 |  | 0.3 x 25 mm/30G | 30G |

### Endotec Luer Duo

|                                            |                                                                                       |                 |     |
|--------------------------------------------|---------------------------------------------------------------------------------------|-----------------|-----|
| <b>Miraject Endotec Duo</b><br>REF 254 230 |    | 0.6 x 25 mm/23G | 23G |
| <b>Miraject Endotec Duo</b><br>REF 254 231 |    | 0.5 x 25 mm/25G | 25G |
| <b>Miraject Endotec Duo</b><br>REF 254 232 |  | 0.4 x 25 mm/27G | 27G |
| <b>Miraject Endotec Duo</b><br>REF 254 233 |  | 0.3 x 25 mm/30G | 30G |

## • INJECTION •

### Carpule

|                                      |                                                                                     |             |     |
|--------------------------------------|-------------------------------------------------------------------------------------|-------------|-----|
| <b>Miraject 17/42</b><br>REF 254 200 |  | 0.5 x 35 mm | 25G |
| <b>Miraject 17/23</b><br>REF 254 201 |  | 0.5 x 21 mm | 25G |
| <b>Miraject 30/42</b><br>REF 254 202 |  | 0.4 x 35 mm | 27G |
| <b>Miraject 30/23</b><br>REF 254 203 |  | 0.4 x 21 mm | 27G |
| <b>Miraject 40/23</b><br>REF 254 204 |  | 0.3 x 21 mm | 30G |
| <b>Miraject 40/10</b><br>REF 254 205 |  | 0.3 x 10 mm | 30G |

### Luer

|                                           |                                                                                       |             |     |
|-------------------------------------------|---------------------------------------------------------------------------------------|-------------|-----|
| <b>Miraject Luer 17/42</b><br>REF 254 210 |  | 0.5 x 42 mm | 25G |
| <b>Miraject Luer 17/23</b><br>REF 254 211 |  | 0.5 x 23 mm | 25G |
| <b>Miraject Luer 30/42</b><br>REF 254 212 |  | 0.4 x 42 mm | 27G |
| <b>Miraject Luer 30/23</b><br>REF 254 213 |  | 0.4 x 23 mm | 27G |

**50 years**  
of highest quality


## Miraject® P / P Super

### Canula for paste-like materials

- Blunt, bendable and individually sterile packed
- For cylinder ampoules
- Color-coded

#### Pack with 100 pieces

| Article description | ø      | Gauge | Length | REF     |
|---------------------|--------|-------|--------|---------|
| Miraject P          | 0.8 mm | 21    | 25 mm  | 254 206 |
| Miraject P Super    | 0.9 mm | 20    | 32 mm  | 254 207 |



## Miraject® PL Super

### Short blunt Luer application canula

- For application of liquids, gels and for rinsing wider root canals
- Sterile packed
- Bendable
- For syringes with Luer tip

#### Pack with 100 pieces

| Article description | ø      | Gauge | Length | REF     |
|---------------------|--------|-------|--------|---------|
| Miraject PL Super   | 0.9 mm | 20    | 22 mm  | 254 214 |



## Miraject® Endo Luer

### Long blunt Luer application canula

- Application of liquids and pastes into the root canal
- Finest controllable measures
- In two sizes: 0.3 mm (30 gauge) and 0.5 mm (25 gauge), easy-to-insert into the root canal
- For syringes with Luer tip
- Color-coded

#### Pack with 100 pieces

| Article description        | ø      | Gauge | Length | REF     |
|----------------------------|--------|-------|--------|---------|
| Miraject Endo Luer, fine   | 0.3 mm | 30    | 40 mm  | 254 222 |
| Miraject Endo Luer, medium | 0.5 mm | 25    | 40 mm  | 254 215 |



## Miraject® Endo Rinse

### Convenient rinsing system for root canals, etc.

- Consisting of disposable syringe and blunt canula, sterile packed
- Syringe: for 3 ml liquid, with 1/10 ml calibration
- Canula: blunt, fine, long and bendable canula (0.5 mm/25 G x 32 mm)

#### Pack with 50 disposable syringes and 50 canulas

| Article description | ø      | Gauge | Length | REF     |
|---------------------|--------|-------|--------|---------|
| Miraject Endo Rinse | 0.5 mm | 25    | 32 mm  | 254 220 |



50  
+ 50