

Traypurol®

Concentrated cleaning agent for impression trays and instruments



REF 2289 Bottle 1 l
REF 2290 Canister 5 l
REF 2304 Dispenser pump – canister 5 l
Not available in Canada

Indications

Self-acting cleaning of impression trays and instruments

Advantages

- Simply diluted in water
- Removes residues of alginates, zinc oxide eugenol, glass ionomer, carboxylate and phosphate cements
- Cleans trays and instruments reliably overnight
- Can be used in ultrasonic devices
- Suited also for non-oxidising metals and plastic
- Biodegradable and gentle

Traypurol® tabs

Active-tabs for cleaning impression trays and instruments



REF 2355 50 pcs.

Indications

Self-acting cleaning of impression trays and instruments

Advantages

- Removes alginates and cements within a short amount of time
- Especially easy to dose
- No vigorous stirring necessary – the tabs automatically dissolve when placed in water
- Gentle to material – suitable for all rustproof metals (including aluminium) and plastics
- Bio-degradable and pH-neutral
- Can be used in ultrasonic baths