

PalaXpress®

The precisely fitting universal denture acrylic.



PalaXpress is the material for a lot of indications with acrylic prosthetics. Its unique formula provides a wide choice of either pourable or injectable processing. Remarkable: You can pour several denture saddles in one procedure thanks to over 10 minutes pouring time.

The advantages in summary:

- Universal denture acrylic for a lot of indications
- High precision with Palajet
- Colour consistency of full and partial dentures or repairs
- Pour several saddles in one procedure
- Completion of opposite saddle without intermediate polymerisation
- Very good mucosal compatibility thanks to low amount of residual monomer
- Colour stability thanks to catalyst system
- Unique feed system of Palajet injection process
- Cadmium free

