



## Contoured Tofflemire Matrices

(0.025mm / .001" - 0.04mm / .0016")



- Easy to burnish, soft stainless steel

Contents	Height	Thickness	REF
<b>Premolars</b>			
30 pcs.	5.5 mm	0.025 mm	6109
30 pcs.	5.5 mm	0.04 mm	6107
<b>Molars</b>			
30 pcs.	6.3 mm	0.025 mm	6108
30 pcs.	6.3 mm	0.04 mm	6106



## Contoured Anatomical Matrices (0.03mm / .0012")



- Anatomically preformed and contoured
- Easy to burnish, soft stainless steel

Contents	Height	Thickness	REF
<b>Premolars</b>			
30 pcs.	4.5-5.5 mm	0.03 mm	6115
30 pcs.	6.5-7.5 mm	0.03 mm	6114
30 pcs.	4.5-7.5 mm	0.03 mm	6113
30 pcs.	4.5-7.5 mm	0.03 mm	6116
<b>Molars</b>			
30 pcs.	4.5-6.3 mm	0.03 mm	6117
30 pcs.	7.0-9.0 mm	0.03 mm	6118
30 pcs.	7.0-9.5 mm	0.03 mm	6119
30 pcs.	7.0-9.5 mm	0.03 mm	6120



## Softrix matrix strip (0.025mm / .001")



- Burnishable soft steel matrix strip

Contents	Height	Thickness	REF
3 m	6 mm	0.025 mm	5906
3 m	8 mm	0.025 mm	5908

## ST matrix strip (0.03mm / .0012")



- Stainless steel matrix strip

Contents	Height	Thickness	REF
3 m	5 mm	0.03 mm	5235
3 m	6 mm	0.03 mm	5236
3 m	7 mm	0.03 mm	5237

## Soft matrix strip (0.05mm / .002")



- Burnishable soft steel matrix strip

Contents	Height	Thickness	REF
3 m	5 mm	0.05 mm	5275
3 m	6 mm	0.05 mm	5276
3 m	7 mm	0.05 mm	5277

## Matrix strip (0.05mm / .002")



- Stainless steel matrix strip

Contents	Height	Thickness	REF
3 m	5 mm	0.05 mm	5255
3 m	6 mm	0.05 mm	5256
3 m	7 mm	0.05 mm	5257