

Gingival Cord Packers

Gingival Cord Packers are used to pack retraction cord into the sulcus. The cord packers with round, non-serrated working ends are used for atraumatic cord placement; serrated cord packers should only be used with braided cord.



Learn more about our [Gingival Cord Packers](#).



BN1 Gingival Cord Packer

#41

[GCPBN1](#)



CSI-1 Gingival Cord Packer, Non-Serrated

#41

[GCPCSI1NS](#)

CSI-1 Gingival Cord Packer, Serrated

#6

[GCPCSI16](#)

#41

[GCPCSI1](#)



113 Gingival Cord Packer, Non-Serrated

#41

[GCP113NS](#)

113 Gingival Cord Packer, Serrated

#6

[GCP1136](#)

#41

[GCP113](#)



7 Guyer Gingival Cord Packer, Non-Serrated

#6

[GCPG7NS6](#)

#41

[GCPG7NS](#)

7 Guyer Gingival Cord Packer, Serrated

#6

[GCPG76](#)

#41

[GCPG7](#)