

TECHNICAL DATA SHEET



Cotton pellets



MATERIAL

- 100 % bleached absorbent cotton
- H2O2 bleached (non chlorine bleached, free of any optical brightener)
- Meeting the European Pharmacopoeia Standards
- In accordance with REACH European Standards

PRODUCT INFORMATION

- Used to dry cavities and to apply medications
- Soft and high stability of form
- EO (Ethylene Oxide) sterilisable
- Shelf life: 5 years from production date
- Class I medical device, single use product
- General tolerance: +/- 10 %

PACKAGING

- Packed in boxes to protect the product from dust
- The goods can also be delivered in bulk
- Other packing and personalized packaging on request
- Refills dimensions: 36x36x43 mm, +/- 3 gr
- Boxes dimensions: 76x76x47 mm
- Export box dimensions: 60x40x60 cm
- Should be stored in a clean and dry place, out of direct sunlight

STANDARDS PACKAGING

Reference	Size	Quantities/box	Boxes/Export box	Export boxes/Pallet
CPA1/4	1 – Large – 7 mm	12 gr ; +/- 420 pellets	420	12
CPA2/4	2 – Medium – 5,5 mm	12 gr ; +/- 780 pellets	420	12
CPA3/4	3 – Small – 4 mm	12 gr ; +/- 1800 pellets	420	12
CPA4/4	4 – Extra small – 3 mm	12 gr ; +/- 2520 pellets	420	12

ALAN & CO

Avenue André Ernst 17 – 4800 Verviers – Belgium
Tel : +32 87 31 10 59 – Fax : +32 87 31 61 69

We are ISO 13485 certified and our products are CE marked

Version 2021/02