


SELF-SEALING STERILIZATION POUCHES



- With internal and external indicators
- Triple adhesive seal to ensure sterile storage
- Blue tinted plastic window helps detect puncture holes that may be caused by jagged instruments
- Easy to read colour steam indicator assures that sterile conditions have been achieved

6000 2 1/4" x 3 1/4" (57 mm x 90 mm) 200 pcs/box £5.95
 6001 2 1/4" x 5" (57 mm x 130 mm) 200 pcs/box £5.95
 6003 3 1/2" x 5 3/4" (90 mm x 145 mm) 200 pcs/box £6.45
 6005 2 3/4" x 9" (70 mm x 230 mm) 200 pcs/box £6.95
 6007 3 1/2" x 9" (90 mm x 230 mm) 200 pcs/box £7.95
 6009 5 1/4" x 10" (135 mm x 260 mm) 200 pcs/box £11.95
 6011 7 1/2" x 13" (190 mm x 330 mm) 200 pcs/box £23.95
 6015 10" x 15" (255 mm x 380 mm) 100 pcs/box £21.95


ACRYLIC POUCH DISPENSERS



- Smart, sturdy, durable construction
- Can be used on table top or wall mounted in either vertical or horizontal position
- Sizes:-
 061-0 (75 mm x 145 mm x 210 mm)
 061-1 (85 mm x 275 mm x 210 mm)
 061-2 (100 mm x 185 mm x 210 mm)
 061-3 (100 mm x 275 mm x 210 mm)
 061-4 (150 mm x 295 mm x 210 mm)

Pouch Dispenser
 061-0 for 2 1/4" x 5" (57 mm x 130 mm) Pouches £36.95
 061-1 for 2 3/4" x 9" (70 mm x 230 mm) Pouches £41.95
 061-2 for 3 1/2" x 5 1/4" (90 mm x 135 mm) Pouches £41.95
 061-3 for 3 1/2" x 9" (90 mm x 230 mm) Pouches £46.95
 061-4 for 5 1/4" x 10" (135 mm x 260 mm) Pouches £62.95

STERILIZATION TAPES WITH INDICATORS



- Colour changing indicators visible before and after sterilization

Indicator Tape
 1255 55 M, 1/2" (12 mm) .. £4.95
 1955 55 M, 3/4" (19 mm) .. £5.95

**Buy 2
get
1 FREE**
(same item)



DISPOSABLE FLUID RESISTANT TRAY COVER PAPERS

- Two ply of facial quality tissue plus a layer of polyethylene film for positive barrier
- Colours: Blue or White

2220 "B" Plastic Tray Covers 29.5 cm X 21.5 cm 200'S £6.45
 2221 "F" Plastic Tray Covers 20 cm X 14 cm 200'S £5.95
 2223 "ISO" Aluminium Tray Covers 27.5 cm X 17.5 cm 200'S .. £6.45

DISPOSABLE ISO TRAY COVER PAPERS for metal tray



3275 ISO Tray Cover Paper
 250 pcs/pack ... £6.95

- Liner is placed on the tray before autoclaving
- Moisture from autoclaving is absorbed by paper
- Diamond burs, finishing burs and small items can not fall through the perforated trays
- Gives the tray an attractive appearance and at the same time reduces the noise created by instruments with metal trays
- ISO size: 27.5 cm x 17.5 cm

PLASTIC STERILIZATION CASSETTES




- Sturdy, durable construction
- Autoclavable up to 134°C
- Not recommended for thermal disinfection
- Sizes:
 3881 190 mm x 90 mm x 25 mm
 3882 190 mm x 135 mm x 25 mm

3881 Small (6 instruments) £17.95
 3882 Medium (8 instruments) .. £20.95

STANDARD COLOURS	Blue	White
	Light Blue	Green
		Lavender

DISPOSABLE NON-WOVEN HEADREST COVERS



- Breathable non-woven fabric
- Soft, provides maximum patient comfort
- Size: Small 28 cm x 26 cm
 Medium 28 cm x 33 cm
- Colours: Light Blue or White

Headrest Covers, 500 pcs/box
 780S Small £27.95
 780M Medium .. £34.95

premium plus™

DISPOSABLE PROPHY BIBS EXTRA LONG

- Ideal for prophylactic, surgical procedures or where extra protection is needed
- Features contour neck design for additional coverage to shoulders and arm
- Two ply of facial quality tissue plus a layer of polyethylene
- Size: W45 cm x L60 cm
- Colour: Light Blue

2208 Prophylactic Bibs Extra Long
 100 pcs/box + Bib Holder XS ... £15.95