

Inverted cone with collar



## 806




Shank	Ring	Length mm	2,2	2,2	3,0	3,5	3,5	3,5		
		ISO	Product Nr.		Ø mm					
FG		806 314 019 534 ...	<b>G806. ...</b>			010	012	014	016	018
		806 314 019 524 ...	<b>806. ...</b>		009	010	012	014	016	018
		806 314 019 514 ...	<b>F806. ...</b>			010	012	014	016	018

Inverted cone



## 807




Shank	Ring	Length mm	3,5	3,5	3,5	4,0	5,0		
		ISO	Product Nr.		Ø mm				
FG		806 314 225 534 ...	<b>G807. ...</b>			012	014	016	018
		806 314 225 524 ...	<b>807. ...</b>		010	012	014	016	018

Barrel



## 811



Shank	Ring	Length mm	4,0	4,2	7,0	7,0		
		ISO	Product Nr.		Ø mm			
FG		806 314 038 534 ...	<b>G811. ...</b>		031	033	037	047
		806 314 038 524 ...	<b>811. ...</b>		031	033	037	047

Amalgam remover



## 813



Shank	Ring	Length mm	1,5	1,5	1,5	2,0	2,0		
		ISO	Product Nr.		Ø mm				
FG		806 314 032 534 ...	<b>G813. ...</b>		010	012	014	016	018
		806 314 032 524 ...	<b>813. ...</b>		010	012	014	016	018

811 Cavity preparation

