

endosurgery

Retro surgery



S12-70D

Retro surgery tip angled at 70°, diamond-coated 30µm, length 5mm, taper 9%

Treatment of posterior areas, in canals that are difficult to access or roots with specific orientations.



P14D

Universal retro surgery tip, diamond-coated 30µm, length 5mm, taper 7%

Preparation of canals in anterior teeth. The micro-retro tips make minimum treatment possible providing fast healing.



P15LD

Left-oriented retro surgery tip, diamond-coated 30µm, length 5mm, taper 7%

Preparation of premolar and molar canals.



P15RD

Right-oriented retro surgery tip, diamond-coated 30µm, length 5mm, taper 7%

Preparation of premolar and molar canals.



S12-70D
F00118



P14D
F00106



P15LD
F00107

P15RD
F00108