


# Dura-Green



Silicon carbide stones with a ceramic binder

Dura-Green stones are made of high-quality silicon carbide and are exceptionally well-suited for fast and splinter proof contouring and finishing of porcelain, composite, precious alloys, and amalgam.



In the dental laboratory the extensive range of shapes helps to optimally produce the individual characteristics of the surface texture of front teeth with their vertical and horizontal structures. In the dental practice Dura-Green stones are used for quick and simple contouring and finishing of restorations made of all modern filling materials.

■ Dura-Green stones can be cleaned or customised to any individual shape with the Abrasives Dresser.



## highlights

- Precise concentricity
- Low vibration performance
- Consistent cutting action
- Binder matched exactly to the grinding particles
- Long-lasting operating life
- Optimal cutting performance

### Indications

#### ■ DENTAL PRACTICE

Contouring and finishing of porcelain, composites, precious metal and amalgam

#### ■ LABORATORY

Contouring and finishing porcelain, indirect composites and metals

### Range of application

Suitable for the following materials:  
porcelain, composite, precious alloys and amalgam



Shape	CN1	CN2	CY2	CY3	FL1	FL2	FL3	FL4	FL5	IC1	IC2	IC3
ISO No.	024	037	025	028	030	025	030	045	065	030	032	050
HP 12 pcs.	0002		0006	0007	0009	0010	0011	0012	0013	0014	0015	0016
HP 72 pcs.						HG10				HG14		HG16
CA 12 pcs.	0052	0053			0059	0060	0061	0062				0066
FG 12 pcs.	0101		0103			0104						



Shape	IC4	IC5	IC6	IC7	IC9	KN1	KN3	KN4	KN5	KN7	PC1	PC2
ISO No.	045	070	065	125	135	050	070	080	060	125	030	050
HP 12 pcs.	0017	0018	0019	0020	0022	0023	0025	0026	0027	0029	0030	0031
HP 72 pcs.		HG18		HG20	HG22					HG29		
CA 12 pcs.											0080	0081



Shape	RD1	RD2	RD3	RD4	RE1	RE2	TC1	TC2	TC3	TC4	WH3	WH4
ISO No.	030	024	036	050	062	062	020	025	030	040	065	060
HP 12 pcs.	0032	0033	0034	0035	0036	0037	0039	0040	0041	0042	0045	0046
HP 72 pcs.										HG42		
CA 12 pcs.	0082	0083	0084		0086		0089	0090		0092	0095	0096
FG 12 pcs.	0107	0108					0109					



Shape	WH5	WH6
ISO No.	065	130
HP 12 pcs.	0047	0048
HP 72 pcs.		HG48
CA 12 pcs.		0098

**Min<sup>-1</sup>** Recom. 5,000-20,000 / Max. 30,000 / CP 3-5 N

PN	Quantity	Product
0002-0098	12 pcs.	Dura-Green stones, HP and CA shank
0101-0109	12 pcs.	Dura-Green stones, FG shank
HG10-HG48	72 pcs.	Dura-Green stones, HP shank

Info: Not suited for the processing of ceramic alloys!