



## Standard Luxating Elevators

Luxating elevators have thin, sharp blades that are used to sever periodontal ligaments and used in a rocking motion vs standard elevators which are used in a lifting motion. Luxating elevators reduce putting pressure on adjacent teeth.



**Luxating Elevator, Curved**  
2 mm

|              |                        |
|--------------|------------------------|
| Black Line   | <a href="#">EL2CX</a>  |
| Standard     | <a href="#">EL2C</a>   |
| Small Handle | <a href="#">EL2CSM</a> |



**Luxating Elevator, Straight**  
2 mm

|              |                        |
|--------------|------------------------|
| Black Line   | <a href="#">EL2SX</a>  |
| Standard     | <a href="#">EL2S</a>   |
| Small Handle | <a href="#">EL2SSM</a> |



**Luxating Elevator, Curved**  
3 mm

|              |                        |
|--------------|------------------------|
| Black Line   | <a href="#">EL3CX</a>  |
| Standard     | <a href="#">EL3C</a>   |
| Small Handle | <a href="#">EL3CSM</a> |



**Luxating Elevator, Straight**  
3 mm

|              |                        |
|--------------|------------------------|
| Black Line   | <a href="#">EL3SX</a>  |
| Standard     | <a href="#">EL3S</a>   |
| Small Handle | <a href="#">EL3SSM</a> |



**Luxating Elevator, Curved**  
4 mm

|              |                        |
|--------------|------------------------|
| Black Line   | <a href="#">EL4CX</a>  |
| Standard     | <a href="#">EL4C</a>   |
| Small Handle | <a href="#">EL4CSM</a> |



**Luxating Elevator, Straight**  
4 mm

|              |                        |
|--------------|------------------------|
| Black Line   | <a href="#">EL4SX</a>  |
| Standard     | <a href="#">EL4S</a>   |
| Small Handle | <a href="#">EL4SSM</a> |



**Luxating Elevator, Curved**  
5 mm

|              |                        |
|--------------|------------------------|
| Black Line   | <a href="#">EL5CX</a>  |
| Standard     | <a href="#">EL5C</a>   |
| Small Handle | <a href="#">EL5CSM</a> |



**Luxating Elevator, Curved**  
5 mm

|            |                        |
|------------|------------------------|
| Black Line | <a href="#">EL5SX</a>  |
| #510       | <a href="#">EL5S</a>   |
| Small      | <a href="#">EL5SSM</a> |

- Refer to [page 15](#) for a detailed
- listing of the [Luxating Elevator Kit](#).