























 50 W6: 40,5 LC: 8,4 WC: 7,6	 80 W6: 41,0 LC: 8,9 WC: 7,7
 51 W6: 42,9 LC: 9,2 WC: 8,1	 82 W6: 41,6 LC: 9,4 WC: 7,8
 52 W6: 43,0 LC: 9,4 WC: 8,3	 83 W6: 42,8 LC: 9,8 WC: 8,1
 53 W6: 43,5 LC: 8,8 WC: 8,6	 84 W6: 44,1 LC: 10,0 WC: 8,2
 54 W6: 45,3 LC: 9,7 WC: 8,5	 85 W6: 43,8 LC: 9,7 WC: 8,1
 55 W6: 45,1 LC: 9,0 WC: 8,4	 86 W6: 44,6 LC: 10,6 WC: 8,5
 56 W6: 44,7 LC: 9,9 WC: 8,4	 87 W6: 46,7 LC: 10,1 WC: 8,7
 57 W6: 49,3 LC: 9,5 WC: 9,1	 74 W6: 42,9 LC: 10,4 WC: 7,8
 58 W6: 45,4 LC: 10,6 WC: 8,4	 75 W6: 46,1 LC: 9,4 WC: 9,0
 69 W6: 51,1 LC: 12,5 WC: 9,7	 78 W6: 45,7 LC: 10,6 WC: 8,7


 U10 W6: 33,8 LC: 8,4 WC: 5,2	 U12 W6: 31,4 LC: 9,0 WC: 4,6
 U13 W6: 34,5 LC: 8,8 WC: 5,1	 U14 W6: 34,3 LC: 9,2 WC: 5,0
 U15 W6: 35,4 LC: 9,5 WC: 5,4	 U16 W6: 35,1 LC: 9,5 WC: 5,0
 U17 W6: 35,8 LC: 9,4 WC: 5,5	 U18 W6: 36,5 LC: 9,3 WC: 5,6
 U19 W6: 36,9 LC: 12,0 WC: 5,8	 U29 W6: 38,0 LC: 10,6 WC: 5,7

 M35 30,9	 M41 32,9
 M23 28,7	 M43 30,6
 M34 30,7	 M44 34,3
 M42 31,2	 S2 33,1
 M45 29,9	 S4 33,2
 M56 33,3	

 M41 30,1	 M43 31,6
 M44 31,2	 M44 33,0
 M44 32,5	 M44 34,2
 S2 30,3	 S4 32,6
 S4 32,0	 S6 33,8
 S6 32,6	

ARTICULACIONES ACONSEJABLES RECOMENDED ARTICULATIONS KOMBINATIONÜBERSICHT ARTICULATIONS CONSEILLÉES						
ANTERIORES ANTERIORS FRONTZÄHNE ANTÉRIEURES			POSTERIORES POSTERIORS BACKENZÄHNE POSTÉRIEURES			
SUP. UPPERS OBERE SUP.	INFERIORES LOWERS UNTERE INFÉRIEURES		M45	M41	M23	
50	U10	U12	M45	M41	M23	
51	U10	U12	M45	M41	M23	
52	U12		S2	M42	M43	
53	U13		M44	M45	S2	
54	U14	U15	M44	M42	S4	
55	U15		M41	M45	M23	
56	U16		S4	M34	M42	
57	U17		M42	M44		
58	U18		M34	S6		
69	U29	U19	S6	M56		
74	U14		M34	M35		
75	U15	U17	S4	M44	M35	
78	U18		M34	M43	S6	
80	U12		M41	M44		
82	U12		M44	M41		
83	U13	U16	M44	M41	M35	
84	U15		M43	S2	M42	
85	U15	U16	M44	S4	M42	
86	U16		M43	S6	M34	
87	U17	U18	S6	M34	M43	

LAS CIFRAS SIGNIFICAN: MEDIDAS (mm)
THE FIGURES ARE MEANING: MEASURES (mm)
DIE ZAHLEN BEDEUTEN: MASSE (mm)
LES NOMBRES SIGNIFIENT: MESURES (mm)

 84 W6: 44,1 LC: 10,0 WC: 8,2
--