

Getting Started with the Versah™ Guided Surgery System



Determine the Following:

- Tooth #
- Implant Diameter



Plan the Case:

- Take a CBCT Scan, Dicom File, and impressions of the site and the opposing jaw.
- Contact a Versah™ Certified Lab to begin the process of printing a custom C-Guide™ Sleeve.



Versah™ Certified Labs

A full list of certified labs available at versah.com/guided-surgery-system/

Dorsey Dental

www.dorseydentallab.com

Cinzara

www.cinzara.com

ROE Dental Lab

www.roedentallab.com

3D Diagnostix

www.3ddx.com

D3 Solutions

www.d3implantsolutions.com

Pro-Art Dental Laboratory

www.pro-artdentallab.com

Versah™ Guided Surgery System Price Breakdown

Item	Unit Price
G-Stop™ Vertical Gauge ...	\$24 each
G-Stop™ Key	\$74 each
C-Guide™ Sleeve.....	\$19 each
G-Stop™ Holder	\$170

I Use CEREC

- Set up your software to produce the guide with a hole that matches the Versah™ C-Guide™ Sleeve outer-diameter.
- Create buccal window to match the Versah™ C-Guide™ Sleeve opening.
- Order the C-Guide™ Sleeve and the compatible G-Stop™ from Versah™.

I Print My Own Guide

- Set up your software to produce the guide with a hole that matches the Versah™ C-Guide™ Sleeve outer-diameter.
- Create buccal window to match the Versah™ C-Guide™ Sleeve opening.
- Order the C-Guide™ Sleeve and the compatible G-Stop™ from Versah™.

I'd Like to Use My Lab

- Versah™ is committed to certify your dental lab to produce the Universal C-Guide™ Surgical Guide.

Please call us at:
517-796-3932



Versah™ Guided Surgery System Compatible Software

A full list compatible software available at versah.com/guided-surgery-system/



3shape.com



BlueSkyBio

blueskybio.com



coDiagnostiX™

codiagnostix.com

ImplaStation

en.implastation.com



R2GATE.COM

Sizing the Versah™ Guided Surgery System

G-Stop™ Vertical Gauge



Available Sizes





S
Placement up to 3.25 mm Implant

M
Placement up to 4.3 mm Implant

L
Placement up to 5.3 mm Implant

XL
Placement up to 6.2 mm Implant

Densah™ Bur Compatibility



G-Stop™ Vertical Gauge + Key



OD 4.2 mm



OD 5 mm



OD 6 mm



OD 7 mm

C-Guide™

ID 4.40 mm



6.25 mm

OD 5.3 mm




6.25 mm

ID 5.20 mm



OD 6.1 mm



ID 6.20 mm



OD 7.1 mm



ID 7.20 mm



OD 8.1 mm

Available G-Stop® Drilling Depths

3 mm | 5 mm | 8 mm | 10 mm | 11.5 mm | 13 mm | 15 mm



C-Guide™ Sleeve



WWW.VERSAH.COM

