




HU-FRIEDY MIRRORS

A BETTER VIEW FOR BETTER PERFORMANCE

How the best perform



HU-FRIEDY MIRRORS

- Front surface reflection provides a sharp, distortion-free image
- High quality stainless steel ensures durability and corrosion resistance
- Double-sided version allows clinicians the ability to maintain indirect vision and illumination while comfortably retracting the patient’s cheek and tongue
- Unique retractor shank and thin profile enhance patient comfort
- Smooth and sleek design allows for efficient cleaning and care




NEW

HIGH DEFINITION MIRRORS

Hu-Friedy’s new High Definition Mirrors feature a proprietary coating that provides superior brilliance and color. The extra brightness increases visibility in all dental procedures. Also scratch-resistant, our HD mirrors are the clear choice for optimal patient treatment.

- 113% reflection factor for exceptional image clarity*
- 38.5% brighter than rhodium coated mirrors*
- 50% brighter than other front surface mirrors*
- Scratch-resistant for product longevity
- Large ergonomic handle options



HIGH DEFINITION VS RHODIUM COATING

RHODIUM-COATED MIRRORS

SINGLE-SIDED

Three Packs

- | | |
|-----------------|--------|
| #5 Mouth Mirror | MIR5/3 |
| #4 Mouth Mirror | MIR4/3 |

Twelve Packs

- | | |
|-----------------|---------|
| #5 Mouth Mirror | MIR5/12 |
| #4 Mouth Mirror | MIR4/12 |

DOUBLE-SIDED

Three Packs

- | | |
|-----------------|----------|
| #5 Mouth Mirror | MIR5DS/3 |
| #4 Mouth Mirror | MIR4DS/3 |

Six Packs

- | | |
|-----------------|----------|
| #5 Mouth Mirror | MIR5DS/6 |
| #4 Mouth Mirror | MIR4DS/6 |

SINGLE-SIDED

Single

- | | |
|--------------------|--------|
| #3 HD Mouth Mirror | MIR3HD |
| #4 HD Mouth Mirror | MIR4HD |
| #5 HD Mouth Mirror | MIR5HD |

Three Packs

- | | |
|--------------------|---------|
| #4 HD Mouth Mirror | MIR4HD3 |
| #5 HD Mouth Mirror | MIR5HD3 |

Twelve Packs

- | | |
|--------------------|----------|
| #4 HD Mouth Mirror | MIR4HD12 |
| #5 HD Mouth Mirror | MIR5HD12 |

DOUBLE-SIDED


Single

- | | |
|--------------------|----------|
| #4 HD Mouth Mirror | MIR4DSHD |
| #5 HD Mouth Mirror | MIR5DSHD |

Six Packs

- | | |
|--------------------|-----------|
| #4 HD Mouth Mirror | MIR4DSHD6 |
| #5 HD Mouth Mirror | MIR5DSHD6 |

*Data on file.



MIRROR HANDLE OPTIONS



#24 Cone Socket
Mirror Handle

| MH24



#1 Cone Socket
Mirror Handle

| MH1



#6 Cone Socket
Mirror Handle

| MH6



#7 Cone Socket
Mirror Handle

| MH7



#8 Cone Socket
Mirror Handle

| MH8

For more information about
Hu-Friedy's Total Solutions visit

HU-FRIEDY.COM


PERIODONTAL




RESTORATIVE



ORTHODONTIC




INSTRUMENT
MANAGEMENT



ENDODONTIC



SURGICAL



ULTRASONIC
SCALING



Hu-Friedy Mfg. Co., LLC: 3232 N. Rockwell St. | Chicago, IL 60618 | USA
Telephone: 1-800-Hu-Friedy (1-773-975-6100) | Website: Hu-Friedy.com

©2014 Hu-Friedy Mfg. Co., LLC. All rights reserved. HF-409/0914

