

## Provicol® · Provicol® QM

Eugenol-free temporary luting cement with calcium hydroxide

### Indications

Temporary luting of inlays, onlays, crowns, partial crowns and bridges  
Temporary filling of small one-surface cavities

### Advantages

- Semi-elastic properties for safe adhesion of the temporary
- Contains calcium hydroxide which is known to support the formation of tertiary dentine and to have an anti-microbial effect
- Ideal before adhesive luting with resin cements, because no inhibition of curing reaction of resin-based filling or luting materials



SC self-curing    auto mix    temporary    contains Ca(OH)<sub>2</sub>

REF 1075	Tube 2 × 25 g (base / catalyst)
REF 1074	QuickMix syringe 5 ml, accessories
REF 2261	Dispenser – QuickMix
REF 2189	Mixing tips type 10, 50 pcs.
<i>REF 1075 is not available in Canada</i>	

## Provicol® C

Proven Provicol in cartridges

### Indications

Temporary luting of inlays, onlays, crowns, partial crowns, bridges, attachments and posts  
Temporary filling of small one-surface cavities

### Advantages

- Semi-elastic properties for safe adhesion of the temporary
- Contains calcium hydroxide which is known to support the formation of tertiary dentine and to have an anti-microbial effect
- Eugenol-free, thus low allergy potential
- Ideal before adhesive luting with resin cements, because no inhibition of curing reaction of resin-based filling or luting materials



SC self-curing    auto mix    temporary    contains Ca(OH)<sub>2</sub>

REF 1076	Cartridge 2 × 65 g, accessories
REF 2158	Dispenser – type 2
REF 2185	Mixing tips type 6, 50 pcs.

## Bifix® Temp

Dual-curing, composite-based temporary luting material

### Indications

Non-permanent luting of temporary restorations

### Advantages

- **Excellent aesthetics:** translucent, tooth-like shade
- **Handling:** film thickness of less than 5 µm for a perfectly accurate fit
- **Dual-curing:** simple removal of excess material after short exposure to light
- **Fluorescent under UV-A light:** for easier detection and removal of even the smallest material residues



DC dual-curing    auto mix    temporary    fluorescent

REF 1055	QuickMix syringe 5 ml universal, accessories
REF 1056	QuickMix syringe 5 ml fluorescent, accessories <b>NEW</b>
REF 2261	Dispenser – QuickMix
REF 2188	Mixing tips type 9, 50 pcs.



CAD / CAM material

Luting

Core build-up

Temporaries

Denture relining

Impression

Repair