



Standard Luxating Elevators

Luxating elevators have thin, sharp blades that are used to sever periodontal ligaments and used in a rocking motion vs standard elevators which are used in a lifting motion. Luxating elevators reduce putting pressure on adjacent teeth.



Luxating Elevator, Curved 2 mm	
Black Line	EL2CX
Standard	EL2C
Small Handle	EL2CSM



Luxating Elevator, Straight 2 mm	
Black Line	EL2SX
Standard	EL2S
Small Handle	EL2SSM



Luxating Elevator, Curved 3 mm	
Black Line	EL3CX
Standard	EL3C
Small Handle	EL3CSM



Luxating Elevator, Straight 3 mm	
Black Line	EL3SX
Standard	EL3S
Small Handle	EL3SSM



Luxating Elevator, Curved 4 mm	
Black Line	EL4CX
Standard	EL4C
Small Handle	EL4CSM



Luxating Elevator, Straight 4 mm	
Black Line	EL4SX
Standard	EL4S
Small Handle	EL4SSM



Luxating Elevator, Curved 5 mm	
Black Line	EL5CX
Standard	EL5C
Small Handle	EL5CSM



Luxating Elevator, Curved 5 mm	
Black Line	EL5SX
#510	EL5S
Small	EL5SSM

• Refer to [page A707](#) for a detailed listing of the [Luxating Elevator Kit](#).