

ORSING

AUTOCLAVABLE LINE

Autoclavable products for operations, implants, root fillings and fine endodontic work.



Hygo® Surge

Hygo Surge, an autoclavable surgical aspirator tip to connect to an 11 mm aspirator tube. The surgical tip has an apex of 4,0 mm, perfect for narrow operations, root fillings and implants. Single use.

- Autoclavable
- For surgical treatments
- Precision tip, 4,0 mm

Hygo® Surge Microtip

Hygo Surge Microtip, a tiny blue autoclavable tip that turns Hygo Surge and Hygo Surge XL into unique suction devices for micro-oral surgery and fine endodontic work. Single use.

- Autoclavable
- For surgical treatments
- Precision tip, 1,7 mm



Hygo® Surge XL

Hygo Surge XL, an autoclavable surgical aspirator tip, fit 11 mm evacuation system. The surgical tip has an apex of 3,6 mm, perfect for narrow operations, root fillings and implants. Single use.

- Autoclavable
- For surgical treatments
- Precision tip, 3,6 mm



Hygovac® Green


Hygovac Green, an autoclavable aspirator tube with one traditional 45° and one S-shaped tip. Fit 11 mm evacuation system. Length 140 mm. Single use.

- Autoclavable
- Dual tips, 45° and S-shaped

Hygovac® Vent Green

Hygovac Vent Green, an autoclavable aspirator tube with one traditional 45° and one S-shaped tip. Ventilated, to make it more pleasant for the patient. Fit 11 mm evacuation system. Length 140 mm. Single use.

- Autoclavable
- Dual tips, 45° and S-shaped
- Ventilated, to avoid getting stuck



ORSING

The future of your Dental Practice

J. H. Orsing AB, a manufacturer of dental supplies focusing on a sustainable future with bio-based fossil-free products. As a leading manufacturer, waiting for a fossil-free autoclavable material for our aspirator tips and tubes, we are proud of our Autoclavable line, in PP.

Important

- Autoclavable line manufactured in autoclavable PP.
- Autoclavable up to 134°C
- No phthalates, no latex
- Products for operations, implants, root fillings and endodontic work.
- Hygo Surge Microtip, for micro oral surgery and fine endodontic work.



ORSING

J.H. Orsing AB
+46 42 295 500
orsing@orsing.se
www.orsing.se

DirectaDentalGroup

DIRECTA Iopdental **ORSING** parkell KOHLER RÖNVIG polydentia
TrollDental Sendoline **physics** CONTACTZ ceramir

