

Bifix® QM

Dual-curing universal composite-based luting system



Also for zirconium dioxide luting



Indications

Adhesive luting of inlays, onlays, veneers, crowns, partial crowns, bridges, root posts and adhesive bridges (Maryland bridges)

Advantages

- For all materials (metal, composite inlays and ceramic, also zirconium dioxide)
- Bubble-free perfect mixing
- Direct application
- Excellent adhesion to dentine, enamel and to metals
- Special one-step Ceramic Bond for optimal adhesion to porcelain
- Radiopaque
- Available in three shades
- Available as an all-in-one set

20 sec
light-curing

DC
dual-curing

auto
mix

3
shades

REF 1217	Set QuickMix syringe 10 g U (universal), Futurabond DC bottle 4 ml each of liquid 1 and 2, Ceramic Bond bottle 5 ml, accessories
REF 1218	QuickMix syringe 10 g universal, accessories
REF 1219	QuickMix syringe 10 g transparent, accessories
REF 1759	QuickMix syringe 10 g white-opaque, accessories
REF 2261	Dispenser – QuickMix

For further accessories see p. 95/96

Bifix® Hybrid Abutment

Universal, dual-curing luting composite for intraoral and extraoral application

With highly opaque shades



Indications

Permanent extraoral and intraoral luting of superstructures (incl. mesostructures) made of metal, ceramics, zirconium dioxide and polymers (composite, hybrid ceramics or PMMA) on

- prefabricated adhesive/titanium bases (extraorally only)
- customised titanium and zirconium dioxide bases (both extraorally and intraorally)
- dental hard tissue

Advantages

- All-in-one – permanent luting of all known dental materials
- Highest flexibility – intraoral application in the practice and extraoral in the laboratory
- High aesthetic – the perfect masking of dark substrates with highly opaque shades
- With Try-In pastes – precise shade simulation for hybrid abutments, abutment crowns and core treatments

DC
dual-curing

auto
mix

3
shades

REF 2410	Implant Set – QuickMix syringe 10 g (universal HO), Try-In syringe 2.5 g, Ceramic Bond bottle 5 ml, Futurabond U <i>SingleDose</i> 5 pcs., accessories
REF 2411	QuickMix syringe 10 g (universal HO), Try-In syringe 2.5 g, accessories
REF 2412	QuickMix syringe 10 g (white HO), Try-In syringe 2.5 g, accessories
REF 2413	QuickMix syringe 10 g (translucent), Try-In syringe 2.5 g, accessories
REF 2261	Dispenser – QuickMix

For further accessories see p. 95/96

CAD / CAM material

Luting

Core build-up

Temporaries

Denture relining

Impression

Repair