



TOP SAFE MAXI PATIENT TOWELS-2 LAYERS

DISPOSABLE DENTAL NAPKINS TO PROTECT THE PATIENT



CONFIGURATION OVERVIEW

Type: Non-Sterile
Latex free

Features: The special features of the materials offer high absorbency and waterproofing, in order to assure a high comfort to the patient.

TECHNICAL SPECIFICATIONS

1 ply of pure wadding of long fiber cellulose 25 g/m² (TISSUE) partially laminated with a 1 ply of high-density polyethylene 8µm . The napkins are embossed and folded.

Dimensions: 40x45 cm (tolerance ± 2%)

PRODUCT VARIATIONS

Art-No.	PATEW15-AKZ-FU-CA2	PATEW15-AKZ-FG-CA2	PATEW15-AKZ-MB-C35	PATEW15-AKZ-MC-CA2
Color	Fuchsia	Fresh Green	Med Blue	Magic Blue
	●	●	●	●

PACKAGING

Pieces/carton: 500
Carton/size: 41x24x25 cm
Carton weight: 3.5 kg
Carton/pallet: 63
Carton/layer: 7
Layers/pallet: 9
Pallet size: 80x120x245 cm
Pallet weight: 235.5 kg

STORAGE

In a dry place at room temperature, away from heat sources, flames, or sparks