

Hi-Tray Light Edentulous

Impression trays in plastic

Impression systems / Accessories

The **Hi-Tray Edentulous** are impression trays in green rigid plastic for edentulous applications.

Characteristics

- Constructed in green rigid plastic
- Available in the following versions:
 - Small upper and lower
 - Medium upper and lower
 - Large upper and lower

Advantages

- Disposable
- Adaptable
- Rigid
- Perforated design for improved mechanical retention

Standard packs 12 pieces available in the following versions:

Size	Small	Medium	Large
Upper	D5ESUP	D5EMUP	D5ELUP
Lower	D5ESLO	D5EMLO	D5ELLO



Hi-Tray Light Edentulous upper



Hi-Tray Light Edentulous lower